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Roti®-Load 3 (LDS)

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Instructions for use

Roti®-Load 3 (LDS)

Protein loading buffer, non reducing
(4x conc.)

- Loading buffer based on LDS
- Optimised for ready-made gels
- Superior solubility at low temperatures

Introduction:

Roti®-Load 3 (LDS) is a modified, non-reducing protein sample-gel loading buffer for gels under denaturing or native conditions. Roti®-Load 3 (LDS) is based on a triethanolamine buffered ready-to-use solution, in which the standard sodium dodecylsulphate (SDS) has been replaced by lithium dodecylsulphate (LDS). At low temperatures, LDS is much better soluble than SDS. Hence, no crystallisation of LDS takes place, and protein stability as well as band resolution is significantly enhanced. Additionally, the formulation of Roti®-Load 3 (LDS) has been optimised for ready-made gels, particularly Roti®-PAGE gels, but also ready-made gels of other standard suppliers.

Roti®-Load 3 (LDS) prevents degradation of proteins during sample preparation and electrophoresis run. The triethanolamine buffer system keeps the sample pH value constant over a wide temperature range, additionally stabilising the sensitive aspartyl-prolyl peptide bonds.

Use:

1. Prepare your polyacrylamide gel.
2. Determine sample volume.
3. Add Roti®-Load 3 (one third of sample volume).
Example:
60 µl sample + 20 µl Roti®-Load 3 (4x conc.)
4. Mix by pipetting.
5. Heat samples for 3-5 mins. to 80-90 °C.
6. Directly load samples to the gel.

Formulation:

LDS (approx. 4 % w/v), glycerol (approx. 40 % v/v), Ficoll® 400 (approx. 4 % w/v), phenol red (approx. 0,025 % w/v), Coomassie® G250 (approx. 0,025 % w/v).

Triethanolamine buffered.

Store at +20 °C

Tip:

For reducing conditions, Roti®-Load 3 (LDS) may be supplemented with 20 µl/ml DTT solution (1 M) or with β-mercaptoethanol.

Safety Information:

 **Warning** H315-H319

P280-P302+P352-P305+P351+P338-P337+P313

Roti®-Load 3 (LDS) 3359.1 10 ml

Loading buffers based on SDS:

Roti®-Load 1 (reducing, Art. No: K929)

Roti®-Load 2 (non reducing, Art. No: K930)

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Postfach 100121

76231 Karlsruhe

Telefon: +49 (0) 721/ 5606-0

Telefax: +49 (0) 721/ 5606-149

E-Mail: info@carlroth.de

Internet: www.carlroth.de

s.t. 04/2017