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## Roti®-Load DNA-small (mit Glycerin)

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
## Roti®-Load DNA small

(with glycerol)

6 x gel loading buffer for DNA-electrophoresis of small fragments

### I. Introduction:

Roti®-Load DNA small is a special 6x concentrated gel loading buffer, optimised for DNA electrophoresis of small fragments. It prevents DNA bands of small fragments from being concealed by bromophenol blue, even in high concentrated agarose gels. Each batch is tested for functionality in electrophoresis and for being DNase-free in order to guarantee the highest safety and reproducible results.

### II. Composition:

Roti®-Load DNA small (with glycerol) is composed as follows:

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The density of your DNA-solution is increased by the proportion of glycerol share in Roti®-Load which guarantees even distribution in the gel pocket base. Xylene cyanol blue alleviates sample application and can be used as a run marker. Bromophenol blue is omitted in order to enable clear depiction of very small DNA fragments. By adding Tris and EDTA, DNA stability is increased.

### III. Application:

Roti®-Load DNA small can either be added to DNA-solution or the DNA can be dissolved directly in Roti®-Load DNA small.

#### Dissolving DNA in 1 x Roti®-Load DNA small:

Dilute Roti®-Load DNA small 1:6 by pipetting 170 µl Roti®-Load DNA into a clean, sterile test tube. Add 830 µl bidistilled, autoclaved water and mix by vortexing briefly. This solution (Roti®-Load DNA small 1x) can now be used to dissolve precipitated DNA and can be applied directly onto the gel.

#### Adding Roti®-Load to a DNA-solution:

Volume DNA-solution		Volume Roti®-Load		Total volume Sample
5 µl	+	1 µl	=	6 µl
10 µl	+	2 µl	=	12 µl
15 µl	+	3 µl	=	18 µl
20 µl	+	4 µl	=	24 µl
25 µl	+	5 µl	=	30 µl

For example, add 1 µl Roti®-Load DNA small to a DNA solution with a volume of 5 µl. Mix the sample by carefully pipetting up and down three times.

#### OPTIONAL:

Irrespective of which method you use, we recommend denaturing the sample for 5 mins at 65 °C and then cooling on ice, otherwise run artefacts may develop (e.g. EcoRI/HindIII-separated Lambda-DNA).

### IV. Storage:

Buffer should be stored at -20 °C for long-term storage (> 3 months). Buffer currently in use can be stored for some weeks at 4 °C (or room temperature).

### V. Quality assurance:

Every Roti®-Load DNA small batch is tested for being DNase-free and for functionality in electrophoresis. For this purpose 1 µg linearized plasmid-DNA is incubated for 48 h at 37 °C in 1x DNA-buffer and the integrity of the DNA electrophoretically compared with a fresh sample.

#### **Roti®-Load DNA small**

**(with glycerol)**      HP03.1      5x1.8 ml

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#### **Carl Roth GmbH+Co. KG**