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Roti®-ImmunoBlock

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Instructions for use

Roti®-ImmunoBlock

Blocking reagent for immunofluorescent applications

A. Introduction

Roti®-ImmunoBlock is an innovative, ready to use blocking reagent made on a polymer basis without detergent and following advantages:

- Your detection system is not contaminated by any foreign proteins. The signal intensity remains constant.
- Specific binding sites are preserved whilst nonspecific reactions are suppressed. Signal intensity increased.
- The ready to use solution saves you time and labour by eliminating the need for dissolving.

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B. Application

Roti®-Block working solution

Dilute Roti®-ImmunoBlock (10 x concentrate) before use 1:10 with H₂O (distilled). The diluted solution can be used directly.

Roti®-ImmunoBlock and samples

When blocking please cover your entire sample with Roti®-ImmunoBlock working solution.

Blocking time and temperature

Blocking should occur under gentle shaking at room temperature. A blocking period from 15 mins to 1 h has proved successful.

Roti®-ImmunoBlock and antibodies

You can also incubate your sample with primary and secondary antibody diluted in Roti®-ImmunoBlock working solution. Alternatively incubation can occur in Roti®-ImmunoBlock working solution diluted 1:10 in PBS or TBS, respectively, or PBST or TBST, respectively.¹⁴

TBST: Roti®-Stock 10 x TBST (Art. No. 1061.1)
(25 mM Tris, pH 7.4, 150 mM NaCl, 3 mM KCl, 0.05% Tween 20)

PBST: Roti®-Stock 10 x PBST (Art. No. 1059.1)
(137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄, pH 7.4, 0.05 % Tween 20)

TBS: Roti®-Stock 10 x TBS (Art. No. 1060.1)
(25 mM Tris, pH 7.4, 150 mM NaCl, 3 mM KCl)

PBS: Roti®-Stock 10 x PBS (Art. No. 1058.1)
(137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄, pH 7.4)

Roti®-ImmunoBlock T144.1 250 ml

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Carl Roth GmbH+Co. KG