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## Roti®-Load DNASTain

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# Instructions for use

## Roti®-Load DNASTain

Ready-to-use, DNase-free, non-toxic  
6x DNA-gel loading buffer with fluorescent  
DNA-dye

### I. Introduction:

Versatile reagent for gel loading and simultaneous staining of double stranded nucleic acid in agarose and polyacrylamide gels. Roti®-Load DNASTain detects up to 0.2 ng per band, and is, therefore, as sensitive as ethidium bromide.

### II. Composition:

Roti®-Load DNASTain is composed as follows:

TRIS/HCl pH 7,5	60 mM
EDTA	60 mM
Glycerol	60 % (w/v)

DNASTain 1: Bromphenolblau, Xylencyanolblau,  
EvaGreen®  
DNASTain 2: Xylencyanolblau, Orange G,  
EvaGreen®  
DNASTain 3: Orange G, EvaGreen®

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
s.t. 04/2017

The density of your DNA-solution is increased by the proportion of glycerol share in Roti®-Load which guarantees even distribution in the gel pocket base. Bromophenol blue, orange G, respectively, and xylene cyanole blue alleviate sample application and can be used as a run marker. By adding Tris and EDTA, DNA stability is increased. EvaGreen® is a non-toxic fluorescence dye that stains nucleic acids during gel run.

### III. Staining:

This innovative reagent contains the dye EvaGreen®, which has been developed specifically for DNA analysis. DNA integrity is not influenced by EvaGreen®; hence, gel eluted DNA may directly be applied to standard down-stream applications, like, e.g., ligations, PCR, transformation, transfection etc. without interfering effects. DNA binding of EvaGreen® is based on the 'release-on-demand' binding technology, where (unstained) dimers are only split into (fluorescent) monomers when binding to DNA. Hence, gel background is very low. Roti®-Load DNASTain can be excited by UV- and blue light, making it fully compatible with standard UV-transilluminators as well as with modern blue-light illumination (see also Rotiphorese®-PROfessional runVIEW, Art. No. 4849.1). When excited, nucleic acid-bound Roti®-Load DNASTain emits a brightly coloured fluorescence that may be documented by standard ethidium bromide- or SYBR® Green foto filters.

Excitation maxima (bound to DNA, black curve):  
ca. 320 nm (UV-light) and 470-500 nm (blue light)  
Emission maximum (bound to DNA, green curve):  
525 nm

Suitable for staining of dsDNA and dsRNA of the following fragment lengths:  
Roti®-Load DNASTain 1: > 500bp  
Roti®-Load DNASTain 2: 100 to 2000 bp  
Roti®-Load DNASTain 3: < 500 bp

### IV. Safe use:

Although it is an intercalating dye, due to the low membrane permeability of EvaGreen®, Roti®-Load DNASTain is non mutagenic, non toxic, and environment friendly. For detailed data see our Technical Info Brochure on EvaGreen® online with the product.

### V. Application:

Roti®-Load DNASTain can either be added to DNA-solution or the DNA can be dissolved directly in Roti®-Load DNASTain.

#### Dissolving DNA in 1 x Roti®-Load DNASTain:

Dilute Roti®-Load DNASTain 1:6 by pipetting 170 µl Roti®-Load DNASTain into a clean, sterile test tube. Add 830 µl bidistilled, autoclaved water and mix by vortexing briefly. This solution (Roti®-Load DNASTain 1x) can now be used to dissolve precipitated DNA and can be applied directly onto the gel.

#### Adding Roti®-Load DNASTain to a DNA-solution:

Add one part of Roti®-Load DNASTain to 5 parts of the DNA sample as shown in the table below.

Vortex gently or pipet a few times, and spin down the tubes to achieve homogeneity.

Volume DNA-solution		Volume Roti®-Load DNASTain		Total volume
5 µl	+	1 µl	=	6 µl
10 µl	+	2 µl	=	12 µl
15 µl	+	3 µl	=	18 µl
20 µl	+	4 µl	=	24 µl
25 µl	+	5 µl	=	30 µl

#### OPTIONAL:

In order to avoid run artefacts, in some cases (e.g. *EcoRI/Hind III*-digested Lambda-DNA) denaturing the sample prior to the run may be recommended. Incubate for 5 mins. at 65 °C and then cool on ice.

#### **VI. Storage:**

For long-term storage, Roti®-Load DNASTain should be stored at -20 °C (> 3 months). Buffer currently in use can be stored for some weeks at +4 °C. Keep protected from light.

#### **VII. Quality assurance:**

All components of the buffer are proved to be DNase-free. The production is carried out under sterile conditions.

#### **VIII. Frequently Asked Questions:**

Q: Are there any samples of Roti®-Load DNASTain?

A: No, sorry. By tube-internal evaporation, the quality of samples of less than 50 µl decreases

rapidly. You may get some interesting test discount, however.

Q: How should I visualise the gels after staining with Roti®-Load DNASTain?

A: Gels can be visualised under UV-light or on a blue-light transilluminator.

Q: On UV-table the Roti®-Load DNASTain stained bands are very bright, but on photos there is virtually no band visible. What is wrong?

A: Some ethidium bromide photo filters are very narrow banded, efficiently excluding other wave lengths like the green range needed here. Use broad-band photo-filters or filters for SYBR® Green.

Q: Can Roti®-Load DNASTain be used for post-run staining or addition to the agarose?

A: No, these applications are not possible.

Q: What is the shelf life of Roti®-Load DNASTain?

A: Roti®-Load DNASTain can be kept for 1 year under cooling.

Q: How should I dispose of Roti®-Load DNASTain?

A: You should follow any existing departmental or company guidance on the disposal of reagents as indicated in the MSDS of Roti®-Load DNASTain. Roti®-Load DNASTain should not be disposed of down the sink.

Q: Does Roti®-Load DNASTain intercalate between base pair stacks?

A: Yes, the molecule intercalates.

Q: Then, why is Roti®-Load DNASTain still non-toxic?

A: Due to the very low membrane permeability, it just doesn't enter the cells. Hence, all respective tests are nearly negative.

Q: May I also use Roti®-Load DNASTain for staining of PAGE gels?

A: Yes, you can use it also for loading of DNA samples into the wells of PAGE gels.

Q: Does Roti®-Load DNASTain also efficiently stain supercoiled plasmids?

A: Yes, staining efficiency of supercoiled plasmids is as good as found for ethidium bromide.

Q: What if the bands on the gel are too faint?

A: Increase the added volume of Roti®-Load DNASTain 2times.

Q: Does Roti®-Load DNASTain also efficiently stain RNA?

A: You need double strands of nucleic acid for binding of the dye. RNA with long hairpins or other double strands is stained; RNA with long single stranded structure may be stained less efficiently.

Q: May gels stained with Roti®-Load DNASTain be used for Southern blotting and detection?

A: Yes, without reservation. During blotting and denaturation the dye simply diffuses from the DNA.

EvaGreen® made by Biotium, Inc., Hayward, CA, USA

Roti®-Load DNASTain is distributed for laboratory research only. Not for diagnostic use.

<b>Roti®-Load DNASTain 1</b>	5783.1	1 ml
	5783.2	5x1 ml
<b>Roti®-Load DNASTain 2</b>	5784.1	1 ml
	5784.2	5x1 ml
<b>Roti®-Load DNASTain 3</b>	6472.1	1 ml
	6472.2	5x1 ml