



# Instructions for use

## Roti®-Bleach Silver Destaining Kit

Easy decolorising of overstained silver gels in foil-supported or free PAA gels. The new staining consistency is firmly preserved.

For film supported and conventional PAA gels  
Sufficient for  
10 gels with a maximum volume of 10 ml or  
5 gels with a maximum volume of 40 ml

### I.1 Technical information

- Please wear gloves.
- Use double distilled, deionized water.
- Please ensure that the temperature of the required solutions is approximately 20 °C.

### I.2 Kit components

10 x sodium thiosulphate in glass vials labelled:

**A** sodium thiosulphate (N767)

10 x potassium hexacyanoferrate in test tubes

labelled: **B** potassium hexacyanoferrate (N768)

*Contents of this Kit may not be bought separately.*

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Postfach 100121

76231 Karlsruhe

Telefon: +49 (0) 721/ 5606-0

Telefax: +49 (0) 721/ 5606-149

E-Mail: info@carlroth.de

Internet: www.carlroth.de

s.t. 04/2017

### I.3 Additional chemicals required

water, double distilled, Art.No. 3478.1

acetic acid 100 % p.a., Art.No. 3738.1

ethanol p.a., Art.No. 9065.1

or ethanol, denaturated, Art.No. K928.1

Rotisol®, Art.No. 7917.1

## II Destaining of silver stained proteins and DNA

### Attention:

Please carry out the various staining steps under gentle shaking. The following procedure is appropriate for both **film supported** and **conventional** gels.

### II.1 With a maximum gel volume of 10 ml (e.g. 10 x 10 cm x 1 mm).

#### A) Washing:

When silver staining is finished (before the drying step) your gel is located in stopping or washing solution. Discard solution and wash your gel for 2 x 10 min in 100 ml H<sub>2</sub>O<sub>dd</sub>.

#### B) Preparation of Solutions:

##### a) Destaining sol.

Add the content of an „A-labelled“ vial to a beaker and dissolve the sodium thiosulphate in 100 ml H<sub>2</sub>O<sub>dd</sub> while stirring. Centrifuge briefly a „B-labelled“ test tube and pipette 1 ml H<sub>2</sub>O<sub>dd</sub> to it. Dissolve the potassium cyanoferrate completely by

vortexing. Now add the content of test tube B to solution A and mix thoroughly.

##### b) Stopping solution (300 ml):

45 ml ethanol abs.

+ 36 ml acetic acid 100 %

Adjust to a volume of 300 ml with H<sub>2</sub>O<sub>dd</sub>.

##### c) Drying solution (100 ml):

10 ml glycerol

+15 ml ethanol

Adjust to a volume of 100 ml with H<sub>2</sub>O<sub>dd</sub>.

### C) Destaining and Washing:

Overlay the gel with the freshly prepared solution and shake gently. Watch the destaining process and discard the solution (right before) just **before** you have reached an optimal result. Immediately wash your gel with plenty of water. Thereafter wash your gel 3 x 10 min in 100 ml stopping solution to stop the reaction completely. While washing the yellow coloured background disappears.

### D) Drying:

Incubate the gel  $\geq 30$  min in the drying solution. The gel can then be dried between two sheets of cellophane film. We recommend using our drying frame (ref. K420.1, K421.1).

## II.2 With a maximum gel volume of 40 ml (e.g. 20 cm x 20 x 1 mm).

### A) Washing:

When silver staining is finished (before the drying step) your gel is located in stopping or washing solution. Discard solution and wash your gel for 2 x 10 min in 100 ml H<sub>2</sub>O<sub>dd</sub>.

### B) Preparation of solutions:

#### a) Destaining sol.

Add the content of two „A-labelled“ glass vial to a beaker and dissolve the sodium thiosulphate in 100 ml H<sub>2</sub>O<sub>dd</sub> while stirring. Centrifuge briefly two „B-labelled“ test tube and pipette 1 ml H<sub>2</sub>O<sub>dd</sub> to each. Dissolve the potassium cyanoferrate completely by vortexing. Now add the content of both test tubes B to solution A and mix thoroughly.

#### b) Stopping solution (600 ml):

90 ml ethanol abs.  
+ 72 ml acetic acid 100 %  
Adjust to a volume of 600 ml with H<sub>2</sub>O<sub>dd</sub>.

#### c) Drying solution (200 ml):

20 ml glycerol  
+30 ml ethanol  
Adjust to a volume of 200 ml with H<sub>2</sub>O<sub>dd</sub>.

### C) Destaining + Washing:

Overlay the gel with the freshly prepared solution and shake gently. Watch the destaining process and discard the solution (right before) just **before** you have reached an optimal result. Immediately wash your gel with plenty of water. Thereafter wash your gel 3 x 10 min in 200 ml stopping solution to stop the reaction completely. While washing the yellow coloured background disappears.

### D) Drying:

Incubate the gel at least 30 mins. in the drying solution. The gel can then be dried between two sheets of cellophane film.

## III. Further products:

### Rotilabo®-staining chambers

(Art. Nos. H591.1, H592.1).

**Gel Drying Frames** (Art. Nos. K420.1, K421.1).

### Roti®-Black P

silver staining kit for **proteins**

L533.1 for 10 gels with a maximum vol. of 10 ml

L533.2 für 10 gels with a maximum vol. of 40 ml

### Roti®-Black N

silver staining kit for **DNA**

N769.1 for 10 gels with a maximum vol. of 10 ml

N769.2 for 10 gels with a maximum vol. of 40 ml

### Roti®-Black NSeq

silver staining kit for **DNA in sequencing gels**

P081.1 1 Kit



**Roti®-Bleach**

**N766.1**