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Instructions for use

ROTI®Load DNA (with glycerol)

1 x gel loading buffer for DNA-electrophoresis

I. Introduction:

ROTI®Load DNA is a 1x concentrated gel loading buffer for DNA electrophoresis. Each batch is tested for functionality in electrophoresis and for being DNase-free in order to guarantee the highest safety and reproducible results.

II. Composition:

ROTI®Load DNA (with glycerol) is composed as follows:

TRIS/HCl pH 7,5	10 mM
Na-acetate	5 mM
EDTA	2 mM
Glycerol	10 % (w/v)
Bromphenol blue	0.02 % (w/v)
Xylene cyanol blue	0.015 % (w/v)

The density of your DNA-solution is increased by the proportion of glycerol share in ROTI®Load which guarantees even distribution in the gel pocket base. Bromophenol and xylene cyanol blue alleviate sample application and can be used as a run marker. By adding Tris and EDTA, DNA stability is increased.

III. Application:

Precipitated DNA can be dissolved directly in ROTI®Load DNA short-run 1x. In order to accelerate solubilisation of DNA, the tube can be heated to 37 °C for 5 minutes. Exceptionally, the sample buffer may be diluted up to 0.6 fold.

OPTIONAL:

Irrespective of which method you use, we recommend denaturing the sample for 5 mins at 65 °C and then cooling on ice, otherwise run artefacts may develop (e.g. *EcoRI/HindIII*-separated Lambda-DNA).

IV. Storage:

The buffer should be stored at –20 °C for long-term storage (> 3 months). Buffer currently in use can be stored for some weeks at 4 °C (or room temperature).

V. Quality assurance:

Every ROTI®Load DNA short-run 1x batch is tested for being DNase-free and for functionality in electrophoresis.

For this purpose 1 µg linearized plasmid-DNA is incubated for 48 h at 37 °C in 1x DNA-buffer and the integrity of the DNA electrophoretically compared with a fresh sample.

1 x ROTI®Load DNA

(with glycerol)

0100.1 5x1 ml

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.