## Mowiol 4-88

Art.-Nr. 0713

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## Mowiol 4-88

Art. No. 0713

Mowiol 4-88 consists of colourless to light yellow granules. It is a partly hydrolysed polyvinyl alcohol with high adhesive power, which is principally used in industry working with adhesives.

In histology Mowiol 4-88 is used as hydrophilic mounting medium. Usually it is applied in combination with the anti-bleaching reagent DABCO (1.4-Diazabicyclo-(2.2.2)octane), which conserves the colours and the fluorescence of durable preparations.

### Instruction for hydrophilic embedding with Mowiol 4-88

#### Ingredients:

- 6.0 g Glycerol (Art. No. 3783)

- 2.4 g Mowiol 4-88
- 6.0 ml distilled water, sterile (Art. No. T143)
- 12.0 ml 0.2M Tris-HCl (pH 8.5) (Art. No. 9090)
- 25 mg per ml solution DABCO (Art. No. 0718) **(optional)**

#### Preparation:


1. Mix 6.0 g Glycerol with 2.4 g Mowiol and dissolve with frequent agitation for one hour at room temperature (RT).
2. Add 6.0 ml water and stir for one more hour at RT.
3. Add 12.0 ml 0.2 M Tris-HCl (pH 8.5) and incubate for two hours at 50 °C under periodical stirring (every 20 minutes for 2 minutes).
4. Note: In many cases Mowiol does not dissolve completely. We recommend centrifugation for 15 minutes at 5000 g. Continue with the supernatant.
5. **Optional:** Subsequently, add 25 mg DABCO per ml solution and dissolve until a clear solution is obtained.
6. Pour the medium into 1.5 ml- or 2 ml-tubes and store at -20 °C.

#### Application:

Thaw a Mowiol- or Mowiol/DABCO-aliquot and keep it on ice until use.

**Note: Mowiol/DABCO-aliquots should undergo 2 cycles of freeze-and-thaw in maximum.**

Apply the tissue specimen onto the slide and remove the excessive buffer. Don't allow the sample to fall dry.

Pipet 50-100 µl Mowiol or Mowiol/DABCO without bubbles onto a cover slip – preferably forming a „T“. Apply the cover slip bubble-free to the slide by lowering with caution (s. Fig.1). Let the samples harden over night.

**Completely hardened specimens can be stored at room temperature.**


Fig. 1

Lower the cover slip wet side down onto the slide starting at the cross of the „T“, in order to disperse the embedding medium optimally.

|               |        |       |
|---------------|--------|-------|
| <b>Mowiol</b> | 0713.1 | 50 g  |
|               | 0713.2 | 100 g |

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
 P.O. Box 100121 • 76231 Karlsruhe  
 Phone: +49 (0) 721/ 5606-0  
 Fax: +49 (0) 721/ 5606-149  
 info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.