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## ROTI®Load DNASTAIN

### SYBR® Green

ready-to-use, 6x konz.

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### SYBR® Green

ready to use, 6x conc.

DNA gel loading buffer with fluorescent dye  
SYBR® Green I

#### I. Introduction

Versatile reagent for gel loading and simultaneous staining of double stranded nucleic acid in agarose and polyacrylamide gels. ROTI®Load DNASTAIN SYBR® Green detects up to 0.01 ng per band, and is, therefore, 20times more sensitive than ethidium bromide.

#### II. Composition

ROTI®Load DNASTAIN SYBR® Green is composed as follows:

TRIS/HCl pH 7,5	60 mM
EDTA	60 mM
Glycerol	60 % (w/v)
DNASTAIN 1: Bromphenol blue, Xylene cyanol FF, SYBR® Green I	
DNASTAIN 2: Xylene cyanol FF, Tartrazine, SYBR® Green I	
DNASTAIN 3: Orange G, SYBR® Green I	

\*SYBR® is a registered trademark of Molecular Probes, Inc., now Invitrogen Corporation, Carlsbad, California, USA

The glycerine content in ROTI®Load DNASTAIN SYBR® Green increases the density of your DNA solution, ensuring that it distributes evenly at the bottom of the gel pocket. The tracking dyes bromophenol blue, xylene cyanol FF, tartrazine, orange G are used to facilitate sample application and as run markers. The addition of Tris and EDTA increases the stability of the DNA. SYBR® Green I is a fluorescent dye that stains the nucleic acid during the gel run.

#### III. Mechanism

The fluorescent dye SYBR® Green I has been developed specifically for DNA analysis. DNA integrity is not influenced by SYBR® Green I; hence, gel eluted DNA may directly be applied to standard down-stream applications, like, e.g., ligations, PCR, transformation, transfection etc. without interfering effects.


SYBR® Green I binds to the minor groove of the DNA and intercalates into the DNA<sup>1,2</sup> but is less mutagenic than ethidium bromide<sup>3</sup>. The specific binding keeps the gel matrix uncoloured and the background signals low.

SYBR® Green I can be excited by UV- and blue light, making it fully compatible with standard UV-transilluminators and modern blue-light illumination. It is also ideally suited for real-time electrophoresis (see also ROTIPHORESE®-PROfessional runVIEW, Art. No. 4849.1).

When excited, nucleic acid-bound SYBR® Green I emits a brightly coloured fluorescence that may be documented by standard ethidium bromide- or SYBR® Green photo filters.

Excitation maxima (bound to DNA, blue curve):  
ca. 254 nm (UV-light) and 495 nm (blue light)

Emission maximum (bound to DNA, green curve):  
521 nm

Suitable for staining of dsDNA and dsRNA of the following fragment lengths:

DNASTAIN 1: > 500bp

DNASTAIN 2: 100 to 2000 bp

DNASTAIN 3: < 500 bp

#### IV. Application

ROTI®Load DNASTAIN SYBR® Green can either be added to DNA-solution or the DNA can be dissolved directly in ROTI®Load DNASTAIN SYBR® Green.

##### 1. Dissolving DNA in 1 x ROTI®Load DNASTAIN SYBR® Green

Dilute 1:6 by pipetting 170 µl ROTI®Load DNASTAIN SYBR® Green into a clean, sterile test tube. Add 830 µl bidistilled, autoclaved water and mix by vortexing briefly. This solution (ROTI®Load DNASTAIN SYBR® Green 1x) can now be used to dissolve precipitated DNA and can be applied directly onto the gel.

## 2. Adding ROTI®Load DNASTAIN SYBR® Green to a DNA-solution

Add one part of ROTI®Load DNASTAIN SYBR® Green to 5 parts of the DNA sample as shown in the table below. Spin down the tubes and vortex or pipet gently a few times to achieve homogeneity.

Volume DNA-solution		Volume loading buffer	=	Total volume
5 µl	+	1 µl	=	6 µl
10 µl	+	2 µl	=	12 µl
15 µl	+	3 µl	=	18 µl
20 µl	+	4 µl	=	24 µl
25 µl	+	5 µl	=	30 µl

## OPTIONAL

In order to avoid run artefacts, in some cases (e.g. *EcoRI/Hind III*-digested Lambda-DNA) denaturing the sample prior to the run may be recommended. Incubate for 5 mins. at 65 °C and then cool on ice.

## V. Storage

For long-term storage, ROTI®Load DNASTAIN SYBR® Green should be stored at -20 °C (> 3 months). Buffer currently in use can be stored for some weeks at +4 °C. Keep protected from light.

## VI. Safety

SYBR® Green DNA dye is classified as non-hazardous and not harmful to the environment. Nevertheless, it is good laboratory practice to always wear gloves when working with the product and not to dispose of the dye solution down the sink (s. security data sheet).

## VII. Quality Assurance

All components of the buffer are proved to be DNase-/RNase-free. The production is carried out under sterile conditions.

## VIII. Frequently Asked Questions

Q: Are there any samples of ROTI®Load DNASTAIN SYBR® Green?

A: No, sorry. By tube-internal evaporation, the quality of samples of less than 50 µl decreases rapidly. You may get some interesting test discount, however.

Q: How should I visualise the gels after staining with ROTI®Load DNASTAIN SYBR® Green?

A: Gels can be visualised under UV-light or on a blue-light transilluminator.

Q: On UV-table, the ROTI®Load DNASTAIN SYBR® Green stained bands are very bright, but on photos, there is virtually no band visible. What is wrong?

A: Some ethidium bromide photo filters are very narrow banded, efficiently excluding other wavelengths as the green range needed here. Use broadband photo-filters or filters for SYBR® Green.

Q: Can ROTI®Load DNASTAIN SYBR® Green be used for post-run staining?

A: No, this application is not possible.

Q: What is the shelf life of ROTI®Load DNASTAIN SYBR® Green?

A: It can be kept for 1 year minimum under cooling.

Q: How should I dispose of ROTI®Load DNASTAIN SYBR® Green?

A: You should follow any existing departmental or company guidance on the disposal of reagents as indicated in the MSDS of ROTI®Load DNASTAIN SYBR® Green. It should not be disposed of down the sink.

Q: Does ROTI®Load DNASTAIN SYBR® Green intercalate between base pair stacks?

A: Yes, the contained dye SYBR® Green I intercalates between base pair stacks<sup>1,2</sup>.

Q: Is ROTI®Load DNASTAIN SYBR® Green mutagenic?

A: The contained dye SYBR® Green I has the potential to be mutagenic when present in high concentrations. However, ROTI®Load DNASTAIN SYBR® Green contains only a very small amount of SYBR® Green I and is therefore harmless and safe.

Q: Is ROTI®Load DNASTAIN SYBR® Green membrane permeable?

A: The contained dye SYBR® Green I can penetrate the cell membrane<sup>4,5</sup>. In the presence of DMSO, this property is particularly pronounced, as DMSO increases the permeability of tissue<sup>6</sup>. However, ROTI®Load DNASTAIN SYBR® Green is an aqueous solution without DMSO, so there is no amplification effect from DMSO.

Q: May I also use ROTI®Load DNASTAIN SYBR® Green for staining of PAGE gels?

A: Yes, you can use it also for loading of DNA samples into the wells of PAGE gels.

Q: Does ROTI®Load DNASTAIN SYBR® Green also efficiently stain supercoiled plasmids?

A: Yes, the staining efficiency for supercoiled plasmids is as good as that of ethidium bromide.

Q: What to do if the bands on the gel are too faint?

A: Increase the added volume of ROTI®Load DNASTAIN SYBR® Green up to 2times.

Q: Does ROTI®Load DNASTAIN SYBR® Green also efficiently stain RNA?

A: You need double strands of nucleic acid for binding of the dye. RNA with long hairpins or other double strands is stained; RNA with long single stranded structure may be stained less efficiently.

Q: May gels stained with ROTI®Load DNastain SYBR® Green be used for Southern blotting and detection?

A: Yes, without reservation. During blotting and denaturation, the dye simply diffuses from the DNA.

For research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or in vitro diagnostics.

Literature:

1. Nucleic Acids Research, Zipper *et al.* 32, 12 (2004).
2. J. Fluorescence, Dragan *et al.* 22, 1189 (2012).
3. Mutat Res, Singer *et al.*, 439, 37 (1999)
4. Acta Histochemica, Briggs & Jones, 107, 301 (2005)
5. J. Micro. Methods, Kirchhoff & Cypionka, 142, 79 (2017)
6. Molecular Probes, Useful Tips (2003)

**ROTI®Load DNastain 1 SYBR® Green**

1.8 ml	1CN5.1
5 x 1.8 ml	1CN5.2

**ROTI®Load DNastain 2 SYBR® Green**

1.8 ml	1CN6.1
5 x 1.8 ml	1CN6.2

**ROTI®Load DNastain 3 SYBR® Green**

1.8 ml	1CN7.1
5 x 1.8 ml	1CN7.2

**Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5 • 76185 Karlsruhe

P.O. Box 100121 • 76231 Karlsruhe

Phone: +49 (0) 721/ 5606-0

Fax: +49 (0) 721/ 5606-149

info@carlroth.com • www.carlroth.com ed 06.2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.