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50 bp-DNA-Leiter SYBR® Green

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50 bp DNA Ladder SYBR® Green

ready-to-use, prestained with SYBR® Green I

I. Introduction

The DNA ladder is used for size determination of DNA fragments pre-stained with SYBR® * Green I, e.g. with the gel loading buffer ROTI®Load DNastain 3 (Art. No. 1CN7). It can be excited with UV and blue light and is ideally suited for real-time electrophoresis (see ROTIPHORESE®-PROfessional runVIEW, Art. No. 4849).


Figure:
5 μl of 50 bp DNA Ladder SYBR® Green on 3.0 % agarose

*SYBR® is a registered trademark of Molecular Probes, Inc., now Invitrogen Corporation, Carlsbad, California, USA

II. Composition

Regular DNA ladder consisting of 10 fragments with 50 to 500 bp, prestained with SYBR®Green I. The 250 bp fragment has been enhanced for easy analysis of the gel.

The marker was generated by PCR amplification based on Lambda phage as template DNA. Following PCR, the DNA was purified, spectroscopically measured and used for marker mixture in the amounts given below. The marker is solubilised in 1x concentrated sample buffer (Tris HCl, pH 7.5, 10 mM; EDTA 10 mM; glycerol 10 %; SYBR® Green I; orange G) and may be applied directly to the gel. The concentration is 105 ng DNA/μl marker solution.

III. Mechanism

The fluorescent dye SYBR® Green I binds to the minor groove of dsDNA and intercalates into DNA^{1,2}, but is less mutagenic than ethidium bromide³ and 20x as sensitive. The specific binding keeps the gel matrix uncoloured and the background signals low. DNA integrity is not influenced by SYBR® Green I; hence, gel eluted DNA may directly be applied to standard down-stream applications.

IV. Application

The DNA ladder is optimised for direct loading on unstained agarose gels. The standard loading of MINI or MIDI gels is 5 μl (500 ng) per lane.

The dye can be excited by UV- and blue light. We recommend use of broadband ethidium bromide or SYBR® Green photo filters.

Excitation maxima (bound to DNA): ca. 254 nm (UV-light) and 495 nm (blue light)

Emission maximum (bound to DNA): 521 nm

The 50 bp DNA ladder SYBR® Green is highly compatible with the 100 bp DNA ladder SYBR® Green (Art. No. 1CN9). The result is a versatile ladder for the small to medium fragment range from 50 bp to 1000 bp.

Fragment sizes** (in bp)	DNA mass (ng/μl marker solution)
500	10
450	10
400	10
350	10
300	10
250	15
200	10
150	10
100	10
50	10

**Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.

V. Storage

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. If necessary aliquote the marker in portions. Short term storage (few weeks) may be done at +4 °C. *Store in the dark!*

VI. Safety

The 50 bp ladder SYBR® Green is classified as non-hazardous and not harmful to the environment. Nevertheless, it is good laboratory practice to always wear gloves when working with the product.

VII. Quality Assurance

All components of the buffer are proved to be DNase-/RNase-free. The production is carried out under sterile conditions.

VIII. Frequently Asked Questions

Q: How can I visualise the bands of the 50 bp DNA ladder SYBR® Green after gel run?

A: Gels can be visualised under UV-light or on a blue-light transilluminator.

Q: On UV-table, the bands are very bright, but on photos, there is virtually no band visible. What is wrong?

A: Some ethidium bromide photo filters are very narrow banded, efficiently excluding other wavelengths as the green range needed here. Use broadband photo-filters or filters for SYBR® Green.

Q: What is the shelf life of the 50 bp DNA ladder SYBR® Green?

A: It can be kept for 1 year minimum under cooling.

Q: How should I dispose of the 50 bp DNA ladder SYBR® Green?

A: You should follow any existing departmental or company guidance on the disposal of reagents as indicated in the MSDS of the 50 bp DNA ladder SYBR® Green. It should not be disposed of down the sink.

Questions about the fluorescent dye SYBR® Green I contained in the ladder:

Q: Does SYBR® Green I intercalate between base pair stacks?

A: Yes, the dye intercalates between base pair stacks^{1,2}.

Q: Is SYBR® Green I mutagenic?

A: SYBR® Green I has the potential to be mutagenic when present in high concentrations. However, the 50 bp DNA ladder SYBR® Green contains only a very small amount of SYBR® Green I and is therefore harmless and safe.

Q: Is SYBR® Green I membrane permeable?

A: SYBR® Green I can penetrate the cell membrane^{4,5}. In the presence of DMSO, this property is particularly pronounced, as DMSO increases the permeability of tissue⁶. However, the 50 bp DNA ladder SYBR® Green is dissolved in an aqueous loading buffer without addition of DMSO, so there is no amplification effect from DMSO.

For research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or in vitro diagnostics.

Literature:

1. Nucleic Acids Research, Zipper *et al.* 32, 12 (2004).
2. J. Fluorescence, Dragan *et al.* 22, 1189 (2012).
3. Mutat Res, Singer *et al.*, 439, 37 (1999)
4. Acta Histochemica, Briggs & Jones, 107, 301 (2005)
5. J. Micro. Methods, Kirchhoff & Cypionka, 142, 79 (2017)
6. Molecular Probes, Useful Tips (2003)

**50 bp DNA Ladder SYBR® Green
1CN8.1**

500 µl

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com

ed 06/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.