

Instructions for use



LuciPac® A3

From Kikkoman

For analysis of contamination of surfaces via ATP/ADP/AMP measuring by Lumitester Smart.

1. Introduction

Using the LuciPac® A3 sampling pen, surfaces can be tested quickly and easily for organic contamination. The organic contamination of surfaces, e.g. by bacteria, is indicated by the presence of ATP / ADP / AMP in the samples. An enzyme/substrate complex in the LuciPac® A3 sampling tube reacts with ATP / ADP / AMP to produce bioluminescence whose intensity correlates directly with the amount of ATP / ADP / AMP present. The Lumitester PD30/Smart then determines the bioluminescence as a direct measure of original contamination.

The unit is easy to operate, even by non-professionals. Just probe the surface with the swab stick. Then put the swab back into the tube and push it down into the reaction chamber. Shake the tube until the reagents are released and measure the bioluminescence in the Lumitester PD30/Smart.

2. Shipment and storage

20 LuciPac® A3 each are packed in a resealable aluminium bag. Kits are to be stored at a low temperature (2 °C to 8 °C (35.6 °F to 46.4 °F)) for long term storage. The kit can be stored below 25 °C (77 °F) for up to 14 days or below 30 °C (86 °F) for up to 5 days before opening an aluminum bag without any adverse effect on the long term stability. Do not freeze the device.

We recommend that you use all 20 sampling devices in a single bag at one time after opening an aluminum bag. If you have leftover sampling devices that you must store after opening a bag, be sure to store them at the recommended low temperature (2 °C to 8 °C (35.6 °F to 46.4 °F)) and use them within two weeks after opening.

Expiry date: printed on the label of the aluminum bag. Store protected from light.


3. Comments on application, detection limits

Do not use this device for any purpose other than what it was designed for as a tester of cleanliness levels of surfaces. Please be aware that this kit cannot be used to test or identify specific types of pathogenic bacteria. Surfaces may carry ATP / ADP / AMP which are caused by improper cleaning of facilities or originate from microorganisms propagating in such environment. Since this device measures ATP / ADP / AMP levels immediately, indications of improper cleaning processes can be promptly detected.

This device shall not provide any guaranty that a given test sample is free of bacterial contamination.


4. Mechanism

This test system uses a cyclic method based on a combination of firefly luciferase, pyruvate orthophosphate dikinase (PPDK) and a pyruvate kinase (PK). The method produces a defined amount of luminescence proportional to the amount of adenosine triphosphate (ATP), adenosine diphosphate (ADP) and adenosine monophosphate (AMP) in a sample.



Firefly luciferase emits light in the presence of ATP and luciferin. The AMP generated from this reaction is converted back to ATP using the PPDK, while existing ADP is converted to ATP by the pyruvate kinase (PK) to produce a luminescent signal that is as good and stable as possible. ATP is a source of energy necessary for various forms of life that is present in organic residues, such as microorganisms, food residues, and biological substances that originate from other living organisms. This method allows you to measure and detect organic residues and micro-organisms at high speed and high sensitivity by detecting ATP using luciferase, to monitor cleaning processes and microbiological control of coolant or other industrially processed water. Furthermore, you can use this technique to determine not only ATP, but also ADP and AMP. AMP is formed from ATP in bacteria when they are in a resting state and shut down the energy balance. ADP and AMP are derived from ATP during processing, such as heat treatment or fermentation.

5. Composition of product



Reagent cocktail:

Luciferin, luciferase, phosphoenolpyruvate (PEP), magnesium acetate, pyrophosphorylase pyruvate/phosphate dikinase

Releasing reagent:

Alkyl dimethylbenzylammoniumchloride

<u>Structural parts</u>	<u>Raw material</u>
Stick holder (green)	Polypropylene
Test tube (body)	Polypropylene, aluminium
Container for releasing reagent	
Measurement tube	

6. Please note

Please make sure to follow the instructions outlined below in order to obtain optimal performance from this device.

1. Do not use products with a shelf life that has already expired. Expired products may not yield accurate results (the expiry date is printed on the label of the aluminum bag holding the sampling devices).
2. Be sure to use only designated products, i.e. type of Lumitester when performing luminescent measurements. This device cannot be used with unqualified products.
3. Residual disinfectants may interfere with measurements. In case there is residual disinfectant on the surface to be probed, clean this surface using water prior to swabbing.
4. Prior to testing the sampling devices should be allowed about 20 min to reach room temperature (20 °C to 35 °C (68 °F to 95 °F)), if they are from a refrigerated stock. Measurement values may deviate from the real value, if the swabs are used without equilibration. Use the sampling devices as soon as possible once they are available at room temperature. Do not keep the device at temperatures exceeding 35 °C (95 °F). High temperatures may cause product performance to drop.
5. Wherever possible, be sure to use up all sampling devices from a single bag that has been opened at one time. However, if you must store leftover sampling devices once finished with a test session, firmly close the aluminum bag and store it in a refrigerated environment (2 °C to 8 °C (35.6 °F to 46.4 °F)). High temperatures may cause product performance to drop.
6. Do not subject the device or any part of it to direct sunlight for long periods of time. Strong light may cause product performance to drop.
7. Do not touch any of the parts inside the device, particularly not any part of the sampling stick (swab) itself with a finger or other object before use. Touching the parts may affect cleanliness levels, making them hard to discriminate.
8. Do not drop the device or any of its parts or allow any parts to be struck or jolted. The interior aluminum sheets and other parts in the kit may become damaged, causing product performance to drop.
9. Do not use the device if any parts become damaged, e.g. the inner aluminum sheets separating the reagent chamber from the swab as well as the luminescent reagent. Such damage will affect product performance. The aluminum sheet is not damaged if the releasing reagent stays in the chamber upon agitation.
10. Ensure that the LuciPac® A3 is not shifted or displaced, discard the item, because accuracy cannot be guaranteed, swabs are not removable from the Lumitester PD30/Smart, or the instrument may cause an error (malfunction).
11. If the test reagent is leaking please don't use the swab!

7. Measurement

Complete the procedures listed below within a temperature range of 20 °C to 35 °C (68 °F to 95 °F). Make sure to always run measurement tests at the same temperature to maintain repeatability for comparison. Remove the LuciPac® A3 from refrigerator and wait until they have reached room temperature (about 20 min, 20 °C to 35 °C (68 °F to 95 °F)). Use the sampling devices then as soon as possible.

Residual disinfectants may interfere with measurements. In case there is residual disinfectant on the surface to be probed, clean this surface using water prior to swabbing.

1. Remove the swab from the main body (casing). Moisten the tip of the swab with tap water
2. Swab the respective surface repeatedly with the swab tip. For an optimal statement, an area of approx. 10 x 10 cm² is recommended.
3. Return the sampling stick to the main body (casing) and push it all the way into the main body. In order to open the reaction chamber, push firmly onto LuciPac® A3 casing.
4. Shake the whole LuciPac® A3 repeatedly so that all remaining liquid in the capsule falls into the reaction tube. Allow the leftover luminescent reagent to thoroughly dissolve. Slightly shake the test tube in order to promote this process.
5. Insert the LuciPac® A3 into the Lumitester to measure the results.
6. Remove the LuciPac® A3 case from the measuring chamber.


8. Safety recommendations

Please observe the following items to ensure safe use of this product.


- Please note safety recommendations given in the Safety Data Sheet!
- Do not attempt to drink the kit reagent or touch it with bare hands or allow it to splash into eyes. Please make sure to read the precautions and instructions in this Instruction Manual before attempting to use the kit and exercise extreme caution when using it.
- Be careful not to get fingers caught when inserting the sampling stick into the main body (casing).
- Handle with care when storing and disposing of the device and its reagents to ensure that none of the substances become mixed with food and other products.
- Please make sure to store this kit and its parts out of the reach of young children.
- Take precaution when disposing of this kit after use to ensure that such substances do not become mixed with food products at food production centers and similar facilities

9. Short instruction

Prior to use, make sure the LuciPac® A3 sticks have come to room temperature. But then use immediately.




1



Pull the cotton swab stick (A) out of the main body (B).

2




When the test object is dry, moisten the swab (A) with tap water*⁴⁾ or wet the object with water*⁴⁾

3



Swab the test object with the swab stick (A).

4




Put the swab stick (A) back to the main body (B) and push it through all the way by putting the tip of the test tube (C) on a palm of hand or table. (Be careful not to get fingers caught when pushing it.)

5



Shake the whole body of the LuciPac a few times so that the liquid in the capsule falls into the test tube (C).

6



Gently shake the whole body of the LuciPac so that the luminescent reagent is entirely dissolved.

7



Insert the whole body of the LuciPac into the measurement chamber of Lumitester. And close the chamber cover.

8



Press the "ENTER" key. Results are obtained in 10 seconds.

LuciPac® A3

20 pieces

1ENH.1

100 pieces

1ENH.2

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021