



# Instructions for use

## ROTI® Prep Gel & PCR

20H6

**Kit for DNA isolation from agarose gel pieces, or from PCR and sequencing reactions.**

### 1. Introduction and product description

- Great flexibility through dual use
- Preparation in the well-established mini-column system
- Fast, easy and reliable
- Extraction time of approx. 3 min for PCR and Sequencing reactions, 20 min for agarose gel pieces
- Fragment length: 60 bp – 30 kbp
- Rate of recovery: 60-95 %
- Average Ratio  $A_{260}:A_{280}$ : 1.7-2.0

ROTI® Prep Gel & PCR was designed to provide the greatest flexibility in the purification of isolated DNA fragments. The kit is used for the purification of DNA from PCR reactions as well as for the purification from gel pieces and requires a minimum of washing steps.

The purification from PCR and sequencing reactions is based on a new two-step process in which the DNA is first bound to the column membrane using the Binding Buffer and then eluted in the second step.

In the purification from agarose gel pieces, the gel piece is first dissolved and the DNA is then washed in several washing and centrifugation steps before it can be eluted.

#### **Suitable source material:**

- TAE and TBE agarose gels (up to 300 mg)
- PCR reactions (up to 50 µl)
- Sequencing reactions (up to 50 µl)"

### 2. Product use and warranty

#### 2.1 Intended Use

- For research use only.
- Please read the instructions for use carefully and follow the instructions accordingly. Reliability of results cannot be guaranteed if there are any deviations in the procedure.

Information in this document is subject to change without notice. Carl Roth GmbH + Co. KG assumes no responsibility for any errors that may appear in this document.


Carl Roth GmbH + Co. KG disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Carl Roth GmbH + Co. KG be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

**Literature Citation:** When describing a procedure for publication using this product, please refer to it as *Carl Roth's ROTI®Prep Gel & PCR*

## 2.2 Safety Precautions

All due care and attention should be exercised in handling the materials and reagents contained in the kit. Always wear gloves while handling these reagents and avoid any skin contact! In case of contact, flush eyes or skin with a large amount of water immediately.

### Caution:

Gel Solubilizer  Danger H314-H302-H312-H332-H412

Binding Enhancer  Danger H226-H314-H315-H319

Binding Buffer  Danger H225-H319-H336-H315

**Attention:** Do not add bleach or acidic components to the waste after sample preparation!

**MSDS:** the appropriate MSDS can be downloaded from our website [www.carlroth.com](http://www.carlroth.com).

## 3. Materials provided in this kit and storage conditions

### 3.1. Included Kit components

| Amount               | Component                    | Storage |
|----------------------|------------------------------|---------|
| 8 / 40 / 180 ml      | Gel Solubilizer              | RT      |
| 1 / 5 / 15 ml        | Binding Enhancer             | RT      |
| 6 / 30 / 140 ml      | Binding Buffer               | RT      |
| 3 / 2x 8 / 2 x 40 ml | Washing Solution WSL (conc.) | RT      |
| 2 / 2x 2 / 15 ml     | Elution Buffer EB            | RT      |
| 10 / 50 / 5x 50      | Mini spin columns            | RT      |
| 10 / 50 / 5x 50      | 1.5 ml Elution tubes         | RT      |
| 10 / 50 / 5x 50      | 2 ml Collection tubes        | RT      |

The ROTI®Prep Gel & PCR Kit should be stored dry, at room temperature (14-25 °C). Before every use make sure that all components have room temperature.

If there are any precipitates within the provided solutions solve these precipitates by careful warming. This kit is stable for at least 6 months after receipt, when stored as directed.

Contents of this Kit may not be bought separately.

The user is responsible to validate the performance of the ROTI®Prep Gel & PCR Kit for any particular use, since the performance characteristics of our kits have not been validated for any specific application.

### 3.2. Required Material and Equipment not included in this kit

- Ethanol 96-99.8 %, e.g. Ethanol p.a. (e.g. 9065.1)

## 4. Application

ROTI®Prep Gel & PCR was developed to provide the greatest flexibility in the purification of isolated DNA fragments. Spin column-based preparation allows elution in a small volume of low-salt buffer, eliminating time-consuming phenol-chloroform extraction or alcohol precipitation. The columns may be used either in micro-centrifuges or on vacuum manifolds.

The kit allows the extraction of high purified DNA from PCR or sequencing reactions, or, alternatively from agarose gel pieces .

Store the extracted DNA at +4 °C. For long time storage placing at -20 °C is recommended.

For centrifugation we recommend a standard microcentrifuge. Centrifugation steps should be carried out at room temperature.

### 4.1. Before start, be sure to...

- **Add Ethanol** (96-99.8 %, not incl. in this kit) to the Washing Solution WSL (conc.) as follows:  
20H6.1 (10 Preps): 12 ml (15 ml final vol.)  
20H6.2 (50 Preps): 32 ml to each bottle (40 ml final vol.)  
20H6.3 (250 Preps): 160 ml to each bottle (200 ml final vol.)  
Mix thoroughly and keep the bottle always firmly closed!
- Heat one thermal mixer or water bath to 56 or 60 °C.
- Pre-fill the needed amount of Elution Buffer EB into a 2.0 ml reaction tube and incubate at 60°C until the elution step.

### 4.2. Workflow

The workflow has been provided on a separate page, in order to use a copy for in-process protocol or filing in your personal lab files.

- DNA Extraction from agarose gel slices (TAE or TBE agarose gels) p. 3
- Purification and Concentration of PCR Products p. 4

#### 4.2.1 DNA Extraction from agarose gel slices (TAE or TBE agarose gels)

| Step (RT = room temperature)  | done |
|---|------|
| <b>1. Sample Preparation</b>  |      |
| Excise the DNA fragment from the agarose gel with a sharp scalpel and transfer it into a 1.5 or 2 ml reaction tube (max. 300 mg weight) |      |
| Add 650 µl <b>Gel Solubilizer</b> and incubate the slice for approx. 10 min at 50 °C until the gel is completely dissolved              |      |
| Add 50 µl <b>Binding Enhancer</b> to the sample. Mix by vortexing or pipetting up and down  |      |
| <b>2. Column Loading</b>  |      |
| Place each Spin Column into a 2 ml collection tube  |      |
| Apply the mix of sample/ Binding Enhancer to the spin column.   |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 1 min/RT and discard the flow-through  |      |
| <b>3. Column Washing</b>  |      |
| Place the Spin Column back into the 2 ml collection tube  |      |
| Add 700 µl of <b>Washing Solution WSL</b> to the Spin Column  |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 1 min/RT and discard the flow-through  |      |
| Place the Spin Column back into the 2 ml collection tube  |      |
| Add 700 µl of <b>Washing Solution WSL</b> to the Spin Column  |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 1 min/RT and discard the flow-through and the collection tube                                |      |
| Place the Spin Column back into the 2 ml collection tube  |      |
| Again centrifuge columns at 12.000-14.000 rpm (or full speed) for 3 mins/RT in order to remove residual ethanol.                        |      |
| <b>4. Elution</b>   |      |
| Place the Spin Column into a clean 1.5 ml elution tube  |      |
| Add 30-50 µl (or 10-20 µl for higher concentration) <b>Elution Buffer EB</b> to the centre of the membrane                              |      |
| Incubate for 2 min at room temperature  |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 1 min/RT to elute DNA  |      |
| Repeat the elution step (by loading the eluted DNA again on the spin column) to increase the yield of extracted DNA                     |      |

#### 4.2.1 Purification and Concentration of PCR Products

| Step (RT = room temperature)   | done |
|--|------|
| <b>1. Binding of PCR fragments</b>   |      |
| Place each Spin Column into a 2 ml collection tube   |      |
| Add 500 $\mu$ l <b>Binding Buffer</b> to the spin filter   |      |
| Add up to 50 $\mu$ l PCR reaction to the spin filter. Mix by vortexing or pipetting up and down                                    |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 3 min/RT and discard the collection tube and the flow-through                           |      |
| <b>4. Elution</b>  |      |
| Place the Spin Column into a clean 1.5 ml elution tube   |      |
| Add 30-50 $\mu$ l (10-20 $\mu$ l for higher concentration of PCR fragments) <b>Elution Buffer EB</b> to the centre of the membrane |      |
| Incubate for 2 min at room temperature   |      |
| Centrifuge at 11.000 g (or 13.000 rpm) for 1 min/RT to elute DNA   |      |
| Repeat the elution step (by loading the eluted DNA again on the spin column) to increase the yield of extracted DNA                |      |

**Ordering information:**  
(for detailed kit content see Table under 1.)



|                                |                  |               |
|--------------------------------|------------------|---------------|
| <b>ROTI®Prep Gel &amp; PCR</b> | <b>10 preps</b>  | <b>20H6.1</b> |
| <b>ROTI®Prep Gel &amp; PCR</b> | <b>50 preps</b>  | <b>20H6.2</b> |
| <b>ROTI®Prep Gel &amp; PCR</b> | <b>250 preps</b> | <b>20H6.3</b> |

### **Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com

LHa 03/2023

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.