## 20 bp-DNA-Leiter (ready-to-use)

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# Instructions for use

## 20 bp DNA Ladder (ready-to-use)


This marker was generated by PCR amplification based on Lambda phage as template DNA. Following PCR, the DNA was purified, spectroscopically measured and used for marker mixture in the amounts given below.

### Delivery includes:

20 bp DNA Ladder, 50 µg in 500 µl  
[100 ng/µl]

Fragment sizes <sup>1</sup> (in bp)	DNA-mass (ng/µl marker solution)
200	10
180	10
160	10
140	10
120	10
<b>100</b>	<b>15</b>
80	10
60	10
40	10
20	10

<sup>1</sup> Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.

### Figure:

1.4 % agarose (High Resolution, Art. No. K297) + 0.8 % Synergel™ (Art. No. 0184), in 1x TAE.

left: 7 µl 20 bp DNA Ladder,  
right: 5 µl each of 20 bp DNA Ladder and 100 bp DNA Ladder, *equimolar*.

### Application:

Our 20 bp DNA Ladder is solubilised in 1x concentrated sample buffer (Tris HCl, pH 7.5, 10 mM; EDTA 10 mM; glycerol 10 %; orange G 0.4 mg/ml) and may be applied directly to the gel. The concentration is 1 µg DNA / 10 µl marker solution.

This gel loading buffer resembles a 6x diluted solution of gel loading buffer ROTI®Load DNA orange 1 (Art.-No. HP04.1) with EDTA concentration enhanced to 10 mM, in order to provide optimal stabilisation of the PCR fragments.

### Sample application:

The standard loading of MINI to MIDI gels per lane is:

- with bands visible in UV-light after Et.-Br.-staining: 0.5 – 1.0 µg
- with detection after Et.-Br.-staining with signal-integrated camera systems: 0.1– 0.5 µg

### Storage:

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. If necessary aliquote the marker in portions. Short term storage (few weeks) may be done at 4 °C.

### Please note:

The 20 bp DNA Ladder is highly compatible with the 100 bp DNA Ladder, *equimolar* (T834, see figure, right lane). Mixed in equal shares (f. e. 5 µl 20 bp DNA Ladder and 5 µl 100 bp DNA Ladder, *equimolar*) they result in a versatile range of DNA marker bands for analysis of all small DNA fragments.

## 20 bp DNA Ladder

2805.1

500 µl

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com

sse 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100420, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.