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### ROTI®Load 3 (LDS)

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## Instructions for use

### ROTI®Load 3 (LDS)

Protein loading buffer, non reducing  
(4x conc.)

- Loading buffer based on LDS
- Optimised for ready-made gels
- Superior solubility at low temperatures

#### Introduction:

ROTI®Load 3 (LDS) is a modified, non-reducing protein sample-gel loading buffer for gels under denaturing or native conditions. ROTI®Load 3 (LDS) is based on a triethanolamine buffered ready-to-use solution, in which the standard sodium dodecylsulphate (SDS) has been replaced by lithium dodecylsulphate (LDS). At low temperatures, LDS is much better soluble than SDS. Hence, no crystallisation of LDS takes place, and protein stability as well as band resolution is significantly enhanced. Additionally, the formulation of ROTI®Load 3 (LDS) has been optimised for ready-made gels, particularly ROTI®-PAGE gels, but also ready-made gels of other standard suppliers.

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

ROTI®Load 3 (LDS) prevents degradation of proteins during sample preparation and electrophoresis run. The triethanolamine buffer system keeps the sample pH value constant over a wide temperature range, additionally stabilising the sensitive aspartyl-prolyl peptide bonds.

#### Use:

1. Prepare your polyacrylamide gel.
2. Determine sample volume.
3. Add ROTI®Load 3 (one third of sample volume).

#### Example:

60 µl sample + 20 µl ROTI®Load 3 (4x conc.)

4. Mix by pipetting.
5. Heat samples for 3-5 mins. to 80-90 °C.
6. Directly load samples to the gel.

#### Formulation:

LDS (approx. 4 % w/v), glycerol (approx. 40 % v/v), Ficoll® 400 (approx. 4 % w/v), phenol red (approx. 0,025 % w/v), Coomassie® G250 (approx. 0,025 % w/v).

Triethanolamine buffered.

#### Store at +20 °C

#### Tip:

For reducing conditions, ROTI®Load 3 (LDS) may be supplemented with 20 µl/ml DTT solution (1 M) or with β-mercaptoethanol.

#### Safety Information:

 **Warning** H315-H319

P280-P302+P352-P305+P351+P338-P337+P313

**ROTI®Load 3 (LDS)** 3359.1 10 ml

#### Loading buffers based on SDS:

ROTI®Load 1 (reducing, Art. No: K929)

ROTI®Load 2 (non reducing, Art. No: K930)