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## Masson-Goldner-Trichrom-Färbekit

Best.-Nr. 3459

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# Instructions for use

## Masson Goldner Trichrome Staining Kit

Art. No. 3459

### Trichrome staining for visualisation of connective tissue

Trichrome staining works with three different stains allowing a differentiated visualisation of tissues.

At first the nuclei are stained with Weigert's iron hematoxylin solution. A simple hematoxylin staining does not suffice for the nuclei would be decolourised by the acid Goldner's stains. Subsequently, the trichrome staining is carried out. The dyes differ in particle size: *Goldner's stain I* contains a fine-particle phase (Ponceau), which infiltrates quickly all structures of tissue, and a coarse-particle phase (Acid Fuchsin), which works more slowly. At first it stains only the coarse structures of tissue by masking the fine-particle phase. You have to stop the staining procedure at that point to avoid an overstaining of the tissue.

During the differentiation step with *Goldner's stain II* (Phosphotungstic acid) it is important to decolour the connective tissue as much as possible. Only then it can be stained with *Goldner's stain III* (Light green).

Goldner's stain I: Muscle fibers and cytoplasm red (Ponceau), connective tissue red (Acid Fuchsin)

Goldner's stain II: Erythrocytes orange (Orange G), connective tissue decolourised (Phosphotungstic acid)

Goldner's stain III: Connective tissue green (Light green SF yellowish)

#### Kit includes:

- Hematoxylin solution A acc. to Weigert (X906.1)
  - Danger** H225-H319-H336 500 ml
- Hematoxylin solution B acc. to Weigert (X907.1)
  - Danger** H290-H318 500 ml
- Goldner's stain I (Ponceau - Fuchsin) (3469.1) 500 ml
- Goldner's stain II (Phosphotungstic acid - Orange G) (3470.1)
  - Danger** H315-H318 500 ml
- Goldner's stain III (Light green SF yellowish) (3473.1) 500 ml

The staining solutions should be filtrated before use! Solutions may be bought separately.

#### Additional chemicals required:

- Acetic acid solution 1% (Acetic acid 100%, Art. No. 3738)
- Clearing Agents: ROTI®Histol (Art. No. 6640) or ROTICLEAR® (Art. No. A538) or Xylene p.a. (Art. No. 4436)
- Appropriate Mounting Media: ROTI®Histokitt (Art. No. 6638), compatible with ROTI®Histol  
 ROTI®Mount (Art. No. HP68), compatible with ROTICLEAR®  
 ROTI®Histokitt II (Art. No. T160), compatible with Xylene

**Instruction:** Use dewaxed, rehydrated sections.

1. Stain with Iron hematoxylin solution acc. to Weigert (Mix solution A+B at a ratio of 1:1). <b>max. 3 min</b>	6. Rinse with Acetic acid solution 1%. <b>30 sec</b>
2. Blue in flowing tap water. <b>10-15 min</b>	7. Counterstain with Goldner's stain III. <b>2-5 min</b> <i>Optionally: Examine by microscope.</i>
3. Stain with Goldner's stain I. <b>5-10 min</b>	8. Wash with Acetic acid solution 1%. <b>2-5 min</b>
4. Rinse with Acetic acid solution 1%. <b>30 sec</b>	9. Dehydrate by ascending alcohol series. 10. Clear with clearing agent. 11. Mount with appropriate mounting medium.
5. Stain with Goldner's stain II until decolouration of connective tissue. <b>1-3 min (normally sufficient) up to 30 min</b> <i>Examine by microscope. Specimens may not dry out.</i>	<b>Please note:</b> Instead of Light green solution it is possible to use a 0.1-0.2% solution of Anilinblue (Art. No. 4002) for counterstaining (step 7).

**Result:** Nuclei: dark brown; Cytoplasm, muscle fibers: red; Erythrocytes: orange; Connective tissue: green

**Please note:** The colour intensity depends on the pre-treatment and the composition of the samples to be stained. It may initially be necessary to adapt the method to the respective conditions.

**Masson-Goldner's trichrome staining kit**

**3459.1**

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021

Revisionen: 10.2021