RNase Inhibitor

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Instructions for use

RNase Inhibitor

40 U/μl, for biochemistry and molecular biology

This RNase inhibitor is a 50 kDa recombinant human placental protein expressed in *Escherichia coli*. The RNase inhibitor inhibits the ribonuclease (RNase) activity of common eukaryotic enzymes such as RNases A, B and C. This inhibition takes place via a non-covalent bond in the ratio of 1:1.

RNase inhibitors are mainly used in applications where the presence of RNases has a negative influence on the quality of RNA or the results of experiments. These include RNA isolation, cDNA synthesis, RT-PCR, *in vitro* transcription and translation or RNase-free monoclonal antibody preparations.

This RNase inhibitor shows no activity towards RNase 1, RNase T1, RNase T2, S1 nuclease and RNase H.

Storage buffer: 20 mM HEPES-KOH (pH 7.6), 50 mM KCl, 8 mM reducing agent, 50 % (v/v) glycerol

Unit Definition: One unit is defined as the amount of enzyme required to inhibit the activity of 5 ng RNase A by 50%.

Features:

- inhibits RNase A, B and C activity
- no polymerase or reverse transcriptase activity
- free of DNase and RNase activity
- Active over a broad pH range (pH 5.5 – 9.0)

Application:

- RNA isolation and purification
- cDNA synthesis, RT-PCR, RT-qPCR
- *in vitro* transcription and translation
- RNase-free monoclonal antibody preparation

Protocol:

- The optimal final concentration of RNase inhibitor in a reaction depends on the level of RNase contamination, the incubation time and the compounds present in the reaction mixture. Usually this is in the range of 1 - 2 U/μl.
- For a standard reverse transcription reaction, use 1 μl (40 U) of RNase inhibitor in a sample volume of 20 μl.
- For optimal activity, a final DTT concentration of 0.5 - 1 mM is essential.
- The RNase inhibitor should be added prior to any components that may be possible sources of RNase contamination.

- The use of the RNase inhibitor does not preclude treatment with RNase H after amplification of the first strand cDNA.

Further information:

- The storage buffer contains 8 mM DTT, but if the ratio of inhibitor to final sample volume is less than 1:8, the addition of DTT to a final concentration of 0.5 - 1 mM DTT is recommended.
- The purity of >90%, was determined using SDS-polyacrylamide gels.

RNase Inhibitor	2000 U	plastic	3727.1
	10000 U	plastic	3727.2

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.