ZelluTrans/Roth Mini Dialyzer (MD)

Vol. 10-100 µl

MWCO	1 MD in microtube	12 MD in microtube	1 strip (= 8 MD)
3500 (3.5 kDa)	4765.1	4766.1	4767.1
6000-8000 (6-8 kDa)	4769.1	4770.1	4771.1
12000-14000 (12-14 kDa)	4772.1	4775.1	4776.1
Schwimmhilfe/Float			4768.1
Abdeckfolie/Sealing Film			4781.1

Technical Data Sheet and Applications

Content

Technical Data Sheet.....	- 2 -
Manual (English Version).....	- 3 -
Manual (German Version).....	- 4 -
Chemical Resistance	- 5 -
Sample Volume Recovery.....	- 6 -
Application Note Protein Binding (BSA).....	- 7 -
Application Note Urea Removal.....	- 8 -
Application Note Salt Removal.....	- 9 -
Application Note Dialysis Speed of different MWCO	- 10 -
Application Note Dialysis of Sucrose	- 11 -
Application Note CHAPS Removal	- 12 -
Application Note SDS Removal	- 13 -
Application Note Volume reduction by evaporation	- 14 -

Technical Data Sheet

ZelluTrans/Roth Mini Dialyzer

Specifications:

- max. sample volume 100 µl
- temperature 1-60 °C
- for aqueous solution
- incubation time up to 24 h
- pH 3-10
- cut off 3.5, 6-8, 12-14 kDa
- storage 4-22 °C
- strips with 8 samples, scalable from 1 to 96 samples
- conform to Microplate Standard (SBS)
- membrane: low binding regenerated cellulose, contain glycerole to prevent embrittlement and traces of elements like sulphides and heavy metals


Applications:

- Protein and peptide sample purification e.g. desalting before mass spectrometry
- Optimization of protein renaturation with different renaturation buffers and steps
- Removal of dyes after protein labeling
- Protein sample rebuffering
- Glycoprotein modification and engineering
- Protein in vitro translation
- Cell culture: studies of a cell line with a virus strain
- Enzyme activity assays
- Plasmid or primer purification

Manual – English Version

ZelluTrans/Roth Mini Dialyzer

Preparation

For a smaller number of samples, break carefully the desired number of single dialyzers from the dialyzer strip or use a single dialyzer in microcentrifuge tube. Please do not touch the membrane!

Buffer Preparation


Pipette dialysis buffer into microcentrifuge tubes (max 1.4 ml/tube) or into a 96-deep well plate (max 1.8 ml/well).

Loading Sample

Load pipette (2-200µl) with sample (max 100µl). Insert the pipette tip upright into the round opening and press on slightly. Inject sample slowly.

Introduction

Place the single dialyzer into a prepared microcentrifuge tube (upright position) and the dialyzer strip into a prepared 96-deep well plate, respectively.

Dialysis

Dialysis can be done in one step. To accelerate the procedure change buffer solution repeatedly. For this purpose simply put the single dialyzer into a tube with fresh dialysis buffer or change the position of the dialyzer strip in the deep well plate.

Dialysis Time

The dialysis time depends on the compounds of the sample and the cut-off of the membrane. Typical dialysis times range from 30 min to 4 h, with repeated replacement of the dialysis buffer (every 15-60 min). For dialysis of salts change buffer every 30 min.

Sample Retrieval

Use pipette 2-200µl and set volume to 140µl. Press the piston to the first stop, hold it, put the pipette upright into the round opening. Aspirate sample carefully.


Storage

Finally, pipette the sample into a microcentrifuge tube or a microtest plate for further analytics.

Materials Required

Single dialyzer or dialyzer strip (consisting of 8 connected single dialyzers).

Deep well plate (2 ml) or microcentrifuge tubes (2 ml)

Pipette 2-200 µl (for sample), pipette 1-2 ml (für dialysis buffer).

Microcentrifuge tubes or microtest plates for sample storage


Forceps to use for disposal of single dialyzers out of microcentrifuge tube.

Suitable protection gloves and goggles

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Test Sample: 100 µl of Congo Red color solution in water dd

Dialysis solution: 1,4 ml of tested chemical

Incubation: 18 h

Determination Method: Optical integrity and leak-tightness to air pressure

Acetonitrile	Good	Acetic acid 25%	Good
Acetone	Good	Acetic acid 96%	Good
Chloroform	Good	Formic acid 25%	Good
Dimethyl sulfoxide	Good	Formic acid 100%	No
Ethanol 70%	Good	Hydrochloric acid 10%	Limited
Ethanol 98%	Good	Hydrochloric acid 25%	No
Ethylacetate	Good	Hydrochloric acid 37%	No
Ethylene glycol	Good	Hydrofluoric acid 50%	No
Glycerol	Good	Nitric acid 25%	No
n-Hexane	Good	Nitric acid 65%	No
<i>iso</i> -Propanol	Good	Phosphoric acid 25%	Limited
Methanol 98%	Good	Phosphoric acid 85%	No
Methylene chloride	Good	Sulfuric acid 98%	No
1-Propanol	Good		
Tetrahydrofuran	Good	Ammonium hydroxide 1N	Limited
Toluene	Good	Ammonium hydroxide 25%	Limited
		Potassium hydroxide 1N	Limited
Hydrogen peroxide 30%	Good	Potassium hydroxide 32%	No
		Sodium hydroxide 1N	Limited
		Sodium hydroxide 32%	No
		Good chemical resistance	Good
		Limited chemical resistance, e.g. pore size cannot be guaranteed	Limited
		No chemical resistance, use not recommended	No

Sample Volume Recovery

ZelluTrans/Roth Mini Dialyzer


Materials and Methods

ZelluTrans/Roth Mini Dialyzer: MWCO 6-8 kDa (Art. No. 4769.1 ff)

Sample: 10 µl, 50 µl, or 100 µl respectively of Aq. dest.

Determination Method: Gravimetry

Sample	Recovery [%]
MD dry 10 µl	84.86
MD dry 50 µl	88.01
MD dry 100 µl	93.60
MD wet 10 µl	100.00
MD wet 50 µl	94.53
MD wet 100 µl	94.37
MD: Mini Dialyzer	
Dry: dry membrane	
Wet: prewetted membrane	

Application Note - Protein Binding (BSA)

ZelluTrans/Roth Mini Dialyzer

Materials and Methods


Sample: 100µl of BSA solution in PBS (20 and 100 µg/ml), n= 5 - 6

Dialysis Buffer: 1,8 ml of Aq. dest., exchange after 1 h, Incubation time: 4h at RT

Determination Method: Protein determination according to Bradford,

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
Remark: Due to the high mobility of urea in solutions and the short diffusion distances in ZelluTrans/Roth Mini Dialyzer, there is only a minor osmotic effect, resulting in a slightly higher volume recovery of about 110 µl.

Application Note - Salt Removal

ZelluTrans/Roth Mini Dialyzer

Materials and Methods


ZelluTrans/Roth Mini Dialyzer: MWCO 6-8 kDa (Art. No. 4769.1 ff)

Sample: 100µl of 50mM H₂NaPO₄/HNa₂PO₄, 150mM NaCl, pH 7.4

Dialysis Buffer: 1.8 ml Aq. dest., Buffer exchange interval 15 min, 20 °C

Determination Method: Wescor Osmometer Vapor 5520

t in min	Osmolarity in mOsm/kg
0	506
15	154
30	77
45	58
60	42
75	35
90	35
105	28
120	27
aqua	15

Application Note - Dialysis Speed of different MWCO

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
Sample: 100µl of 1mM *p*-Nitrophenol in PBS, pH 7.4, n=3

Dialysis buffer: 1.8 ml of PBS, pH7.4, exchange interval 30 min.

Determination Method: Tecan Sunrise Photometer, 420 nm

Time / min	mM <i>p</i> -NP		
	3.5 kD	6-8 kD	12-14 kD
0	1.06	1.06	1.06
15	0.66	0.61	0.64
30	0.51	0.43	0.47
60	0.30	0.21	0.26
120	0.10	0.07	0.05

Application Note - Dialysis of Sucrose

ZelluTrans/Roth Mini Dialyzer

Materials and Methods


ZelluTrans/Roth Mini Dialyzer: MWCO 6-8 kDa (Art. No. 4769.1 ff)

Sample: 100µl of 1M, 0,5 M or 0,1 M Sucrose, respectively, in Aq. dest

Dialysis Buffer: 1.8 ml Aq. dest., Buffer exchange interval 30 min, 20 °C

Determination Method: Wescor Osmometer Vapor 5520

Time min	Residual Sucrose in % of original concentration		
	0.1 M	0.5M	1M
0	100.0	100.0	100.0
30	83.2	55.5	42.1
60	53.1	35.3	20.0
120	25.6	23.2	16.1
240	7.8	11.8	

Remark: With higher molarities of sucrose (e.g. 0,5 M and 1 M), an increasing osmotic effect appears resulting in dilution of the sample, which is simulating a higher speed of dialysis. To prevent this, it is best to keep the concentration difference around 100 mM. If this is not possible or desired, reduce the sample volume to 50 µl and close the ZelluTrans/Roth Mini Dialyzer with sealing for MD. For salt, this osmotic effect was not as pronounced as for sugars.

Application Note - CHAPS Removal

ZelluTrans/Roth Mini Dialyzer

Materials and Methods


ZelluTrans/Roth Mini Dialyzer: MWCO 6-8 kDa (Art. No. 4769.1 ff)

Sample: 100µl of 1 mg/ml BSA in Aq. dest., 1% (w/v) CHAPS

Dialysis Buffer: 1.8 ml Aq. dest., Buffer exchange interval 2 h for 8 h, 20 °C

Determination Method: Agilent HP1100/MSD

t in h	Residual CHAPS in %
0	100
2	53,4
4	19,8
8	2,5
24	0,0

Application Note - SDS Removal

ZelluTrans/Roth Mini Dialyzer

Materials and Methods


ZelluTrans/Roth Mini Dialyzer: MWCO 6-8 kDa (Art. No. 4769.1 ff)

Sample: 100µl of 1 mg/ml BSA and 0.1 % (w/v) SDS in PBS, pH 7.4

Dialysis Buffer: 1.8 ml Aq. dest., Buffer exchange interval 2 h for 8 h, 20 °C

Determination Method: Agilent HP1100/MSD, ESI negative mode

t in h	Residual SDS in %
0	100
2	63,8
4	41,2
8	25,3
24	12,6

Application Note - Volume reduction by evaporation

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
Sample: 100 µl of PBS (Phosphate buffered saline: 10mM H₂NaPO₄/ HNa₂PO₄, 150mM NaCl, pH 7.4)

Determination Method: Gravimetric analysis of the weight of the residual solution in the Mini Dialyzer

Experimental: MD were filled with 100 µl PBS and placed in an "Egg Crate" (from Corning) which was attached to a stand to allow air to pass the Mini Dialyzer. Ambient temperature, n=2.

Time / min	Weight of residual solution in Mini Dialyzer/ mg		
	Sample 1	Sample 2	Mean
0	100	100	100
15	81	78	79.5
30	64	62	63
45	55	52	53.5
60	46	47	46.5
75	38	45	41.5
90	30	22	26
120	26	11	18.5

Remark: Results may vary with temperature and air flow