Art.-No. 6478.1 125 ml
Art.-No. 6478.2 3 x 125 ml

OCT Embedding Matrix

Instructions for Use


www.cellpath.com | T: +44 (0)1686 611333

OCT Embedding Matrix

KMA-0100-00A

- 3 **EN** Instructions for Use
- 4 **DA** Brugsanvisning
- 5 **DE** Gebrauchsanweisung
- 6 **EL** Οδηγίες χρήσης
- 7 **ES** Instrucciones de uso
- 8 **FI** Käyttöohjeet
- 9 **FR** Instructions d'utilisation
- 10 **HU** Használati utasítás
- 11 **IT** Istruzioni per l'uso
- 12 **LV** Lietošanas instrukcija
- 13 **NL** Gebruiksaanwijzing
- 14 **NO** Bruksanvisning
- 15 **PL** Instrukcja użytkowania
- 16 **PT** Instruções de utilização
- 17 **RO** Instrucțiuni de utilizare
- 18 **SV** Användar instruktioner

EC	REP	EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands
----	-----	---

CellPath 
INNOVATION IN CELLULAR PATHOLOGY

 **CellPath Ltd**
80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com

OCT Embedding Matrix

KMA-0100-00A

EN

A tissue freezing and embedding material intended to be used as an embedding matrix for the preparation and processing of frozen sections from biological tissues or clinical specimens.

Ready to use solution.

For laboratory use.

Keep in a cool, dry, well ventilated area. Keep containers tightly closed. Store in correctly labelled containers.

Please see expiry on product label.

Instructions for Use

Important: Tissue stored at -80°C should be placed in the chamber for at least 30 minutes prior to sectioning, along with any instruments which are required for use when sectioning, such as: brushes, forceps, etc. Leave microscope slides to be used to collect tissue section at room temperature.

- 0 Ensure the cryostat to be used is cold, with a working temperature set between -10°C to -20°C. Place the tissue specimen to be sectioned in the cryostat chamber for 20 minutes, in order for it to reach the same temperature as the chamber.
- 6 Remove cap and cut the nozzle to enable dispensing.
- @ Place a specimen disk inside the chamber and apply CellPath OCT Embedding Matrix to the specimen disk. When the OCT matrix begins to freeze (starts to turn a white colour), place the tissue specimen in the OCT matrix and freeze completely, in order to mount the specimen on the disk. Additional mounting medium may be applied around the tissue for extra support to the specimen.
- 8 Replace cap.
- 0 Once the specimen is securely embedded in frozen OCT matrix, place the specimen disk in the cryostat specimen head, orient your tissue as desired and adjust the plane of the specimen.
- 0 Begin trimming the specimen at section thickness of 5-150 µm. Once you are ready to collect a tissue section, set the cryostat to the desired section thickness ~5µm and place the anti-roll plate in position. It is recommended to

cut several practice sections, so that the user can make any adjustments required to prepare smooth tissue sections.

Warning: For user safety, please ensure the microtome/knife blade is always guarded when not in use and the cryostat section wheel is locked.

Note: If sections curl or shred, adjust the anti-roll plate position in order to achieve a high-quality tissue section. Only minor adjustments of the anti-roll plate are necessary to make a difference in tissue section quality.

& Once all the necessary adjustments are made, lift the anti-roll plate away and use a fine section brush to carefully position the section prior to collecting the section on a room temperature positively charged microscope slide.

Important: Once the OCT section is mounted on a slide it needs to be fixed as quickly as possible by immersion in 95% ethanol for 45 to 60 seconds. It is very important that the section is fixed as quickly as possible after it is mounted on the slide, this is paramount for a high-quality stained specimen. Then wash slides gently in running tap water then rinse in distilled water to remove the OCT.

@ Once all desired sections have been collected, remove the remaining tissue specimen from the specimen head. The residual remaining tissue can be easily removed from the specimen disk, by removing the disk from the chamber to allow the OCT matrix to thaw.

0 Use a large section brush to sweep up discarded tissue sections and clean the cryostat chamber with 100% Ethanol.

Important: DO NOT use water as it will cause the cryostat chamber to become frosty.

(ID) Dispose of empty bottle in accordance with local or government guidelines.

! Any complaints or adverse incidents need to be reported to CellPath immediately.

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CellPath Ltd, 80 Mochdre Enterprise Park, Newtown, Powys, SY16 4LE, UK