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200 bp-DNA-Leiter (ready-to-use)

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Instructions for use

200 bp DNA Ladder (ready-to-use)

This marker was generated by PCR amplification based on Lambda phage as template DNA. Following PCR, the DNA was purified, spectroscopically measured and used for marker mixture in the amounts given below.

Delivery includes:

200 bp DNA Ladder, 50 µg in 500 µl
[100 ng/µl]

Fragment sizes ¹ (in bp)	DNA mass (ng/µl marker solution)
2000	10
1800	10
1600	10
1400	10
1200	10
1000	15
800	10
600	10
400	10
200	10

¹Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.


Figure:

5 µl 200 bp DNA Ladder on 1.8 % agarose.

Application:

Our 200 bp DNA Ladder is solubilised in 1x concentrated sample buffer (Tris HCl, pH 7.5, 10 mM; EDTA 10 mM; glycerol 10 %; bromophenol blue, xylene cyanole) and may be applied directly to the gel. The concentration is 1 µg DNA / 10 µl marker solution.

This gel loading buffer resembles a 6x diluted solution of gel loading buffer ROTI®Load DNA with glycerol (Art. No. X904.1) with EDTA concentration enhanced to 10 mM, in order to provide optimal stabilisation of the PCR fragments.

Sample application:

The standard loading of MINI to MIDI gels per lane is:

- With bands visible in UV-light after ethidium bromide staining: 0.5-1.0 µg (5-10 µl).
- With detection after ethidium bromide staining with signal-integrated camera systems: 0.1- 0.5 µg (1-5 µl).

Storage:

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. If necessary aliquote the marker in portions. Short term storage (few weeks) may be done at +4 °C.

Please note:

The 200 bp DNA Ladder is highly compatible with the 20 bp DNA Ladder (Art. No. 2805). The two markers can easily be combined, resulting in a highly versatile ladder from 20 bp up to 2000 bp. We recommend use of this marker on the following special agaroses: HR-PLUS (Art. No. HP30), LM/PCR (Art. No. HP31), or an agarose with additional Synergel™ (Art. No. 0184).

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Carl Roth GmbH + Co. KG