

# Instructions for use



## ROTI®Pol TaqUltra

DNA-free recombinant Taq DNA polymerase for PCR amplifications, particularly of bacterial DNA  
9347

### 1. Description

Recombinant full-length form of the heat stable Taq DNA polymerase from the thermophilic bacterium *Thermus aquaticus* in storage buffer, and tested on the absence of bacterial DNA.

For research use only. Not approved for use in clinical or *in vitro* diagnostics.

### 2. Applications

This ROTI®Pol TaqUltra is recommended for use in PCR applications in general, or if bacterial DNA, or in RT-PCR, bacterial 16S rRNA shall be detected. The polymerase is appropriate for use in the amplification of DNA from genomic, viral, and plasmid templates, for high-yield PCR, for colony PCR and TA-cloning.

ROTI®Pol TaqUltra is purified using a multiple-step process that minimises contaminating bacterial DNA to a none detectable level. Each lot of the polymerase undergoes strict quality control testing in order to ensure the absence of detectable amounts of contaminating bacterial DNA.

ROTI®Pol TaqUltra is able to amplify PCR products up to 3 kb and is appropriate for use in the amplification of DNA from eukaryotic as well as prokaryotic templates. The TaqUltra DNA polymerase possesses a 5' → 3' polymerase- as well as a 5'-3' exonuclease activity, and generates a 3'dA (adenine)-overhang which may well be used for TA-cloning purposes.

### 3. Content

ROTI®Pol TaqUltra polymerase (Art. No. 9347) in storage buffer containing 50 % glycerol  
Filled in orange-capped tubes.

| Reagent             | Lid colour | 9347.1 | 9347.2  |
|---------------------|------------|--------|---------|
| TaqUltra polymerase | orange     | 1 tube | 5 tubes |

This DNA polymerase is provided without reaction buffer.

Since the ROTI®Pol TaqUltra is a particularly high purified version of the TaqS polymerase, the buffer included in the ROTI®Pol TaqS products (as well as in the ROTI®Pol TaqHY products) is perfectly suited for use as reaction buffer for the TaqUltra as well (see also 12. Additionally recommended products). In addition to our specifically optimised reaction buffers, general Tris-based Taq polymerase reaction buffers as, for instance, given in *Molecular Cloning* (Sambrook and Russell, Cold Spring Harbour Laboratory Press), work well for the ROTI®Pol TaqUltra enzymes. We recommend use of 20 mM MgCl<sub>2</sub> and 0.1 % Tween®-20.

**Please note**, however, that the batches of our reaction buffers *have not been tested for presence of DNA or residual bacteria* which might lead to false-positive results. For each buffer that shall be used for the TaqUltra enzyme, either taken from a ROTI®Pol product or self-made, we recommend to do preliminary studies in order to ensure the freeness of DNA functioning as target for the particular primer pairs to be used.

#### 4. Storage Buffer

50 mM Tris-HCl (pH 8.0), 100 mM KCl, 0.1 mM EDTA, 0.5 % IGEPAL CA-630, 0.5 % Tween-20, 1 mM DTT, 50 % glycerol

#### 5. Enzyme activity

5 units/ $\mu$ l enzyme solution

#### 6. Unit definition

One unit of activity is defined as the amount of enzyme required to incorporate 10 nmoles of dNTP into an acid-insoluble DNA fraction in 30 minutes at 72 °C.

#### 7. Suggested pipetting scheme

Prepare on ice:

| Components                                  | Apply for PCR reaction of 20 $\mu$ l volume | Final concentration (recommended) |
|---|---|-----------------------------------|
| PCR buffer (10x)*                           | 2 $\mu$ l                                   | 1x                                |
| dNTP-Mix (2 mM)                             | 2 $\mu$ l                                   | 800 $\mu$ M (200 $\mu$ M each)    |
| Forward primer (e.g. 5 pmol/ $\mu$ l)       | variable (e.g. 1 $\mu$ l)                   | 0.1-0.5 $\mu$ M                   |
| Reverse primer (e.g. 5 pmol/ $\mu$ l)       | variable (e.g. 1 $\mu$ l)                   | 0.1-0.5 $\mu$ M                   |
| Template DNA                                | variable                                    | 0.01-10 ng / reaction             |
| ROTI®Pol TaqUltra polymerase (5 U/ $\mu$ l) | variable (i.e. 0.2 $\mu$ l)                 | 0.5-1.5 U                         |
| Sterile dest. water                         | adjust to 20 $\mu$ l final volume           |                                   |

\*also see 3. Content

#### 8. Basic amplification protocol

| Step                 | Time                                       | Temperature |
|----------------------|--|-------------|
| Initial denaturation | 2 minutes                                  | 92-95 °C    |
| <b>25-40 cycles</b>  |  |             |
| Denaturation         | 2-10 seconds                               | 92-95 °C    |
| Annealing            | 2-10 seconds                               | 55-68 °C    |
| Extension            | variable, depends on the length of product | 72 °C       |

#### 9. Notes

**IMPORTANT:** For specific species-DNA detection and to avoid false positive results make sure your PCR buffer and other PCR reagents are free of DNA contaminations.

For maximum yield and specificity, annealing temperatures and annealing time as well as extension time and cycle numbers should be optimised for each template target and primer pair. Usually the optimal annealing temperature is 2-5 °C below the melting temperature of the primers. Recommended elongation time is 30-60 secs. per 1 kb of target. Elongation times of 30 secs. per 1 kb may be sufficient but longer elongation times may be necessary depending on the complexity of the template DNA.

Product is not covered by pending or issued patents or may have certain limitations. To our best knowledge, however, this product does not provide any conflict with pending or issued patents.

## 10. Recommended MgCl<sub>2</sub> concentration

2-4 mM

In case the MgCl<sub>2</sub> concentration has to be adjusted, use a separate MgCl<sub>2</sub> solution (100 mM) in PCR quality and add in appropriate amounts according to the scheme below. We recommend doing PCR with a MgCl<sub>2</sub> gradient in order to find the optimal concentration.

### Pipetting scheme for additional MgCl<sub>2</sub>

| Final MgCl <sub>2</sub> conc. in mM   | 2.5    | 3      | 3.5    | 4      |
|---|--------|--------|--------|--------|
| Add 100 mM MgCl <sub>2</sub> solution in following amounts to 20 µl reaction volume | 0.1 µl | 0.2 µl | 0.3 µl | 0.4 µl |

## 11. Storage conditions

Store at -20 °C. Infrequent short term storage (few hours) of the enzyme may be done at +4 °C.

## 12. Additionally recommended products

For our Thermal cyclers please contact us under 0721 / 5606 - 0

|  |               |
|--|---------------|
| <b>ROTI®Mix PCR 3</b> (10 mM per dNTP dATP, dTTP, dGTP, dCTP)                          | Art. No. L785 |
| <b>ROTI®Mix PCR 3 (pH 7)</b> (10 mM per dNTP dATP, dTTP, dGTP, dCTP)                   | Art. No. 0179 |
| <b>dNTP-Set 1</b> (≥99 %, 100 mM pure solutions dATP, dTTP, dGTP, dCTP)                | Art. No. K039 |
| <b>dNTP-Set 1 (pH 7)</b> (≥99 %, 100 mM pure solutions dATP, dTTP, dGTP, dCTP)         | Art. No. 0178 |
| <b>PCR water</b> for molecular biology, sterile, ready-to-use                          | Art. No. 1HPE |
| <b>Magnesium chloride solution</b> 25 mM, for PCR, for molecular biology               | Art. No. 1HY7 |
| <b>Mineral oil</b> (for or overlaying PCR and other enzymatic reactions)               | Art. No. HP50 |
| <b>ROTI®Nucleic acid-free</b> (ready-to-use solution for removal of DNA from surfaces) | Art. No. HP69 |
| <b>ROTI®Nucleic acid-free eXtra</b> (ready-to-use, gentle solution for DNA removal)    | Art. No. 1312 |
| <b>DNA AWAY®</b> (ready-to-use solution for removal of DNA from surfaces)              | Art. No. X996 |

Please note our full range of DNA polymerases and MasterMixes:

|                                  |               |
|----------------------------------|---------------|
| <b>ROTI®Pol TaqS</b>             | Art. No. 9223 |
| <b>ROTI®Pol TaqS Mix</b>         | Art. No. 9239 |
| <b>ROTI®Pol TaqS Red-Mix</b>     | Art. No. 9241 |
| <b>ROTI®Pol Hot-TaqS</b>         | Art. No. 9245 |
| <b>ROTI®Pol Hot-TaqS Mix</b>     | Art. No. 9248 |
| <b>ROTI®Pol Hot-TaqS Red-Mix</b> | Art. No. 9256 |
| <b>ROTI®Pol TaqHY</b>            | Art. No. 9345 |
| <b>ROTI®Pol TaqHY Mix</b>        | Art. No. 1K33 |
| <b>ROTI®Pol TaqHY Red-Mix</b>    | Art. No. 1K34 |
| <b>ROTI®Pol Hot-TaqHY</b>        | Art. No. 9346 |
| <b>ROTI®Pol ProofRead</b>        | Art. No. 9344 |
| <b>ROTI®Pol TaqUltra</b>         | Art. No. 9347 |
| <b>ROTI®Pol Hot-TaqUltra</b>     | Art. No. 9350 |

|                          |         |               |
|--------------------------|---------|---------------|
| <b>ROTI®Pol TaqUltra</b> | 200 U   | <b>9347.1</b> |
|                          | 1.000 U | <b>9347.2</b> |

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

jh 01/2022

