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## ROTI®Cell Myco-free Lösung

9830

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## ROTI®Cell Myco-free Solution

9830

**For the elimination of mycoplasma contamination in cell culture. Mycoplasma eradication solution, antibiotic mixture.**

### I. Introduction

The contamination of cell cultures with mycoplasmas is a well-known phenomenon in cell biology. Mycoplasmas can produce biochemical and immunological changes in cells and thereby radically alter results. Mycoplasma-infected cells cannot always be discarded, and they multiply rapidly. Thus, it is important to eliminate mycoplasmas from cell cultures. The source of the contamination is often the user, or cell culture reagents.

Mycoplasmas lack a cell wall. Hence, they are resistant to many antibiotics that direct their effects towards the bacterial cell wall. Subsequently, mycoplasma contamination may often not be detected in the cell culture medium, since, while the antibiotics added to the medium prevent growth of other bacteria, mycoplasmas can proliferate unrestricted.

### II Product

ROTI®Cell Myco-Free is a reliable, fast-acting reagent for the elimination of contamination by mycoplasmas and other bacteria in cell culture.

ROTI®Cell Myco-Free is a sterile mixture of specific antibiotics in PBS, which are already effective at very low concentrations. The antibiotic action is based on the specific disruption of bacterial protein biosynthesis by impeding translation and inhibiting the gyrase of the transcription apparatus. Eukaryotic cells are not affected.

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5  
76185 Karlsruhe  
Postfach 100121  
76231 Karlsruhe  
Telefon: +49 (0) 721/ 5606-0  
Telefax: +49 (0) 721/ 5606-149  
E-Mail: [info@carlroth.de](mailto:info@carlroth.de)  
Internet: [www.carlroth.de](http://www.carlroth.de)

gh 01/2020

Advantages:

- Simple application
- Reliable eradication
- Effective on all mycoplasmas and bacteria
- No cytotoxicity with standard application

### III. Concentration and pH

50x stock solution in PBS (0,5 mg/ml).  
pH 4,0-4,5.

### IV. Use

Only for research purposes. Not for diagnostic, therapeutic, or clinical use on humans or animals.

A 100ml bottle is enough for 300 applications.

### V. Storage

Store in the dark (ROTI®Cell Myco-Free is slightly light sensitive).  
Store at -20 °C.  
We recommend you prepare small aliquots for storage at -20 °C.  
Short term storage: 4 °C for 1 week is possible.

### VI. Method

1. Remove the culture medium from the flask and wash with PBS. Treat cells with trypsin or PBS/EDTA.
2. Wash cells thoroughly when passaging.
3. Count cells and add the appropriate cell number, required for culturing the cells, into the cell culture flask.
4. Add 2 ml ROTI®Cell Myco-free per 100 ml cell culture medium (1/50 dilution; working concentration: 10 µg/ml).
5. Cultivate the cells for 2-3 days.
6. If the cell layer or suspension grows too dense, passage cells as usual (wash well) and subsequently add medium plus ROTI®Cell Myco-free to the cells.
7. Test for mycoplasmas using standard methods e.g. a PCR-based method.

If the test shows positive for mycoplasma contamination, repeat steps 1-7 and increase the concentration of ROTI®Cell Myco-free. You may need to increase the number of passages.

The concentration of ROTI®Cell Myco-free can be increased to 1/25 (2x), if required (s. VII. Suggestions)

### VII. Suggestions

The overall time for the initial eradication time in medium plus ROTI®Cell Myco-free should be 7 to 10 days.

If you can still detect Mycoplasma contamination following this incubation time, incubate for a further 7-10 days in medium plus ROTI®Cell Myco-

free at a concentration of (1,25x), 1/30 (1,67x) or 1/25 (2x), depending on the magnitude of residual contamination (2nd eradication time).

You can keep ROTI®Cell Myco-free in the medium for extended time periods, but in that case, please use the lowest concentration of ROTI®Cell Myco-free. Most mycoplasmas will be removed after the first eradication time. Generally, all mycoplasmas will usually have been successfully eliminated after the second eradication time.

Please note:

1. It is important to always test for mycoplasmas after each eradication time, to ensure that they have been completely removed. No eradication solution can guarantee 100% successful removal. In addition, the contamination can potentially reappear after some time, as single cells can be below the mycoplasma test's detection limits.
2. Complete mycoplasma removal should be ensured to prevent the development of antibiotic resistances. Please check for successful mycoplasma eradication using a mycoplasma detection method after each eradication time.
3. Very high concentrations of ROTI®Cell Myco-free can have toxic effects on sensitive cells. Please observe your cells to detect if there is a possible decrease in cell proliferation, or any possible changes in cell morphology. If you can see any toxic effects on the cells after increasing the concentration, use the previously used concentration of ROTI®Cell Myco-free.
4. Cell cultures can be cross-contaminated by other cell cultures. Hence, we recommend regular testing of all cells for mycoplasma contamination. Exchange contaminated media, FBS, trypsin and buffers with fresh cell culture reagents.

**ROTI®Cell Myco-free** 100 ml **9830.1**