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## ROTI®Blue

### Kolloidale Coomassie-Färbung

#### A. Einleitung

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### Colloidal Coomassie staining

#### A. Introduction

ROTI®Blue is a colloidal CBBG-250 staining solution for proteins in polyacrylamide gels. Its colloidal properties allow the dye to accurately bind to the proteins and only minimally to the gel matrix. This means that the extensive washing procedure normally required in background destaining is no longer necessary. The sensitivity of ROTI®Blue solution can be compared to that of silver staining (<30 ng/protein). ROTI®Blue is also suitable for staining of membranes after transfer (blot).

#### B. Staining with ROTI®Blue in general

Proteins can be detected after 30 mins. They are stained on the surface. Depth staining occurs in the course of the next few hours (depending on gel thickness). Background staining remains minimal even after overnight incubation (15 h). Gently shake the gel during the staining process. Depending on shaking and temperature conditions, greater colloid complexes are precipitated during longer staining procedures. The quality of the staining is not impaired.

We recommend that upon completion of staining, the gel be transferred to another clean staining box prior to washing. Precipitated dye complexes on the gel can be dissolved in the washing solution. In case gels are to be blotted after staining we recommend omitting fixation.

#### C. Staining Procedure

The formulas for the fixing solution, ROTI®Blue staining solution, washing solution, stabilizer and drying solution are listed below.

##### A) Staining without fixation

Recommended for gels of up to 12 ml volume prior to blotting

1. Following electrophoresis incubate gel for 2-15 h in 100-200 ml ROTI®Blue staining solution. Proceed as described in B.2. Washing.

Alternatively:

##### B) Staining with fixation

Suitable for all gels that shall not be blotted

1.a. Following electrophoresis incubate gel for 60 mins. in 100-300 ml fixing solution

1.b. Subsequently incubate gel for 2-15 h in 100-300 ml ROTI®Blue staining solution.

2. Washing: Transfer gel to a clean staining box. Incubate gel in 100 ml washing solution whilst shaking for 5 mins. Any possible dye complexes on the gel can be dissolved in the washing solution by jerking and jolting. Gels can now be stored (3.) or dried (4.).

3. Storage: Gels may be stored in stabilizing solution over a longer period. Use an air-tight container as gel shrinks upon loss of water.

4. Drying: Incubate gel whilst shaking in 100 ml drying solution for 30 mins. Change solution to fresh Drying solution after 15 mins. The gel can then be dried between cellophane sheets as described there.

##### C) Staining of membranes prior to Western Blot reactions

1. Following blotting, incubate membrane for 15 mins. to 1 h directly in 100-200 ml 1x ROTI®Blue staining solution. Bands should be stained only weakly.

2. Wash membrane for approx. 5 mins. in washing solution.

3. Prior to detection reactions: Destain membrane for approx. 30 mins. in PBS or TBS.

#### D. Buffers and solutions

Fixing solution (100ml):

79 ml H<sub>2</sub>O + 1 ml o-phosphoric acid 85 % + 20 ml methanol

ROTI®Blue staining solution (100 ml):

**Thoroughly shake the ROTI®Blue (5 x conc.) solution several times before use!**

1. Mix 60 ml H<sub>2</sub>O + 20 ml methanol.
2. Stir in 20 ml ROTI®Blue (5 x concentrate).

**The prepared staining solution is not a clear solution.** It contains dye in colloidal form. The dye beads only dissolve gradually during the dyeing process. Overstaining of the gel is thus avoided as the dye binds immediately to the proteins and only minimally to the gel matrix.

Washing solution (100 ml):

25 ml MeOH 99.8 % + 75 ml H<sub>2</sub>O

Stabilizer solution (100ml):

Dissolve 20 g ammonium sulphate in H<sub>2</sub>O to make 100 ml total volume

Drying solution (100 ml):

10 ml glycerine 99.5 % + 20 ml ethanol + 70 ml H<sub>2</sub>O

#### We recommend use of the following Roth reagents:

- Methanol, p.A., Art.no. 4627.1
- Phosphorsäure, p.A., Art.no. 6366.1
- Ethanol, p.A., Art.no. 9065.1
- Glycerin p.A., Art.no. 3783.1
- Ammoniumsulfat, Art.no. 9212.1
- ROTI@Stock 10x PBS, Art.no. 1058.1
- ROTI@Stock 10x TBS, Art.no. 1060.1
- Gel drying frames MINI Art. No. K420.1
- Replacement cellophane sheets MINI Art. No. K422.1
- Gel drying frames MINI Art. No. K421.1
- Replacement cellophane sheets MINI Art. No. K423.1

#### Hazard and Precautionary Statements

Please note safety data given on label and MSDS.

**Warning** H290 P234-P390-P406

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

<b>ROTI®Blue</b>	<b>A152.2</b>	<b>250 ml</b>
	<b>A152.1</b>	<b>1 l</b>