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ROTI® Quick-Kit

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Literatur

¹ Chomczynski P. and Sacchi N. (1987) Single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction. *Anal. Biochem.* 162:156

ROTI®Quick-Kit

A979.1

Instructions for use

ROTI® Quick-Kit

For Isolation of total RNA from Cells and Tissues

Column-free, very versatile isolation system for reliable isolation of total RNA from nearly every type of tissue.

- Easy, most versatile application
- High recovery rate of intact, very pure total RNA
- Column-free one-step isolation
- Applicable on nearly each tissue type

Suitable for application on:

- Mammalian and non-mammalian tissue
- Cell culture material
- Biopsies
- Fat- or protein rich tissue
- Very hard plant material, and so on.

May be used on fresh as well as on frozen material.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

Kit Contents

ROTI® Quick 1 (A976)

Danger H302+H332-H314-H412-EUH032-EUH208

ROTI® Quick 2 (A977)

Danger H301+H311+H331-H314-H341-H351-H361d-H372-H411

ROTI® Quick 3 (A978)


Danger H225-H319-H336

Contents of this Kit may not be bought separately.

Storage Temperature: Store all solutions at 4 °C.

Method:

This ROTI® Quick-Kit is based on the well-known GITC-method published by Chomczynski und Sacchi¹. Following homogenisation of tissues, RNA is purified in one single step and is then isolated by precipitation.

Time necessary for isolation:

Approx. 2.5 hours; several isolations may easily be performed in parallel.

Up-/Downscaling:

Optimised for RNA isolation from approx. 0.2 g tissue or cell material. In case other amounts shall be used for isolation, please recalculate amounts proportionally.

Isolated RNA:

Isolated total RNA is very pure and may be applied directly to all standard down-stream applications, like e.g. RT-PCR, Northern-Blotting, Reverse Transcription etc.

Protocol

Volumes of solutions given below apply to 0.2 g cell material. In case you have other amounts of cell material please up- or downscale solutions in relation to the amount used.

Please make sure that all materials used are made RNase free (by DEPC-treatment or by heating to 180 °C).

1. **Solid tissue:** Homogenize N₂-frozen tissue in 2.0 ml Quick solution 1
Cell culture cells: Scrape cells from the bottom of culture bottles and suspend cells in 2.0 ml Quick solution 1
Shear cell suspension with a needle using a syringe.
2. Add 2.6 ml 4 °C chilled Quick solution 2 and vortex thoroughly.
Note: Immediately before use, shake Quick solution 2 until mixture is homogenous (2 phases separate when stored).
3. Incubate on ice for 10 mins.
4. Centrifuge at 10.000 rpm for 15 mins. at 4 °C.
5. Pipette the upper aqueous layer into a fresh reaction tube and keep the lower phase for control.

Add 2.0 ml of solution 3 to the upper phase.

Note: In case phases do not separate efficiently (e.g. during extraction from tissue containing fat), steps 2-4 should be repeated using 1 vol. ROTI®Quick 2. Clear phase separation is necessary for pure, isolated RNA.

6. Incubate at -20 °C for 40 mins. or at -80 °C for 10 mins.
7. Centrifuge at 10.000 rpm for 15 mins. at 4 °C.
8. Remove supernatant, resuspend RNA pellet in 600 µl Quick sol. 1 plus 600 µl Quick sol. 3.
9. Incubate at -20 °C for 30 mins. or at -80 °C for 10 mins.
10. Sediment RNA by centrifugation at 4 °C (13.000 rpm, 20 mins.).
11. Add 1 ml 70 % ethanol to the RNA pellet (first wash).
12. Centrifuge at 4 °C (13.000 rpm, 5 mins.).
13. Repeat steps 11 and 12 (second wash).
14. Remove ethanol **completely** from the pellet. Centrifuge shortly (4 °C) and remove all residual ethanol.
15. Dissolve RNA pellet in ca. 200 µl autoclaved DEPC-water.
16. Optional: Quantitation at 260 nm
Standard results are $OD_{260/280} \geq 2,0$.

Additionally Required Reagents

70 % Ethanol (Art. No. T868)
DEPC-treated Water (Art. No. T143)
DEPC (Art. No. K028)

Literature

¹ Chomczynski P. and Sacchi N. (1987) Single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction. *Anal. Biochem.* 162:156)

ROTI®Quick-Kit

A979.1