ROTI®Hybriquick

I. Einleitung

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ROTI®Hybriquick

I. Introduction

ROTI®Hybriquick is a sodium phosphate buffered hybridisation buffer. The buffer is to be used directly as a 1 x hybridisation buffer and is suitable for all Southern- and Northern-Blot hybridisation processes from 50 °C onwards. Depending on specific assay, ROTI®Hybriquick should also be used for the preceding prehybridisation. In many cases hybridisation can be carried out by adding the denatured probe directly to the hot prehybridisation buffer and then returning the hybridisation buffer to the membrane filters. Diluted at 1:2 to 1:10 ROTI®Hybriquick is also suitable, with increasing stringency, for subsequent wash steps. Suitable for all membrane types and detection methods. We recommend using an uncharged nylon membrane (AE50.1) for Northern-Blot hybridisation in order to reduce the background.

II. Technical Information

In order to ensure highest possible flexibility and stability, the buffer contains no additional blocking substances such as salmon or herring sperm DNA. Please add these yourself after denaturing. Standard concentrations are app. 100-500 µg/ml. As the buffer contains SDS, SDS flakes may form at lower temperatures. These are of no consequence with regard to the quality of the buffer and can be dissolved easily by heating.

III. Protocol of a typical Southern-Blot hybridisation with ROTI®Hybriquick

III.I Prehybridisation

- Wet the filter membranes in 50 mM sodium phosphate / 1 mM EDTA (pH 8.0) and equilibrate in warm ROTI®Hybriquick (app. 40-50 °C).

- Prehybridization in app. 10 ml ROTI®Hybriquick per 500 cm² filter surface. We recommend adding 100 µg/ml sheared and denatured herring sperm DNA. Prehybridize for at least 2 hours at 65 °C.

III.II Hybridisation

- Optional: Competition of repetitive DNA fragments by adding 100 µg sheared and denatured genomic DNA of the same species.

- Heat-denaturize probe-DNA by boiling for 5 to 10 minutes and transferring immediately onto ice.

- Cast the prehybridisation buffer from the filters. The probe can be added directly to the prehybridisation buffer or to preheated, fresh ROTI®Hybriquick.

- Immediately add the hybridisation solution to the filters. As a rule, hybridisation should be carried out overnight at 65 °C under gentle shaking or in the rotation oven.

III.III Washing the membranes

- Wash the filter membranes in preheated ROTI®Hybriquick at 65 °C under vigorous shaking:

A) at a dilution of 1:2 for 10 mins

B) at a dilution of 1:5 for 10 mins

C) at a dilution of 1:10 for 10 mins

(Step C can be repeated until the background is reduced).

- Detection of signals can then be carried out via direct X-ray film exposition or chemoluminescence detection.

IV. Protocol of a typical Northern-Blot hybridisation with ROTI®Hybriquick

All solutions should have been freed from RNAses by treating them with 0.1 % DEPC (Art. No. K028.1). Work process should also be RNase-free. We recommend using uncharged nylon membranes.

IV.I Prehybridisation

- Wet the filter membranes in 50 mM sodium phosphate / 1 mM EDTA (pH 8.0) and equilibrate in warm ROTI®Hybriquick (app. 40-50 °C).

- Prehybridisation in app. 10 ml ROTI®Hybriquick per 500 cm² filter surface. We recommend adding 100 µg/ml chopped and denatured herring sperm DNA. Prehybridize for at least 2 hours at 65 °C.

IV. II Hybridisation

- Optional: Competition of repetitive RNAs by adding 100 µg sheared and denatured genomic DNA or denatured RNA of the same species.

- **ds-probes:** heat-denaturize the probe by boiling for 5 to 10 minutes and immediately transferring onto ice. Alternatively: denaturize by adding 0.1 volume 3 M NaOH, incubate for 5 minutes at room temperature. Neutralize by transferring onto icy water and adding 0.05 volume 1 M Tris-Cl (pH 7.2) and 0.1 volume 3 M HCl.

- **ss-probes:** do not have to be denaturized. However, in the event of probes which form a strong secondary structure, denaturation may improve hybridisation.

Following is applicable for subsequent hybridisation: Highly stringent hybridisation (high homology of probe target) at 65-68 °C overnight in 1 x ROTI®Hybriquick. For low-stringent hybridisation (weak homology of probe target), we recommend using following hybridisation buffers:

3 ml ROTI®Hybriquick
+ 5 ml formamide
+ 1 ml SDS (20 %)
+ 1 ml NaCl (2 M)

- Casting the prehybridisation buffer from the filters. The probe can be added directly to the prehybridisation buffer ROTI®Hybriquick or to preheated, fresh formamide hybridisation buffer.

- Add the hybridisation solution to the filters immediately. As a rule, hybridisation should be carried out overnight at 65 -68 °C (highly stringent hybridisation) or 37-42 °C (low-stringent hybridisation) under gentle shaking or in the rotation oven.

IV. III Washing the membranes

- Wash the filter membranes in preheated washing buffer at 65-68 °C (highly stringent hybridisation) / 23 °C (low-stringent hybridisation) under vigorous shaking:

A) 2 x SSC / 0.1 % SDS for 15 mins

B) 1 x SSC / 0.1 % SDS for 15 mins

C) 0.5 x SSC / 0.1 % SDS for 15 mins

Step C can be repeated until the background is reduced. Increasing the SDS-concentration to 1 % can further reduce the background.

- Detection of signals can then be carried out via direct X-ray film exposition or chemoluminescence detection.

ROTI®Hybriquick A981.1 1 |

 Danger H290-H318

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe

P.O. Box 100121 • 76231 Karlsruhe

Phone: +49 (0) 721/ 5606-0

Fax: +49 (0) 721/ 5606-149

info@carlroth.com • www.carlroth.com

ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.