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## ROTI® Liquid Barrier Marker

AN91 (farblos)

AN92 (rot)

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## ROTI® Liquid Barrier Marker

AN91 (colourless)  
AN92 (red)

ROTI® Liquid Barrier Marker consists of hydrophobic substances dissolved in an organic solvent. The solution is applied to slides or other glass surfaces. After drying, a resistant hydrophobic film remains providing a barrier to aqueous liquids.

ROTI® Liquid Barrier Marker is delivered in a glass bottle with screw cap and integrated brush. The solution is colourless or – for better differentiation – stained red.

### Areas of application

ROTI® Liquid Barrier Marker is ideal for immunological, histological and molecular biological tests. Aqueous liquids such as antibody solutions or chromogenic substrates can be enclosed on limited areas of the slide. This allows a more effective working: The required volume of the reagents is reduced; an intermingling of different liquids can be prevented. It is also possible to process 2 different specimens side by side on the same slide.

*We do not recommend using the marker with slides coated with chrome alum gelatin.*

### Application

A thin barrier film is applied to the glass surface, e.g. by forming a circle around the specimen. You can use the integrated brush or other convenient tools such as pipettors. After drying the aqueous liquid can be added. The drying process is completed after ca. 30 seconds depending on the thickness of the film.

It is necessary to close the bottle tightly after use.

### Formulation

ROTI® Liquid Barrier Marker contains different petrol ethers, paraffins and gum arabic. Additionally, AN92.1 is stained with a red dye.

### Chemical compatibility

The barrier film is resistant to ethanol and acetone, but can be wiped off the glass surface with a dry cloth or a cloth impregnated with xylene.

### Storage

Store tightly sealed at room temperature.

**Danger** H225-H304-H315-H319-H336-H411

P210-P280-P301+P310-P302+P352-P304+P340-P403+P235

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com

sse 06/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

### ROTI® Liquid Barrier Marker (colourless)

10 ml	glass bottle with brush	AN91.1
5 x 10 ml	glass bottle with brush	AN91.2

### ROTI® Liquid Barrier Marker (red)

10 ml	glass bottle with brush	AN92.1
5 x 10 ml	glass bottle with brush	AN92.2