Bürkle GmbH | Basler Straße 83 | D - 79540 Lörrach  
 Telefon +49 (0) 7621 9331-0  
 Telefax +49 (0) 7621 9331-31

info@buerkle.de | www.buerkle.de

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Filling Machines | Sampling | Plastic Labware  
 for Laboratory, Industry, Science

Tablette Nr. Pellet No.	Gewicht in g Weight g	Tablette Nr. Pellet No.	Gewicht in g Weight g	Tablette Nr. Pellet No.	Gewicht in g Weight g
1.	40,8	15.	51,9	29.	51,3
2.	47,8	16.	48,3	30.	51,4
3.	49,4	17.	46	31.	52,5
4.	50,8	18.	50,5	32.	50,7
5.	50,3	19.	47	33.	47
6.	46,4	20.	46,3	34.	48,6
7.	49,1	21.	51,8	35.	44,8
8.	44	22.	53	36.	40,1
9.	43,9	23.	52,2	37.	49,6
10.	42,3	24.	49,4	38.	47,6
11.	41,1	25.	51,7	39.	17,1
12.	48,1	26.	50,8		Flasche leer
13.	50,7	27.	52,2		Empty cylinder
14.	47,6	28.	49,5		

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The performance test of SnowPack 50 reached the result of 40 pellets/10 kg bottle which was established under optimal conditions. Practical experiences are showing a result of about 30 pellets (21.6 % below optimum result). The causes which decrease the quantity of CO<sub>2</sub> are different (room-temperature, filling time, etc.).

## Gebrauchsanleitung · Instructions for use SnowPack Trockeneisgerät · Dry Ice Machine Art.-Nr. 9805-0030

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### Description

SnowPack enables you to obtain a -79 °C dry ice tablet weighing 50 g within half a minute, without complicated technology, without electricity! You thus have the possibility for low-cost cooling to very low temperatures at any time you like. Dry ice is odorless and tasteless and easy to handle. The SnowPack dry ice machine can be installed on any commercial CO<sub>2</sub>-bottle with a standpipe (W 21,80x1/14" acc. to DIN 477 No. 6), these are often already present in many production facilities and laboratories. If this is not the case, a bottle can be rented for a small price from the local specialised dealer.

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### Areas of Use

SnowPack is often used in hospitals, universities, schools, institutes, laboratories, in chemical industry, in materials testing facilities. Examples for use: cooling in cooling baths or cooling traps, pre-cooling of tissue samples, the safe transportation or the embedding of biological samples, freezing of blood and tissue samples, analyses in chemistry, physics and materials technology, etc.

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### Safety

Solid CO<sub>2</sub> is intensely cold and can cause cold burns and frostbite to bare skin. Wear protective gloves therefore. Gaseous CO<sub>2</sub> is heavier than air. Care should be taken to disperse the gas formed if the SnowPack is used in a low lying confined space. CO<sub>2</sub> cylinder valves are fitted with a safety disc to guard against overpressure. To avoid the risk of rupture, always store cylinders indoors, in the shade and away from any source of heat. (Cold cylinders make more pellets.)

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
### Montage

Screw SnowPack onto a CO<sub>2</sub> bottle. The bonded fibre sealing washer ensures a gas-tight seal with only light tightening torque of the loose nut. The cylinder must be upright, but the SnowPack handle can be pointed down or horizontal.

With the SnowPack open, turn on the cylinder valve briefly to ensure that the jet is clear and the cylinder is full of liquid CO<sub>2</sub>. Turn the cylinder off and close SnowPack by fastening the clamping belt at the hooks.

### Producing of Dry Ice

Open the cylinder valve fully and inject continuously gas for 30 seconds (for a 50 gramme pellet), then turn off.

When CO<sub>2</sub> injection has ceased, unfasten the clamping belt. The dry ice pellet now lies in the cup attached to the handle. Remove the pellet from the SnowPack, if necessary rocking the handle gently to free the pellet.

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### Problems – What Now?

- No producing of dry ice:
  - you are using a CO<sub>2</sub> bottle without feeding tube. Device only operates thanks to a bottle with feeding tube (please check remarks on the bottle).
- Bad producing of dry ice:
  - bottle is too hot. When the room temperature is too high, SnowPack cannot work correctly.
  - never let the bottle in the sunshine.
- If no dry ice is still being produced and if the mentioned problems have already been checked, the internal feeding tube may be faulty.
  - Please let it replace from your retailer.

Important for you: a CO<sub>2</sub> bottle cannot be totally emptied. About 20 % of the contents always remains in the bottle.

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### Test Specification and Results:

- Fill a 22 lb (10 kg) cylinder and leave it for 12 hours at shop temperature (20 °C).
- Make pellets by injecting for 30 seconds with a 120 second gap between the removal of one pellet to the start of the next.
- Weight each pellet immediately after removal and tabulate weight against sequence number of pellets made.
- Continue until 22 lb cylinder is exhausted. Please note: a CO<sub>2</sub> bottle cannot be totally emptied. About 20 % of the contents always remains in the bottle.

### Results:

Ambient temp. = 20 °C  
 Injection time = 30 secs.  
 Repetition time = 120 secs.