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Kolbenpumpe

CR-KP292

2015-10-30

Carl Roth GmbH + Co. KG

Schoemperlenstr. 3-5
D-76185 Karlsruhe

T: +49 721 5606-0
F: +49 721 5606-149
info@carlroth.de

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Schoemperlenstr. 3-5
D-76185 Karlsruhe

T: +49 721 5606-0
F: +49 721 5606-149
info@carlroth.de

Reprinting or reproduction of this manual, including extracts, is not allowed without the prior written permission.

All rights under the copyright laws are expressly reserved by manufacturer.

We reserve the right to make changes and amendments.

Contents

1	Important Information	4
1.1	General Information	4
1.2	Target Groups	4
1.3	Intended Use.....	4
1.4	Use for an Unauthorized Purpose	4
1.5	Safety Devices	5
1.6	Meaning of the Warning notes	5
1.7	Product Standards, Safety Regulations	5
2	Basic Safety Instructions	6
2.1	General Information	6
2.2	Electricity.....	6
2.3	Mechanical Systems	6
2.4	Hazardous Substances.....	7
2.5	High Temperatures	7
3	Description	8
3.1	Design	8
3.1.1	Connections - Suction-/Pressure side	8
3.1.2	Connecting to the electricity supply	8
3.2	Protection measure against liquids in the pump	9
3.2.1	Condensate separator (Suction side)	9
3.3	Areas of Application	9
4	Technical Data	10
5	Installation and Operation	11
5.1	Unpacking	11
5.2	Setting up and connecting	11
5.2.1	Setting up	11
5.2.2	Connecting.....	11
5.2.2.1	Electrical Connection	11
5.2.2.2	Vacuum connection (Suction side)	11
5.2.2.3	Exhaust connection (Pressure side)	11
5.3	Operation	12
5.3.1	Start-up	12
5.3.2	Decommissioning.....	12
5.4	Storage.....	12
5.5	Scrap Disposal	12
6	Maintenance and Servicing	13
6.1	Maintenance Performed by the User	13
6.2	Maintenance by the Manufacturer	13
7	Troubleshooting	14

EC Declaration of Conformity

Important Information

1 Important Information

1.1 General Information

The **Piston Pump CR-KP292** conforms to the following directives:

2006 / 42 / EC	Machinery Directive
2004 / 108 / EC	Electromagnetic Compatibility Directive
97 / 23 / EC	Pressure Equipment Directive
2014 / 29 EU	Simple Pressure Vessels Directive

The CE sign is located on the rating plate.

Observe the binding national and local regulations when fitting the pump into installations!

1.2 Target Groups

This Operating Manual is intended for the personnel planning, operating and maintaining Piston Pump.

This group of people includes:

- Designers and fitters of vacuum apparatus
- Employees working on commercial laboratory and industrial vacuum technology applications
- Service personnel for Piston Pumps

The personnel operating and maintaining the Piston Pump must have the technical competence required to perform the work that has to be done.

The user must authorize the operating personnel to do the work that has to be done.

The personnel must have read and understood the complete Operating Manual before using the Piston Pump.

The Operating Manual must be kept at the place of use and be available to the personnel when required.

1.3 Intended Use

- The layout of the Piston Pump must be appropriate for the conditions of use. The user bears the sole responsibility for this.
- The Piston Pump may only be operated under the conditions stated
 - in the "Technical Data" section,
 - on the type plate, and
 - in the technical specification for the order concerned.
- The Piston Pump is approved for extracting, pumping and compressing gases and vapours. If these gases and vapours are toxic or explosive, then the user must observe the currently valid safety regulations for this application.

1.4 Use for an Unauthorized Purpose

It is forbidden to use the pump for applications deviating from the technical data stated on the type plate or the conditions stated in the supply contract, or to operate it with missing or defective protective devices.

1.5 Safety Devices

Measures such as the following are for the safety of the operating personnel:

- electrical connection with a protective conductor (operating mode S1) and an earthing plug
- Motor protection device (thermal)
- "Hot Surface" label on the pump - warning notice 

The Piston Pump must not be operated without these elements.

1.6 Meaning of the Warning notes

Take note of the warning notices. They are each in the following box:

	CAUTION ! / WARNING !
Hazard which may lead to serious injuries or material damage.	

1.7 Product Standards, Safety Regulations

The Piston Pump meets the following product

DIN EN ISO 12100-2010	Safety of machinery - General principles for design - Risk assessment and risk reduction
DIN EN ISO 13857:2008-06	Safety of machinery - Safety distances to prevent hazard zones being reached by upper and lower limbs
DIN EN 1012-1 DIN EN 1012-2	Compressors and Vacuum pumps - Safety requirements - Part 1: Compressors Part 2: Vacuum pumps
DIN EN ISO 2151	Acoustics - Noise test code for compressors and vacuum pumps - Engineering method (grade 2)
DIN EN 60204-1	Safety of machinery - Electrical equipment of machines - Part 1: General requirements
DIN EN 61000-6-2 DIN EN 61000-6-4	Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments Part 6-4: Generic standards - Emission standard for industrial environments
DIN EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements
DIN EN 50110-1	Operation of electrical installations
Directive 2012/19/EC	Electrical and electronics - old devices (WEEE)
Directive 2011/65/EC	Dangerous materials in electrical and electronics devices (RoHS)
China - RoHS	Environment protection law - China 2007-03

The following additional safety regulations apply in the FR Germany:

BGV A3	Electrical equipment and operating materials
VBG 5	Power-driven machines
BGR 120	Guidelines for laboratories
BGI 798	Hazard assessment in the laboratory
BGG 919 (VBG 16)	Accident prevention regulations for "compressors"
BGR 189 (BGR 195;192;197)	Use of protective working clothes

Observe the standards and regulations applying in your country when you use the Piston Pump.

Basic Safety Instructions

2 Basic Safety Instructions

2.1 General Information

Warning notices must be observed. Disregarding them may lead to damage to health and property.

The Piston Pump must be operated by personnel who can detect impending dangers and take action to prevent them from materialising.

The manufacturer or authorized workshops will only service or maintain the Pis-ton Pump if it is accompanied by a fully completed damage report. Precise information about the contamination (also negative information if necessary) and thorough cleaning of the Pis-ton Pump are legally binding parts of the contract.

Contaminated Piston Pumps and their individual parts must be disposed of in accordance with the legal regulations.

The local regulations apply in foreign countries.

2.2 Electricity

The Piston Pump of operation mode S1 is supplied. When the location of operation mode S1 devices is changed, please note that the testing must be repeated in accordance with DIN EN 0105, DIN EN 0702 and BGV A2.

The local regulations apply in foreign countries.

Please note the following when connecting to the electrical power supply system:

- The electrical power supply system must have a protective connector according to DIN VDE 0100-410 (IEC 60364-4-41).
- The protective connector must not have any breaks.
- The connecting cable must not be damaged.

2.3 Mechanical Systems

Improper use can lead to injuries or material damage. Observe the following instructions:

- Only operate the pump with hose of the specified dimensions.
- The maximum permissible pressure of 1 bar at the suction connection must not be exceeded.
- Hazardous substances must be separated out as far as this is technically possible before they reach the pump.
- External mechanical stresses and vibrations must not be transmitted to the pump. Only use flexible laboratory hoses for connecting pumps.
- The overpressure generated at the pressure port must not exceed 3 bar.
- The pump must not be used to suck up fluids. Lay the exhaust pipe so that it slopes downwards, so allowing condensate to flow out of the pump. Collect the condensate and dispose of it in an environmentally compatible manner.
- Prevent dyes exuding.
- Maintain a space of least 20 mm between the pump and adjacent parts in order to enable the pump to cool.

	CAUTION !
Solid particles in the pumping medium impair the pumping action and can lead to damage. Prevent solid particles penetrating into the pump.	

2.4 Hazardous Substances

	ACHTUNG !
<p>The operating company bears the responsibility for the use of the Piston Pump. Hazardous substances in the gases to be pumped can cause personal injuries and property damage. Pay attention to the warning notices for handling hazardous substances.</p>	

The local regulations apply in foreign countries.

Combustible, aggressive and explosive Gases

Don't pump combustible, especially aggressive or explosive gases or vapors or operate this pump in an atmosphere containing combustible or explosive gases or vapors.

Examine before switching on whether that can form gas combustible gas/air mixtures which can be promoted! Consider the regulations of the guideline 1999/92/EC.

The Piston Pump is not certified according to ATEX guidelines 94/9/EC.

The Piston Pump is not recommended for pumping acidic, basic or organic vapors.

Poisonous Gases

Use a condensate separator (eg. Woulff bottle) when pumping poisonous or harmful gases. Prevent such substances from leaking out of the appliance or pump. Treat these substances according to the applicable environmental protection regulations.

Test the strength and leak-tightness of the connecting lines and the connected apparatus. Prevent environmental poisons, e.g. mercury, getting into the pump.

Fulfil the requirements, for example:

- German Hazardous Substances Regulation (GefStoffV) of 01. December 2010
- Regulation (EC) 1272/2008
(Classification, Packaging and Labelling of hazardous substances),
- Manufacturer's safety data sheets on hazardous substances.

2.5 High Temperatures

The Piston Pump may heat up as a result of the temperature of the gas being pumped and through compression heat.

Prevent the following maximum permissible temperatures from being exceeded.

- + 40 °C for the environment, and
- + 60 °C for the gas to be pumped.

The electric motor has a thermal overload protection.

Description

3 Description

3.1 Design

3.1.1 Connections - Suction-/Pressure side

3.1.2 Connecting to the electricity supply

The Piston Pump is supplied with complete electrical wiring. It is connected via a mains connection cable and a power plug. Mains connection cable and plug must comply with the requirements of the line disconnection devices (current, output).

Motor is provided with a thermal overload protection ex works, protecting the vacuum pump from damage or destruction, respectively.

	<p>WARNING !</p>
<p>Should the user change the electrical connection, for example for fitting into a system, then this may only be performed by a electrical specialist under observance of the accident prevention regulations.</p>	

	<p>WARNUNG !</p>
<p>Please keep in mind that the barometric pressure changes from day to day. This also changes the displayed pressure values because of latent air pressure fluctuations.</p>	

3.2 Protection measure against liquids in the pump

3.2.1 Condensate separator (Suction side)

If the possibility of penetration of larger liquid quantities we recommend the installation of a larger collection vessel (eg. Woulff bottle).

3.3 Areas of Application

The Piston Pump is intended to:

- Pumping and compressing neutral and aggressive gases and vapours.
- Generating a vacuum down to a minimum ultimate pressure 292 mbar.
- Use in physical and chemical laboratories in trade and industry.
- Use for vacuum filtration, vacuum drying and other vacuum technology applications.

Technical Data

4 Technical Data

Parameter	Unit	CR-KP292
Pumping speed at atmospheric pressure	m ³ / h	0.55
	l / min	9.2
Ultimate pressure	mbar	292
Overpressure, max.	bar (ü)	3.3
Suction-/ Pressure connection	-	Hose nozzle DN 6 for hose inside diameter 6 mm
Ambient temperature	°C	+ 10 to + 40
Max. Operating gas temperature		+ 60
Noise level DIN EN ISO 2151	dB (A)	≤ 45
Voltage / Frequency	V / Hz	230 / 50
Motor power	W	25
Type of protection DIN EN 60529	-	IP 20
Weight	kg	2.0
Dimensions (W/D/H)	mm	194 / 114 / 191

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of products.

It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

Carl Roth GmbH + Co. KG does not warrant, guarantee or assume any obligation or liability in connection with this information.

5 Installation and Operation

5.1 Unpacking

Carefully unpack the **Piston Pump**.

Check the pump for:

- Transport damage,
- Conformity with the specifications of the supply contract (type, electrical supply data),
- Completeness of the delivery.

Please inform us without delay if there are discrepancies between the delivery and the contractually agreed scope of delivery, or if damage is detected.

Please take note of the general terms of business of the manufacturing firm.

In case of a claim under warranty, the device must be returned in packaging that is suitable for protecting it during transport.

5.2 Setting up and connecting

5.2.1 Setting up

- Set the pump on a flat and horizontal surface.
- Remove the protective caps on the suction and pressure ports.
- Note that the cooling of the pump is guaranteed, *see chapter 2.3*.
- Note that on the installation location no moisture acts on the pump.

5.2.2 Connecting

5.2.2.1 Electrical Connection

Before the electrical connection of the pump review the specifications on the nameplate with the existing electrical connection conditions, *see chapter 2.2*.

CAUTION !

**The electric motor has a thermal overload protection.
After triggering the protective fuse after a certain cooling occurs, an auto restart (Reset) the motor.**

5.2.2.2 Vacuum connection (Suction side)

The vacuum connection (suction port) consists of a hose nozzle DN 6.

Suitable vacuum line:

- Vacuum hose 6 x 3 mm and hose clamp

Make sure that the vacuum line is kept as short as possible from the pump to the apparatus.

5.2.2.3 Exhaust connection (Pressure side)

The exhaust connection (pressure port) also consists of a hose nozzle DN 6.

Exhaust line for vacuum- and compressor application:

- Fabric hose 6 x 3 mm and hose clamp

Installation and Operation

	CAUTION !
Use only for compressed air applications suitable conduit systems!	

5.3 Operation

	CAUTION !
Observe the basic safety instructions when using the Piston Pump, <i>chapter 2</i>.	

5.3.1 Start-up

- Make sure that when you start the pump, which may arise in the connection lines condensates, cannot penetrate into the pump.
- Always try to avoid aspiration of foreign body particles!
- The Piston pump is switched on using the power switch.
- It is recommended to let the pump run for a few minutes before use. The warm-up improves the eligibility of the pump.

5.3.2 Decommissioning

- If vacuum applications with special steam load, the pump needs to run after the process with an open vacuum port about 2 minutes and rinsed with atmospheric air. This measure is used to cleaning the pump chamber and minimizes the corrosive attack wetted material pump parts.
- The Piston pump is switched off using the power switch.

5.4 Storage

The pumps are to be stored in a low-dust, interior room within the temperature range from + 5 to + 40 °C and at a relative air humidity < 90%.

Leave the protective elements on the suction and pressure ports. Another equally good protection may be used.

5.5 Scrap Disposal

	CAUTION !
The Piston Pump must be disposed of in accordance with the 2012/19/EC guideline and the specific national regulations. Contaminated Piston Pumps must be decontaminated according to the laws.	

6 Maintenance and Servicing

Piston pump is 100% oil-free.

The piston and cylinder of the pump are lubricated maintenance free.

All bearings are sealed and permanently lubricated.

6.1 Maintenance Performed by the User

	WARNING !
<p>Only perform the work that is described here, and that which is permitted to be done by the user.</p> <p>All other maintenance and service work may only be performed by the manufacturer or a dealer authorized by him.</p> <p>Beware of the pump parts being possibly contaminated by hazardous substances. Wear protective clothing if there is contamination.</p> <p>Renew defective parts, if necessary! Do not clean with compressed air!</p>	

- Check the pump daily for unusual running noises and heat building up on the surface of the pump.
- Check the electrical and vacuum connections daily.

6.2 Maintenance by the Manufacturer

Repairs and maintenance going beyond the extent of the work described *in chapter 6.1* or reconditioning or modification may only be performed by the manufacturer or authorized workshops.

	WARNING !
<p>The user shall be liable for the consequences of an incorrect damage report or a contaminated pump. The statements in the damage report are legally binding.</p>	

Troubleshooting

7 Troubleshooting

During the warranty period, intervention in the Piston Pump and accessory components may only be made by manufacturing firm.

Trouble	Cause	Remedy	
		by:	with:
Pump does not start	No power supply	Qualified electrician	Check electrical installation
	Motor defective	Service workshop	Exchange
Pump does not generate a vacuum or overpressure	Connected apparatus and/or connecting elements leaking	User or Service workshop	Identify and seal the leak, replace the seals and/or hoses if necessary.
	Pump leaking		Exchange the hoses and/or fittings if necessary.
	Pump dirty		General maintenance / Cleaning
	loose screws		Tighten the screws
Running noise	loose screws or connectors		Tighten the screws and / or connectors
Cable	defective and/or brittle	Qualified electrician	Exchange of the cable

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