# Gebrauchsanleitung Instructions for Use / Instruction d'emploi

## Anti-Statik-Set Best.-Nr. 0213.1

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### Safety Regulations

Observe general references to danger and safety regulations, e.g. wearing of protective clothing, safety goggles, and protective gloves during the transfer operation. Use connecting cables only for the purpose of electrostatic discharging, under no circumstances for other electrical connections!

Use bare copper cable only for liquids which do not attack copper chemically. Transfer substances which generate dangerous vapours only in smoke funnel or rooms with good ventilation or outside.

Observe safety advices of the manufacturers of reagents. Observe with combustible liquids the safety regulations/Explosion-protection. The transfer of hydrocarbons into non-conductive plastic containers is limited to a maximum of 2 litres!

### Packaging/Content

1 Plastic case

1 Connecting cable 3.00 m, red

1 Connecting cable 3.00 m, black

1 Connecting cable 1.40 m, bare copper

1 Instruction for Use

Please check the content for completeness.

Damage to the package which can be seen on the outside should be immediately shown to the deliverer/forwarding agent who in turn must confirm it.

### Function

When transferring liquids, electrostatic charges may result, if the transfer equipment and containers are not sufficiently earthed. This can lead with many substances to ignition or explosion.

The operating instructions B210 to B236 and the safety advice S33 of the Ordinance concerning Hazardous Substances and not in the least one's own interest, therefore, require that during the transfer of ether, benzole, benzine, carbon bisulphide and similar liquids a complete earthing must be given.

This can be easily and reliably achieved by the **Anit-Statics Set**.

### Handling

The **Anti-Statics Set** consists of three connecting cables. Two insulated colour coded connecting cables (red and black) are equipped with clamping nippers at each end. The third connecting cable of bare copper possesses one end with clamping nippers and one free end. When transferring liquids, the two conductive containers are being connected by the red cable by clamping; one of the two containers is being earthed by the black cable, e.g. to an earth connected conductor, e.g. water pipe (no plastic tubes).

If it is not possible to fasten the clamping nippers at the container which is to be filled or if this consists of poorly conductive material, e.g. glass, the cable of bare copper with only one clamp is being used for the connection of both containers. In these cases the free end of the cable is immersed in the liquid of the container to be filled. At this point the chemical resistance of copper against the transferred substance and the electric conductivity of the transferred substance have to be checked (e.g. pure hydrocarbons are considered to be non-conductors).

**Important:** The clamps must always be fixed to the bare metal of containers, pumps and earthing connections. If needed, they have to be freed of coatings.

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