# Nitrat

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# Nitrate

## Test kit for performing colorimetric tests on nitrate ions in surface water and sewage

**Method:**

Nitrate ions are reduced to nitrite ions in an acidic medium. Combined with a suitable aromatic amine, these form an orange-yellow azo dye.

**Measurement range:**

1–120 mg/L NO<sub>3</sub><sup>-</sup>

**Contents of test kit (\*refill pack):**

sufficient for 110 tests

- 30 mL NO<sub>3</sub>-1\*
- 5 g NO<sub>3</sub>-2\*
- 1 measuring spoon 70 mm\*
- 2 screw-plug measuring glasses
- 1 slide comparator
- 1 color chart
- 1 plastic syringe 5 mL
- 1 instructions for use\*

**Hazard warning:**

This test does not contain any hazardous substances which must be specially labelled as hazardous.

**Instructions for use:**

also refer to the pictogram on the back of the color chart

1. Pour a **5 mL water sample** into each of the measuring glasses using the plastic syringe.  
Place a measuring glass on position A in the comparator.

**Only add the reagent to measuring glass B.**

2. Add **5 drops of NO<sub>3</sub>-1**, seal the glass and mix.
3. Add **1 level measuring spoonful of NO<sub>3</sub>-2**, seal the glass and **immediately shake the mixture well for 1 min.**
4. Open the glass after **5 min** and place it on position B in the comparator.
5. Slide the comparator until the colours match in the inspection hole on top. Check the measurement reading in the recess on the comparator reed. Mid-values can be estimated.
6. After use, rinse out both measuring glasses thoroughly and seal them.

The reagents can be used for the **photometric evaluation** with photometer PF-12 / PF-12<sup>Plus</sup>.

This technique can be used also for analyzing sea water (see „Conversion table“).

**Disposing of the samples:**

The used analysis specimens can be flushed down the drain with tap water and channelled off to the local sewage treatment works.

**Interferences:**

Depending on their concentration, oxidizing substances may reduce the measurement reading or suppress the reaction totally. Chlorine ≤ 10 mg/L does not interfere.

Nitrite interferes (same reaction). This can be circumvented by addition of amido sulphonic acid (REF 918 973).

The water sample should be between 18 and 30 °C. At lower temperatures the reaction takes place at a significantly slower rate, and the results are limited.

**Conversion table:**

mg/L NO <sub>3</sub> <sup>-</sup>	mg/L NO <sub>3</sub> -N (Nitrate nitrogen)	mmol/m <sup>3</sup>	mg/L NO <sub>3</sub> <sup>-</sup> in sea water
1	0.2	16	1
3	0.7	48	3
5	1.1	81	5
10	2.3	160	12
20	4.5	320	25
30	6.8	480	40
50	11	810	65
70	16	1130	95
90	20	1450	120
120	27	1940	160

**Storage:**

Store the test kit in a cool (< 25 °C) and dry place.

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