# Mangan

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# Manganese

**Test kit for performing colorimetric tests  
on manganese ions in surface water and sewage**

**Method:**

Determination of total manganese with formaldoxime

**Measurement range:**

0.1–1.5 mg/L Mn

**Contents of test kit (\*refill pack):**

sufficient for 70 tests

30 mL Mn-1\*

30 mL Mn-2\*

4 g Mn-3\*

1 measuring spoon 70 mm\*

2 screw-plug measuring glasses

1 slide comparator

1 color chart

1 plastic syringe 5 mL

1 instructions for use\*

**Hazard warning:**

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

**Instructions for use:**

*also refer to the pictogram on the back of the color chart*

1. Pour a **5 mL water sample** (the pH value of the sample must be between pH 1 and 13) into each of the measuring glasses using the plastic syringe.

Place a measuring glass on position A in the comparator.

**Only add the reagent to measuring glass B.**

2. Add **5 drops of Mn-1**, seal the glass and mix.
3. Add **7 drops of Mn-2**, seal the glass and mix. Wait **1 min**.
4. Add **1 level measuring spoonful of Mn-3**, seal the glass and shake well.
5. Open the glass after **5 min** and place it on position B in the comparator.
6. Slide the comparator until the colors match in the inspection hole on top. Check the measurement reading in the recess on the comparator reed. Mid-values can be estimated.
7. After use, rinse out both measuring glasses thoroughly and seal them.

The reagents can be used for the **photometric evaluation** with photometer PF-12/ PF-12<sup>Plus</sup>.

The method can be applied also for the analysis of sea water.

**Disposing of the samples:**

Information regarding disposal can be found in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

**Interferences:**

The following ions will not interfere:

≤ 0,1 mg/L Ni<sup>2+</sup>; ≤ 2 mg/L Cr, Co<sup>2+</sup>; ≤ 20 mg/L Cu<sup>2+</sup>, Fe<sup>3+</sup>;

≤ 50 mg/L Al<sup>3+</sup>, PO<sub>4</sub><sup>3-</sup>

**Storage:**

Store the test kit in a cool (< 25 °C) and dry place.

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