**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

## **ZX4 Advanced Vortex Mixer**

F202A0280

### **General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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### **This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Mixing by spinning can be carried out quickly and safely by placing the tube containing the sample on the mixing cup head. The Vortex Mixer can be used in automatic (SENSOR) or continuous mode (CONT.): in both cases the mixing speed can be adjusted using the speed selector knob (1).

When the instrument is switch on, it is always set in SENSOR mode. Mixing starts when the test-tube crosses the infrared sensor (6) field near the mixing cup head (4). For continuous mode press the mode select button (2), mixing starts automatically at the selected speed, no pressure is required on the cup head.

The led (5) indicates that the instrument is running in SENSOR mode.

The press-on fitting means that the mixing cup head can easily be replaced by any of the various accessories available for test-tubes of different sizes.


Fig. 1

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- ZX4 Advanced IR Vortex Mixer with mixing cup head
- EU adapter plug
- Power supply 100-240V/12V no plug
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up


Check that the mixing cup head is correctly positioned. Also check that the speed regulation knob is turned to "0 rpm". Use the on-off switch on the front of the unit to turn the unit on (position "ON"). The start-up and the control of the speed are effected with the knob. When switched on the unit starts up in "SENSOR" mode with the indicator lights on.

### 2.3 Information about construction materials

Housing	Zinc alloy/Polymer	Feet	MQ/MVQ
Paint	PET	Platform	PP
Frontal control panel	PET	Foamed part for platforms	PE
Cup head	SEBS		

## 3. Operating controls

Connect the unit to mains and turn it on using the on-off switch. Select the operating mode using the mode select button on the front of the unit.

When the **Touch** mode is selected, mixing starts automatically when the test-tube crosses the sensor field (infrared optical system) shown in Figure 2.

In **Continuous** mode the mixing cup head runs continuously at the set speed without pressure on the mixing cup head. However, it is important to grip the tube firmly during mixing.

In both operating modes there is no need for any pressure on the mixing cup head. The sensor field is purpose-designed to detect the presence of any kind of test-tube.

Fig. 2

## CONTINUOUS-SENSOR

Pushing the mode select button, the respective modes are activated:

- **Sensor** (indicator light on) allows automatic operation mode even with different accessories.
- **Continuous** (indicator light off). The instrument runs in continuous mode with various accessories.

## SPEED REGULATION KNOB

The knob on the front of the unit can be used for quick precision regulation of mixing speed from 0 to 3000 rpm.

## ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.

Always turn the unit off after use.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

To replace the cup head, pull it upwards and off. The other push-fit accessories can then be inserted.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Dimensions (WxHxD)	150x130x165 mm (5.9x5.1x6.5 in)	
Weight	2.7 Kg (6.2 lb)	
Speed range	0-3000 rpm	
Spin diameter	4.5 mm (0.177 in)	
Max. weight on the mixing cup head	0.5 kg (1.1 lb)	
Permitted operation	Continuous	
Operation mode	Continuous or Sensor	
Protection degree	IP 42-EN 60529	
Pollution degree	2- EN61010-1	
Admitted temperature	Storage: -10...+60 °C (+14...+140 F°)	Working: 5°...+40 °C (+41...+104 °F)
Admitted humidity	Max. 80%	
Max. altitude	4000 m	
Power	15 W	
Power supply input	AC 100-240V ; 50-60 Hz ; 0.5A	
Power supply output	DC 12V ; 1.25A	

## 6. Accessories / Spare parts

A00000012	Foam stand for 19 microvials 1.5 ml - Eppendorf®
A00000013	Customizable soft foam top
A00000014	Foam stand 5 test tubes Ø16mm
A00000015	Foam stand for microtiter
A00000016	Small rubber supporting plate Ø 50mm
A00000019	Foam stand 4 test tubes Ø29mm

Do not exceed 800 rpm when using one of the above accessories. High speeds may cause the sample to spill.

A00000012	A00000013	A00000014
A00000015	A00000016	A00000019

10005195	Mixing cup head TX, ZX, RX
10005213	Knob 24D blue

10002564	Foot 28Dx5H antivibration
40001099	Switching VELP 100-240V/12V without plug

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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## **Thank you for having chosen VELP!**

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment.

VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

### **Food&Feed Line**

DK and DKL Digesters

JP Recirculating water pump for fumes aspiration

SMS Scrubber

UDK Distillation units

SER Solvent extractors

FIWE Raw fiber extractors

Dietary fiber extractors

Oxitest Oxidation test reactor

NDA 701 Elemental analyzer

Consumables

### **Environment Line**

ECO Thermoreactors

BOD. Determination systems

Refrigerated thermostats and incubators

Flocculators

Overhead mixer

Mineralization unit for trace heavy metals determination

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