InLab®710, 717, 720, 730 and 740 Conductivity Sensors with Automatic Temperature Compensation, with 1 metre cable.

Intended Use and Principles

These products combine a state of the art conductivity sensor with an Automatic Temperature Compensating device (a thermistor in the sensor body) to provide an easy to use, reliable sensor. They utilise a potentiometric method for measuring conductivity using a series of plates or rings formed in or attached to the sensor shaft. The potentiometric method for measuring conductivity is the most accurate method available today.

The InLab®710 is a 4-plate (platinum) sensor designed for measurement in aggressive media.

The InLab®717 is a 4-ring (platinum) sensor designed for use solely with DL series titrators.

The InLab®720 is a 2-plate (platinum) sensor designed for the measurement of low conductivity samples.

The InLab®730 is a 4-ring (carbon) sensor ideal for general purpose use.

The InLab®740 is a 2-pole (stainless steel) sensor designed for the measurement of low conductivity samples.

Operating Instructions

For optimum performance, use the following operating procedures:

1. Remove all packaging material from sensor before use.
2. To prevent carry over from solution to solution, rinse the sensor with distilled water between measurements.
3. When measuring, make sure the solution is above the cell chamber slot, plates or rings.
4. Ensure the cell chamber is free of bubbles when measuring. To reduce air bubbles immerse the sensor at an angle and then raise to the vertical position.
5. Allow sufficient time for the sensor to stabilise when measuring samples of different temperatures. Manual endpointing is advised with these samples.

6. The InLab®710, 717 and 730 sensors are not recommended for low conductivity solutions (<10µS/cm), for these applications use InLab®720 or 740.
A flow cell (ME-51302257) may also be needed for accurate work on low conductivity solutions.
7. Clean the sensor with distilled water after use.
8. If there is a build up of solids inside the cell chamber, carefully remove it with a cotton wool bud soaked in detergent solution, then rinse the sensor with distilled water.
9. Ensure that the cable connectors are kept free of moisture and dirt.
10. For calibration and measurement procedures, refer to the Instruction Manual supplied with the conductivity meter. A booklet, "A Guide to Dissolved Oxygen & Conductivity Measurement", is also available, see the Instruction Manual for details.
11. In unstirred solutions drift effects are possible. In general, measurement using stirred samples is preferred.

Reagents and Supplies

To achieve optimum performance of this Conductivity Sensor, the correct calibration solutions must be used. See the Instruction Manual supplied with the conductivity meter for details.

⚠ Precautions and Limitations

1. To prevent static damage to the meter avoid touching the cell chamber area of the sensor. Always disconnect the sensor from the meter before cleaning.
2. Do not expose the sensor to organic solvents, either when cleaning or when taking measurements.
3. Do not use the sensor outside the recommended temperature range.
4. Specification and performance of this sensor are guaranteed only when used with the appropriate conductivity meter.

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