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## 1. On/Off function **GB**

- Press Cal for 1 sec.  
Calibration is performed
- CAL 2 blinks in the display
- Press briefly
- If 1 appears, deactivate with
- Automatic shut-off after 8 min. of inactivity
- Press and hold for 3 sec. to switch off the device
- Display lighting on/off

## 2. Measuring the pH value

- Hold the device without protective cap in the sample water to be measured so that the device is immersed in the sample water no higher than the seal ring.
- The measurement is displayed.
- Briefly press 1 to freeze/release the measurement
- Press and hold for 3 sec. to save the measurement  
The storage location (e.g. 01) appears in the display while the button is pressed

## 3. Calibration

- The electrode should be cleaned thoroughly with distilled water before it is immersed in another buffer solution.  
Recalibrate as necessary.  
For use of Lovibond buffer solutions or pH buffer tablets, refer to the instructions (see 4.)
- Remove the protective cap, switch on the device, rinse the electrode with distilled water, then immerse in the pH 7 Puffer solution.
- Press Cal for 3 sec.  
CAL 1 blinks in the display
- Wait until goes out
- Calibration is performed
- CAL 3 blinks in the display
- End the 2-point calibration or continue with
- Rinsing of the electrode with distilled water and immersing it in the pH 10 Puffer solution.
- CAL 2 blinks in the display
- Wait until goes out
- Press Cal for 1 sec.  
CAL 3 blinks in the display
- 3-point calibration is finished

- 1. A clean container is filled with 20 ml of deionised / distilled water.
- 2. A buffer tablet is added and dissolved in the measured value.
- 3. After a wait time of approximately two minutes, the active tablet components have dissolved and the insoluble components have settled to the bottom of the container.
- 4. The electrode is immersed in the buffer solution and the calibration is performed (see 3.).
- 5. The electrode is removed from the container and thoroughly cleaned.

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- Selection of the buffer system (7.00 or 6.86)
- Save and close

## 9. Type of display of date and time

- Switch on
- Press for 3 sec.
- 4 x drücken
- Date appears
- Confirm
- Choose between 12-hour clock (AM/PM) and 24-hour clock
- Save and close

## 6. Viewing saved data

- Switch on
- Press for 3 sec.
- The last storage location blinks
- Confirm
- The last storage location is shown
- Switch between storage locations
- Close

## 7. Deleting saved data

- View the saved values (see 6.)
- Press for 3 sec.
- CLEAR blinks
- Confirm deletion of values  
All values are deleted
- Process closes automatically

## 8. Select buffer system

- Switch on
- Press for 3 sec.
- 1 x drücken  
Adjusted buffer system blinks (see 3.).
- Confirm

## Save and switch to year

- To increase
- To decrease

## Save and switch to month/day

- To increase
- To decrease

## Save and close

- Save and close

## 11. Adjusting the temperature unit

- Switch on
- Press for 3 sec.
- 3 x drücken
- Temperature unit blinks
- Confirm
- Selection of temperature unit (°F/°C)
- Save and close

## 10. Adjusting date and time

- Switch on
- Type of display must be selected beforehand (see 9.)
- Press for 3 sec.
- 2 x drücken
- Time blinks
- Confirm

## Switch on

- Confirm

## Selection of temperature unit (°F/°C)

- Save and close

## View the saved values (see 6.)

- Press for 3 sec.

## CLEAR blinks

- Confirm deletion of values  
All values are deleted

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- 
- 
- 
- 
- 
- 

## Error messages and procedure

Description	Measuring range	Procedure
Err 1 Electrode locked incorrectly or defective	Thermistor: 0.5–100 kΩ	Connect correctly or replace electrode
Err 2 pH outside of the measuring range	0–14 pH	Use new buffer solution or replace electrode place electrode in the solution in solution
Err 3 Temperature outside of the measuring range	0–60 °C	Bring measuring solution specified temperature-range
Err 4 Offset outside of the range	-60–60 mV	Use new solution or replace electrode
Err 5 Drift outside of the range	85%–115%, 50–68 (mV/pH)	Use new solution or replace electrode
Err 6 mV outside of the range	-1800–1800 mV	Use new mV solution

## Accessories

- 19 50 021-2 AAA batteries in 2-pack (2 batteries required)
- 19 48 20 pH spare electrode
- 38 48 01 Measuring beaker, 100 ml
- 241 86 09 KCl solution, 3 M, 100 ml
- 19 48 10-16 Plastic box
- 45 70 22 VE water for rinsing the electrode

## Warranty

The warranty for the SD 50 is valid for a period of 2 years from the date of purchase; the warranty for the electrode is 6 months. This does not cover damages resulting from improper handling or based on mechanical damage.

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## Storing the electrode

The electrode must be stored in a 3 M KCl (Order no. 418609) solution for at least 5 hours before the initial use. For this purpose, fill a suitable glass until the glass portion of the electrode is completely immersed. Observe the 'Important information for pH electrodes' for exact details on storage and safekeeping.

## SD 50 pH Specifications

- Temperature range 0-60 °C, 0-14 pH, -1800 mV to +1800 mV
- pH resolution 0.01 pH, Accuracy  $\pm 0.05$  pH
- Buffer system selection (pH 7.00 or pH 6.86): 1-, 2-point or 3-point calibration with automatic recognition
- Automatic temperature compensation
- mV resolution 0.1 mV within  $\pm 1000$  mV and 1 mV outside of  $\pm 1000$  mV, Accuracy  $\pm 20$  mV
- Temperature resolution 0.1 °C, selection of °C or °F
- Display of time and date with 25 storage locations (not permanent)
- 22x22 mm LCD display with yellow-brown display background lighting
- 2 x AAA batteries 1.5 V
- 10Battery capacity > 350 h with continuous use (without display lighting)
- 10Battery capacity > 350 h with continuous use (without display lighting)
- 11 Automatic device shut-off after 8 minutes of inactivity

## Accessories

- 19 50 021-2 AAA batteries in 2-pack (2 batteries required)
- 19 48 20 pH spare electrode
- 38 48 01 Measuring beaker, 100 ml
- 241 86 09 KCl solution, 3 M, 100 ml
- 19 48 10-16 Plastic box
- 45 70 22 VE water for rinsing the electrode

## Warranty

The warranty for the SD 50 is valid for a period of 2 years from the date of purchase; the warranty for the electrode is 6 months. This does not cover damages resulting from improper handling or based on mechanical damage.

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