Instruction Manual Manuale di istruzioni Manuel d'instructions Manual de instrucciones Bedienungsanleitung

Rotilabo®-Mini-Magnetrührer M 3

Roth-# AAN1.1

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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1. Introduction

The M 3 magnetic stirrer is a small, simple and effective stirrer. Thanks to the white cover surface it is specially designed for microtitrations, and it is extremely useful where reliable and small instrument is needed. The stirrer stays cold even after several days of continuous use. This feature makes it highly appreciated in microbiology and biochemistry.

The technopolymer structure is studied and tested in order to give to the instrument an excellent resistance to the attack of chemical agents and to corrosion in general. The stirring speed can be set within a range of 0 - 1100 rpm using the knob on the front panel. The Stirrer led shows when the stirring is running.

NOTE: it's important to choose the most suitable magnetic stirring bar related to the quantity and the liquid to be stirred as well as to the type of the beaker. The magnetic stirring bar that satisfies most of the application is A00001056 (6x35mm).

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- Magnetic stirrer M 3
- Plug for power supply
- Switching VELP 100-240V/12V without plug
- Instruction manual

2.1 Electrical connection

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.


2.2 Start-up

Rotate the speed knob completely to the left. Place the flasks containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed by turning the speed knob.

2.3 Information about construction materials

Structure PA66/6 Feet MQ/MVQ

3. Operating controls

Turn the stirring speed knob on the front panel in order to start the stirring function.
It is possible to select speeds between 0 and 1100 rpm.

When the led is lighted it means that the stirring is on

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Admitted power supply	100-240 V / 50-60 Hz
Power	0.6 W
Dimensions (WxHxD)	130x50x150 mm (5.11x1.96x5.9 in)
Weight (with plug)	0.4 Kg
Speed range	0 - 1.100 rpm
Operating mode	Continuous
Construction material	ABS
Environmental ambient temperature	5 - 40 °C (41 - 104°F)
Storage temperature range	-10... + 60 °C (14 - 140°F)
Max humidity	80%
Pollution degree	2

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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Thank you for having chosen VELP!

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment.

VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

Food&Feed Line

DK and DKL Digesters

JP Recirculating water pump for fumes aspiration

SMS Scrubber

UDK Distillation units

SER Solvent extractors

FIWE Raw fiber extractors

Dietary fiber extractors

Oxitest Oxidation test reactor

NDA 701 Elemental analyzer

Consumables

Environment Line

ECO Thermoreactors

BOD. Determination systems

Refrigerated thermostats and incubators

Flocculators

Overhead mixer

Mineralization unit for trace heavy metals determination

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Carl Roth GmbH + Co. KG