**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

## **REC Heating Plate**

F20700081, F20710081

### **General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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### **This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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The REC heating plate with ceramic plate, is used in the laboratory for general use and for all those applications that require a precise regulation of the heating plate temperature.

The ceramic heating plate (1) means that higher temperatures can be obtained compared to the traditional aluminium heating plate and it takes less time to heat the sample. Ceramic is an inert and very hard material and is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The instrument is turned on using the ON/OFF button (2). The temperature (4) of the heating plate is adjusted using the relative knob on the front panel, the value is indicated on the digital display (3).

**NOTE:** If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

Fig. 1

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- REC Heating Plate
- Power supply
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Rotate the temperature knob completely to the left. Place the flask containing the sample. Then, set the temperature by turning the knob.

## 3. Operating controls

Connect the unit to mains and turn it on using the on-off button. The green led indicates that the instrument is On.

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures of from 5 to 550 °C can be selected at 5 °C intervals. The temperature setting is shown on the display. When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50 °C, the message "Hot" flashes on the display until the temperature of the heating plate falls to below 50 °C (or until the instrument is turned off using the ON/OFF).

**NOTE:** this warning is not active if the instrument is not connected to the power supply.

The knob on the front panel is easily accessible and is well away from the heating plate in order to ensure maximum operator safety as well as safeguarding the electronic components inside the unit.

### 3.1 Error messages

If an error message appears on the display, please contact VELP Scientifica's technical service department.

When the display shows an error message, heating stops automatically.

**NOTE:** To remove the error message disconnect the instrument from the power supply.

- AL1 Thermocouple is out of range
- AL2 Excessive heating time

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Dimensions WxHxD	203x93.5x344 mm (8x3.68x13.5 in)
Weight	3.5 Kg (7.7 lbs)
Overall power	800 W
Dimension of the heating plate	180 x 180 mm
Programmable temperature range	5 – 550 °C
Type of temperature control	Digital
Overtemperature protection	Yes
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Pollution degree CEI EN61010-1	2

## 6. Accessories / Spare parts

A00001071	Hemispheric bowl for 250ml flasks
A00001072	Hemispheric bowl for 500ml flasks
A00001073	Hemispheric bowl for 1000ml flasks
A00001069	Support rod

10002097	Knob
10000239	Foot

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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## **Thank you for having chosen VELP!**

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment.

VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

### **Food&Feed Line**

DK and DKL Digesters

JP Recirculating water pump for fumes aspiration

SMS Scrubber

UDK Distillation units

SER Solvent extractors

FIWE Raw fiber extractors

Dietary fiber extractors

Oxitest Oxidation test reactor

NDA 701 Elemental analyzer

Consumables

### **Environment Line**

ECO Thermoreactors

BOD. Determination systems

Refrigerated thermostats and incubators

Flocculators

Overhead mixer

Mineralization unit for trace heavy metals determination

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