***Instruction Manual  
Manuale di istruzioni  
Manuel d'instructions  
Manual de instrucciones  
Bedienungsanleitung***

**PW - LH - LS - ES Overhead stirrers**

F201A0150, F201A0156, F201A0151, F201A0152

**General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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**This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

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The values indicated on the rating plate of the instrument must correspond to those of the power supply.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Safe working conditions are ensured only when the accessories described in the dedicated chapter are used.

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The working speed set on the instrument must be such as to avoid wobbling and/or splashes.

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Do not use the instrument in explosive environment or with dangerous substances. The unit must not be immersed in water.

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It is dangerous to run the unit with the stirring blade turning in free air. Always place the stirring shaft in the receptacle before turning the unit on.

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## Contents / Indice / Index / Índice / Inhalt

1.	INTRODUCTION.....	5
2.	ASSEMBLY AND INSTALLATION.....	5
	2.1 ELECTRICAL CONNECTION.....	5
	2.2 START-UP.....	5
	2.3 INFORMATION ABOUT CONSTRUCTION MATERIALS.....	5
3.	OPERATING CONTROLS.....	5
	3.1 OVERLOAD PROTECTION.....	5
4.	MAINTENANCE.....	6
	4.1 CLEANING.....	6
5.	TECHNICAL DATA.....	6
6.	ACCESSORIES / SPARE PARTS.....	6

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
Overhead stirrers are generally used for laboratory applications where the solubilization, mixing, suspension, emulsifying and homogenization of samples is required.

The PW, LH, LS and ES overhead stirrers with electronic speed control and advanced safety features are able to satisfy the most diverse needs of the laboratory in terms of both viscosity and volume.

The chuck is keyless and ensures a secure and easy-to-fit mount for stirring shafts of up to 10mm in diameter. The stirring speed is electronically controlled. The shaft bearings and belt drive ensure silent running.

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Overhead stirrer
- Instruction manual
- Allen key
- Allen head screw
- Shaft
- Power supply cable

### 2.1 Electrical connection

After unpacking the instrument, place the unit on the laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Assemble the shaft using the Allen head screw. Fasten the instrument securely to the support rod (A00001300) with the double clamp (A00001301). Secure the receptacle using the ribbon clamp (A00001302) to avoid damage.

### 2.3 Information about construction materials

Housing	Die casting Polymer	Structure	Aluminum	Control panel	PET
---------	---------------------	-----------	----------	---------------	-----

## 3. Operating controls

Turn the unit on using the mains switch. Turn the potentiometer knob to be select the stirring speed.

Stirring speeds are indicated on the analogue scale around the knob.

**NOTE:** The speed setting is saved when the instrument is turned off or if there is a power failure.

When the instrument is turned on again, stirring starts at the previously set speed.

### 3.1 Overload protection

The current and voltage are electronically limited in order to safeguard the motor and the electronic parts against electric surges.

In the event of abnormal conditions arising the unit shuts down automatically.

**Over Load** Current surges

The stirrer can compensate for eventual current surges for a limited period of time.

If the situation persists, the speed of the motor is gradually reduced, when a speed of 0 is reached the unit shuts down.

**Over Curr** Rpm error

This error occurs if the chuck is accidentally blocked or if the number of revs exceeds the max. permitted value.

**Over Temp** Temperature error

This error occurs if the unit is used for heavy-duty work for prolonged periods of time and at an ambient temperature that exceeds the max. permitted temperature.

Before resetting the instrument allow it to cool down.

**Over Volt** Voltage surge

The software limits the power supply of the motor to the nominal value of the motor itself (60V).

A high speed setting or a heavy load on the chuck may prompt a higher voltage. In this case the voltage is limited to the nominal value of the motor (60 Volt) and the number of revs is reduced accordingly.

If for any reason the power supply of the motor exceeds the maximum threshold value of 70V, "OVERLOAD" the unit will shut down.

**NOTE:** To reset the unit remove the load from the chuck and turn the unit off and on again using the mains switch.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.

## 5. Technical data

Code Model	F201A0150 PW	F201A0156 LH	F201A0151 LS	F201A0152 ES
Admitted power supply	From 80 to 260Volt 50/60Hz	From 80 to 260Volt 50/60Hz	From 80 to 260Volt 50/60Hz	From 80 to 260Volt 50/60Hz
Max. input / output power	190 W / 435 W	190 W / 280 W	120 W / 180 W	30 W / 45 W
Dimensions (WxHxD)	80x230x196 mm (3.1x9.0x7.7 in)	80x230x196 mm (3.1x9.0x7.7 in)	80x215x196 mm (3.1x8.5x7.7 in)	80x160x200 mm (3.1x6.3x7.9 in)
Weight	2,9 Kg (6,4 lbs)	2,9 Kg (6,4 lbs)	2,3 Kg (5 lbs)	1,3 Kg (2,86 lbs)
Speed range at nominal load	20 – 1200rpm	50 – 2000rpm	50 – 2000rpm	50 – 1300rpm
Speed indicator	Analogue	Analogue	Analogue	Analogue
Max torque stirrer shaft	120Ncm	80Ncm	40Ncm	15Ncm
Max stirring quantity H <sub>2</sub> O	70 litres	40 litres	25 litres	15 litres
Clamping chuck clamping range	From 1 to 10mm	From 1 to 10mm	From 1 to 10mm	From 1 to 10mm
Overload and overtemperature protection	Automatic cut-off	Automatic cut-off	Automatic cut-off	Automatic cut-off
Enviromental ambient temperature	5 - 40°C (41 - 104°F)	5 - 40°C (41 - 104°F)	5 - 40°C (41 - 104°F)	5 - 40°C (41 - 104°F)
Storage temperature range	-10... + 60°C (14 - 140°F)	-10... + 60°C (14 - 140°F)	-10... + 60°C (14 - 140°F)	-10... + 60°C (14 - 140°F)
Max humidity	80%	80%	80%	80%
Operating mode	Continuous	Continuous	Continuous	Continuous
Electr. safety level CEI-EN60529	IP40	IP40	IP40	IP40
Pollution degree / installation cat. CEI-EN61010-1	2	2	2	2
Max altitude installation	4000m (13123ft) asl	4000m (13123ft) asl	4000m (13123ft) asl	4000m (13123ft) asl
Fuse diam. 5x20mm	2x3A	2x3A	2x3A	2x3A

## 6. Accessories / Spare parts

A00001300 Support rod and base  
 A00001301 Double clamp  
 A00001302 Ribbon clamp  
 A00001304 Stirring shaft with floating blade  
 A00001305 Stirring shaft with folding blade  
 A00001306 Stirring shaft with fixed blade  
 A00001307 Stirring shaft with propeller  
 A00001308 Stirring shaft with paddle, 6 holes  
 A00001309 Stirring shaft with turbine  
 A00001310 Stirring shaft with turbo propeller  
 A00001311 Stirring shaft with anchor

			
A00001300	A00001301	A00001302	A00001304
			
A00001305	A00001306	A00001307	A00001308
			
	A00001309	A00001310	A00001311

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)

2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2006/95/EC / Bassa tensione 2006/95/CE / Basse tension 2006/95/CE / Baja tensión 2006/95/CE / Niederspannung 2006/95/EG

- Electromagnetic compatibility directive 2004/108/EC / Compatibilità elettromagnetica 2004/108/CE / Compatibilité électromagnétique 2004/108/CE / Compatibilidad electromagnética 2004/108/CE / Elektromagnetische Verträglichkeit 2004/108/EG

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment.

VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

### Food&Feed Line

DK and DKL Digesters  
JP Recirculating water pump for fumes aspiration  
SMS Scrubber  
UDK Distillation units  
SER Solvent extractors  
FIWE Raw fiber extractors  
Dietary fiber extractors  
Oxitest Oxidation test reactor  
NDA 701 Elemental analyzer  
Consumables

### Environment Line

ECO Thermoreactors  
BOD. Determination systems  
Refrigerated thermostats and incubators  
Flocculators  
Overhead mixer  
Mineralization unit for trace heavy metals determination  
Turbidimeter  
Radiation detector

### Stirring Line

Heating magnetic stirrers  
Vertex digital thermoregulator  
Ultraflat magnetic stirrer  
Magnetic stirrers  
Overhead stirrers  
Heating plates  
Vortex mixers  
Homogenizer

### Other Lab Solution

Open circulating bath  
Recirculating water vacuum pump  
Peristaltic pump

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