Bleiacetatpapier

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Lead acetate paper

Lead acetate paper allows the quick and easy detection of hydrogen sulfide (H₂S). The gas occurs in the processing of raw oil. It is toxic even in low concentrations. Sulfide (S²⁻) containing solutions also give a positive reaction.

Safety precautions:

Lead acetate paper contains lead(II) acetate 0.47–3.9% CAS 6080-56-4. DANGER May damage the unborn child. Do not handle until all safety precautions have been read and understood. IF exposed or concerned: Get medical advice/attention. Store locked up. For further information, please ask for safety data sheet (see www.mn-net.com/MSDS).

Color reaction:

The white paper reacts with H₂S to form black lead sulfide. With lower concentrations, however, the white paper assumes only a brown color.

Limit of sensitivity:

1 drop of a solution with 5 mg/L sulfide (S²⁻) gives a brown ring which is just visible.

Disposal:

Please observe local regulations for collection and disposal of hazardous waste and contact waste disposal company, where you will obtain information on laboratory waste disposal (waste code number 16 05 06). Close container tightly. Dispose of contents / container to regulated waste treatment.

Storage:

Avoid exposing test paper to sunlight and moisture. Keep cool and dry (storage temperature 4–30°C).

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