# IKA

designed for scientists

**IKA Plate (RCT digital)**

8

/// DEUTSCH

62

/// ENGLISH

116

/// ESPAÑOL

170

/// FRANÇAIS

224

/// ITALIANO

278

/// PORTUGUÊS

332

/// POLISH

386

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440

/// TÜRKÇE

494

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ENGLISH

### Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the standards or standardized documents: EN 61010-1, EN 61010-2-010, EN 61010-2-051 and EN 61326-1.

## 68

### /// DEVICE SETUP

- › IKA Plate (RCT digital) . . . . . 68
- › Display . . . . . 69

## 70

### /// SAFETY INSTRUCTIONS

- › About our warning symbols . . . . . 70
- › General information . . . . . 71
- › Device design . . . . . 71
- › Permissible medium / contaminants / side reactions . . . . . 72
- › Procedures during sample runs . . . . . 73
- › Accessories . . . . . 74
- › Power supply / switching off . . . . . 74
- › About protecting your device . . . . . 75

## 76

### /// PROPER USAGE

- › Basics . . . . . 76

## 77

### /// UNPACKING

- › Scope of delivery . . . . . 77

## 78

### /// ASSEMBLING THE SUPPORT ROD

- › Process . . . . . 78

## 80

### /// START-UP

- › Operation . . . . . 80
- › Setting the safety temperature limit . . . . . 82
- › Controlling the medium temperature limit using a contact thermometer . . . . . 84

## 86

### /// MENU

- › Menu structure . . . . . 86
- › Operating mode . . . . . 90
- › Temperature control mode . . . . . 91
- › Intermittent mode / rotation direction . . . . . 92
- › Beep . . . . . 92
- › Temperature unit . . . . . 92
- › Error 5 . . . . . 93
- › Expiration of timer . . . . . 93
- › Safe / Visc. . . . . 94
- › Factory reset . . . . . 95
- › External probe calibration . . . . . 95
- › Software version . . . . . 97
- › Timer / Counter (F, 4) . . . . . 97
- › Lock (G, 15) . . . . . 100

## 101

### /// INTERFACES AND OUTPUTS

- › USB and RS 232 interfaces . . . . . 101

## 106

### /// MAINTENANCE AND CLEANING

- › Proper handling . . . . . 106

## 107

### /// ACCESSORIES

- › An overview . . . . . 107

## 109

### /// ERROR CODES

- › Troubleshooting correctly . . . . . 109

## 113

### /// WARRANTY

- › Regulations . . . . . 113

## 114

### /// TECHNICAL DATA

- › IKA Plate (RCT digital) in detail . . . . . 114

# Device setup

/// IKA Plate (RCT digital)

I  
H  
C  
E  
F  
G  
A  
D

N

O

B A


J K L M

## Explication of symbols

<b>A</b>	Power button (left » on, right » off)	<b>I</b>	Standby LED
<b>B</b>	Adjustable safety circuit	<b>J</b>	Power socket
<b>C</b>	Rotating knob „temperature“	<b>K</b>	USB interface
<b>D</b>	Rotating knob „speed“	<b>L</b>	RS 232 interface
<b>E</b>	„Menu“ button	<b>M</b>	Connection for PT1000 temperature sensor series, temperature probes or contact plugs
<b>F</b>	„Timer“ button	<b>N</b>	Heating plate
<b>G</b>	„Lock“ button	<b>O</b>	Threaded hole for stand
<b>H</b>	Display		

/// Display

## Display

<b>01</b>	Set speed value	<b>10</b>	Operating mode
<b>02</b>	Set temperature value	<b>11</b>	Safety circuit temperature
<b>03</b>	Motor running / rotation direction	<b>12</b>	Viscosity trend value
<b>04</b>	Timer / counter	<b>13  </b>	PT1000 temperature sensor connected
<b>05</b>	Temperature value	<b>14</b>	ETS-D5 / ETS-D6 connected
<b>06</b>	Heating function activated	<b>15</b>	All parameters locked
<b>07</b>	Actual hotplate temperature	<b>16</b>	USB connection with PC
<b>08</b>	Actual speed value	<b>17</b>	Connection with PC established
<b>09</b>	Temperature control mode	<b>18</b>	Intermittent mode activated

# Safety instructions

## /// About our warning symbols

Indicates an imminently hazardous situation, which, if not avoided, **will result in death, serious injury.**

Indicates a potentially hazardous situation, which, if not avoided, **can result in death, serious injury.**

Indicates a potentially hazardous situation, which, if not avoided, **can result in injury.**

Indicates handling, which, if not avoided, **can result in equipment damage.**

**Attention!**  
Indicates the risks of magnetism.

**Danger!**  
Indicates the exposure to a hot surface.

## /// General information

- › **Read the operating manual in its entirety before using the IKA Plate (RCT digital) and follow the safety instructions.**
- › Keep the operating manual in a location that is easily accessible by everyone in the laboratory.
- › Ensure that only trained staff operates the device.
- › Follow the safety instructions, (laboratory) guidelines, occupational health and safety and accident prevention regulations (OSHA).
- › Electrical outlet must be grounded (protective ground contact).

### **Attention – Magnetism!** (📵)

Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).

### **Risk of burns!** (🔥)

Exercise caution when touching parts of the housing and the heating plate. The heating plate can reach dangerous temperatures. Pay attention to the residual heat on the heating plate after switching off the stirrer.

**The device must only be transported when the heating plate has cooled off and is no longer hot to the touch!**

## /// Device design

### **Danger!** (⚡)

**Do not** operate the device in explosive atmospheres, with hazardous substances or under water.

- › Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- › The feet of the device must be clean and undamaged.
- › Ensure that the power cord does not touch the heating plate.
- › Check the device and accessories for damage before each use. Do not use damaged components.

**Caution! (⚠)**

Device surface is partially made of glass!

- › Glass surface can be damaged by impact.
- › If glass surface is damaged it could cause injury, don't use the device anymore and contact your IKA Service team.

### /// Permissible medium / contaminants / side reactions

**Caution! (!)**

Only process and heat media that has a flash point higher than the adjusted safe temperature limit that has been set. The safety temperature limit must always be set to at least +25 °C lower than the fire point of the media used.

**Warning! (!)**

Beware of hazards due to:

- › flammable materials.
  - › combustible media with a low boiling temperature.
  - › glass breakage.
  - › incorrect container size.
  - › overfilling of media.
  - › unsafe condition of container.
- 
- › Process pathogenic materials only in closed vessels under a suitable fume hood.

Please contact IKA if you have any questions.

**Danger! (☒)**

Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.

- › The heating plate can heat up due to the action of the magnets at high motor speeds, even if the heater is switched off.
- › Please consider any possible contaminations and unwanted chemical reactions.
- › It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- › When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C – 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

*(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)*

### /// Procedures during sample runs

**Warning! (!)**

Wear your personal protective gear in accordance with the hazard category of the media to be processed. There may be a risk from:

- › splashing and evaporation of liquids.
- › ejection of parts.
- › release of toxic or combustible gases.

Reduce speed if:

- › medium splashes out of vessel because the speed is too high.
- › device is not running smoothly.
- › container moves on the heating plate.
- › an error message is displayed.

### /// Accessories

- › Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- › Always disconnect the plug before attaching accessories.
- › Observe the operating instructions of the accessories.
- › Ensure that the external temperature sensor is inserted into the medium to a depth of at least 20 mm when connected.
- › Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.

### /// Power supply / Switching off


- › The voltage stated on the type plate must correspond to the mains voltage.
- › The outlet for the mains plug must be easily accessible.
- › The device can only be disconnected from the mains outlet by pulling out the power cord or the connector plug.

#### **Warning!** (▲)

- › The device will automatically restart in **mode B** following any interruption to the power supply.

### /// About protecting your device

- › The stirrer may only be opened by qualified and IKA approved experts.
- › Do not cover the device, even partially e.g. with metallic plates or film. This may result in overheating.
- › Protect the device and accessories from bumps and impacts.
- › Keep base plate clean.
- › Observe the minimum distances:  
between devices min. 100 mm,  
between device and wall min. 100 mm,  
above the device min. 800 mm.

# Proper usage

## /// Basics

### Usage

- › The magnetic stirrer has been designed for mixing and / or heating liquids.

### Area of use (only indoors)

- › laboratories
- › pharmacies
- › universities
- › schools

### This device is suitable for use in all areas except

- › residential areas.
- › areas that are connected directly to a low-voltage supply network that also supplies residential areas.

### The safety of the user cannot be guaranteed

- › if the device is operated with accessories that are not supplied or recommended by the manufacturer.
- › if the device is operated improperly or contrary to the manufacturer's specifications.
- › if the device or the printed circuit board are modified by third parties.



# Unpacking

## /// Scope of delivery

Unpack the device carefully. **Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).**

### Scope of delivery

- › IKA Plate (RCT digital)
- › power cord
- › operating manual
- › warranty card
- › screwdriver (safety circuit)
- › USB cord
- › magnetic stir bar **IKAFLOX 30** and **40 mm**
- › temperature sensor **PT1000.60**

# Assembling the support rod

/// Process


## Note! (▲)

Review the mounting and safety instructions of the IKA boss head clamp prior to using it. (see "Accessories")

**The device must not be suspended from the support rod!**


## Procedure

1. Remove the threaded plug in the threaded hole (O).
2. Remove the protective cap from the support rod.
3. Place the washer between the housing and the nut.
4. Screw in the support rod manually until it no longer can be tightened.
5. Tighten the nut using an SW17 open-end spanner.
6. Assemble the accessories using a boss head clamp.


## Note! (▲)

When using bath attachments with a diameter over 180 mm, use support rod in conjunction with an extension. (see "Accessories")

## Risk of tipping! (▲)

Make sure that the center of mass of the attached device does not protrude over the safe area indicated below with a dotted rectangle.

# Start-up

## /// Operation

Remove the protective film from the plate before using device!

### Power-up

1. Press the power button **(A)** to the right.
2. Insert the power cord into the power outlet **(J)**.
3. Press the power button **(A)** to the left.
  - › *Standard mode A is selected.*

### Stirring

1. Set the speed using the rotating knob **(D)**.
  - › *The set speed value (01) is shown on the display (H).*
2. Press the rotating knob **(D)** to start the stirring function.

### Heating

1. Set the safety temperature limit **(11)**.
2. Set the target temperature using the rotating knob **(C)**.
  - › *The set temperature value (02) is shown on the display (H).*
3. Set the temperature control mode **(09)**.
4. Press the rotating knob **(C)** to start the heating function.

### General heating information

- › The target and current temperature **(02, 07)** are shown permanently on the display **(H)**.
- › When the heating function is switched on, the heater activated symbol **(06)** is displayed.
- › When the device is switched off while the heating plate temperature is higher than +50 °C, the display **(H)** shows **HOT!** and the current temperature **(07)** is shown, even if the device is switched off.

### Connecting an external thermometer

1. Press the power button **(A)** to the right.
2. Connect the safety contact thermometer compliant with DIN 12878 Class 2 or temperature sensor PT1000 (individual sensor) to the connection **(M)**.
3. Press the power button **(A)** to the left.

### Temperature sensor PT1000

- › The actual temperature **(07)** shown on the display **(H)** corresponds to the temperature of the medium. The symbol **(13)** is shown on the display.

### Contact thermometer ETS-D5 / ETS-D6

- › Follow the instructions for the contact thermometer. The actual temperature **(07)** is shown on the display **(H)**. The symbol **(14)** is shown on the display.

### Menu

1. Press the "Menu" key **(E)** for to enter the menu mode.
2. Set the menu items by turning the rotating knob **(D)** and confirm by pushing it.
3. Press the "Menu" key **(E)** to enter the menu only possible if stirring and heating is switched off.
4. Press the "Menu" key **(E)** to exit the menu from each state of the Settings menu.

### /// Setting the safety temperature limit

The maximum achievable heating plate temperature is restricted by an adjustable safety temperature limit. Once this limit has been reached, the device stops heating.

#### Caution! (▲)

The denominated temperatures should always refer to the center of the hotplate.


#### Warning! (▲)

The safe temperature limit must always be set at least +25 °C lower than the fire point of the media to be processed!

The adjustable maximum heating plate temperature must always be set at least +15 °C under the set safety temperature limit.

**Setting range:** [+50 °C] to [max. SET temperature +60 °C]

**Factory setting:** [max. SET temperature +60 °C]

#### Setting the safety temperature limit

After switching on the device, the safety temperature limit (**11**) can be adjusted using a screwdriver. Do not turn the setting screw beyond the clockwise or anticlockwise stop. This will cause irreparable damage to the potentiometer.

1. Press the power button (**A**) to the left.
2. Using the screwdriver supplied, turn the "Safe Temp" setting screw (**B**) to the clockwise stop.
3. Turn the temperature rotary knob (**C**) to set the target temperature to the desired "Safe Temperature" and push it to start the heating function, then wait until the temperature is achieved.
4. Turn the "Safe Temperature" setting screw (**B**) slowly anticlockwise until the heating function switches off and the display shows „Err. 25“.
5. The "Safe Temp" value (**11**) is displayed on the screen.

#### Warning! (▲)

Setting the safety temperature limit has to be done as described. The value "Safe Temperature" shown on the display only serves for visualization.

#### Functional test: safety circuit switch-off

- › Heat the device to a temperature of over +100 °C.
- › Set the safety temperature limit as far to the left as possible.
- › The display will show: Err. 25

## /// Controlling the medium temperature limit using a contact thermometer

The preferable method for controlling the average temperature is with contact thermometer. After the set point temperature has been adjusted, this results in a short heating-up time, practically no temperature drift and only minor fluctuation in temperature.

A 6-pin jack is located on the rear side of the device for connecting the PT1000 series, contact thermometer or the contact plug. The electronics of the devices returns a test current that must flow via connector pins 3 and 5 for the heating plate to heat up.

### Safety contact thermometers

Acc. to DIN 12 878 class 2 or acc. Gerstel are connected with a three-wire cable, the test current flows through the contact thermometer.

### Safety function

If the test current is interrupted because of e.g. breakage of contact thermometer or falling out of the cable plug, the heating cuts off.

### Contact thermometer without safety circuit

Acc. to DIN 12 878 class 0. The device only heats if the test current circuit is closed by an electrical connection of the plug pins 3 and 5.

### 2-wire connecting cables

Connect plug pins 3 and 5 of the device plug.

### 3-wire connecting cables

Here the test current circuit can also be produced in the terminal head of the contact thermometer (connect plug pins 2 and 3).

A 3-wire cable with the required bridge is available (accessories).

### Settings

For detailed instructions for settings and limit values, please refer to the operating manual of the instrument you are connecting.


**The desired average temperature can be adjusted on the contact thermometer. The required surface temperature of the heating plate can be selected with the rotating knob (temperature) or button.**

Adjusting the temperature of device to the maximum adjustable temperature will result in the fastest possible heating time, but the average temperature may fluctuate to values above the set-point temperature on the contact thermometer. By adjusting the rotating knob (temperature) or button to approximately twice the set-point value of contact thermometer (with a set-point of +60 °C, the temperature of device would be set to +120 °C), you will reach a good compromise between a fast heating time and over-shooting the set point. If you adjust the temperature of the device to exactly the set-point temperature, the medium will not reach the set-point temperature because some loss of the heat will always occur between the heating plate and the medium.

**If a fault occurs in the control circuit, the maximum heating plate temperature can be adjusted with the safety temperature setting screw.**

# Menu

/// Menu structure

### /// Operating mode

Operating the device in mode A, B or D.

#### Mode A

All settings will be stored if the device is switched off or disconnected from the power supply. The stirring and heating functions will be set to OFF when the device is powered on. The safety circuit can be set or modified. When the power supply switch is switched on, A is displayed.

#### Mode B

All settings will be stored if the device is switched off or disconnected from the power supply. The stirring and heating functions will be set to ON or OFF when the device is powered on, depending on the previous status of the device. The safety circuit can be set or modified. When the power supply switch is switched on, B is displayed.

#### Mode D

##### Switch the device on

The safe temperature has to be confirmed by pushing one of the rotating knobs (**C or D**).

##### During operation

In Mode D the device behaves the same as in Mode A with the exception that:


- › The temperature / speed setting must be confirmed by pressing the temperature / speed knob. To change the temperature / speed setting, turn the temperature / speed knob until the desired value appears.
- › The new value flashes for 5 seconds on the display. Confirm the new temperature / speed setting by pressing the temperature / speed knob, otherwise the temperature / speed setting value will skip back to its previous value.

### /// Temperature control mode

When using an external PT1000 temperature sensor, you can choose between two types of control:


#### PID mode

Good control results, minimized overshooting, slow rise in temperature.

#### 2P-mode (two-point controller)

Maximum heating rate, increased overshooting.

## /// Intermittent mode / rotation direction

### Intermittent mode

Intermittent mode consists of regular, cyclical interruption of the stirring function. The run time (0 / 10 sec – 10 min, 10 sec steps) and pause time (3 sec – 5 min, 1 sec steps) should both be set. When the intermittent mode is activated, a decimal point in front of the first letter appears on the display **(18)**.

### Rotation direction

You can select the continuous rotation. Therefore the time has to be set to 00:00:00 by pushing the rotating knob. **(D)**. Afterwards the rotation direction (clockwise or counterclockwise) can be selected by turning the rotating knob **(D)**. Alternative to the continuous rotation, you can set the automatic reverse rotation function. Then the rotation direction changes after each pause. In this way the magnetic rod is captured and any liquid still rotating is gradually slowed down. Then the drive accelerates to the target speed.

## /// Beep

The menu item **“Beep”** allows the user to turn on / off the sound that occurs when the timer reaches 00:00:00 or an error message appears.

## /// Temperature unit

The menu item **“Temperature unit”** allows the user to select the desired unit between “°C” or “°F” for displaying the temperature.

## /// Error 5

Providing safety protection, Error 5 is taken to recognize that the sensor has not been placed in the medium when the heating is switched on.

### Note! (▲)

You can set a value from 1 to 30 min for this time limit depending on the application. If time limit is set to 0 min, the Error 5 detection is disabled.

### Caution! (▲)

This function will only be active if:

- › sensor temperature is  $< +50\text{ °C}$
- › difference target temperature / sensor temperature  $> 5\text{ K}$

## /// Expiration of timer

You can specify the procedure of heating and stirring function if the timer achieves 00:00:00. Next to the visually / audible information you can decide independently if the heating and stirring should be stopped automatically if the timer is expired or should continue.

### Caution! (▲)

If heating and stirring function is set to stop after timer is expired, you have to ensure that there is no heat distortion for this application. If heat distortion is not desirable, the stirring function should continue.

### /// Safe / Visc

Once the "VISC / SAFE" has been set to "VISC", the viscosity trend value (12) will be shown on the display. The torque trend measurement is used to deduce the change in viscosity of the reaction medium. The devices are not designed to measure absolute viscosity. They only measure and display the relative change in the viscosity of the medium from a starting point specified by the user.

Once the speed of the motor and magnetic rod in the medium have stabilized to the target speed, the viscosity measurement is started at 100 %. Consequently, an increase in the measured torque trend is indicated by a value > 100 % and a reduction in values < 100%.

The value can always be reset to 100 % by pushing the „Menu“ button (E) for 2 seconds.

#### **Note!** (▲)

**Torque trend measurement only works for a constant set speed for the duration of the measurement. As a result, intermittent mode cannot be used in conjunction with torque trend measurement.**

The current control variable is saved as the reference 100 % ΔP and shown on the digital display. The change in the viscosity is then shown in %.

*When on the "VISC" screen, the screen will switch to "SAFE" if the adjustment of the safety temperature is activated, during which you may set the safety temperature and the screen will switch back to "VISC" once the setting is done, otherwise it will switch back to "VISC" after 2 seconds.*

### /// Factory reset

The menu item "**Factory reset**" allows the user to reset all the system settings to the original standard values set at dispatch from the factory (see "Menu structure").

### /// External probe calibration


In order to use tolerances to limit temperature deviations, the user can calibrate the temperature sensor in conjunction with the device. If the calibration is not set at default (customer-specific settings), the user can simply reset the calibration. If the calibration is set at default (factory setting), the user can use below method to calibrate the temperature sensor.

#### **2-point calibration**

Calibration using two temperatures.


#### **Note!** (▲)

Before starting the calibration process, please select an appropriate speed value according to the application environment. Compare the actual temperature with a calibrated reference measuring instrument.

### /// Software version

Turn the knob **(D)** to change the menu item to "Software version".

### /// Timer / Counter (F, 4)

Press the key **(F)** for 2 seconds to activate the timer / counter function.

#### Timer mode (count down)

Turn the knob **(D)** to select "hh" value. Press the knob **(D)** to accept "hh" value and the segment jumps automatically to "mm" value.

Turn the knob **(D)** to select "mm" value. Press the knob **(D)** to accept "mm" value and the segment jumps automatically to "ss" value.

Turn the knob **(D)** to select “ss” value. Press the knob **(D)** to accept “ss” value and the system is now in “Timer” mode. Timer is now ready to start.

Press the key **(F)** to start the timer, and then the timer is running and the time is counting down.

Press the key **(F)** to stop the timer (take a break).

Press the key **(F)** to re-start the timer.

If the timer achieves 00:00:00, the display starts to flash and beep (depending on the menu setting).

### Counter mode (count up)

To select the counter mode, you must set all the values to 00:00:00.

Press the knob **(D)** to accept “hh” value and the segment jumps automatically to “mm” value.

Press the knob **(D)** to accept “mm” value and the segment jumps automatically to “ss” value.

Press the knob **(D)** to accept “ss” value and the system is now in “Counter” mode. Counter is now ready to start.

Press the key **(F)** to start the counter, and then the counter is running and time is counting up.

Press the key **(F)** to stop the counter (take a break).

Press the key **(F)** to re-start the counter.

If the elapsed time exceeds the value of 100 hours, the display switches from **hour:minute:second mode** to **day:hour mode**.

If the elapsed time exceeds the value of 100 days, the counter will reset to 00:00:00.

### Note! (▲)

From each state of the timer / counter, press the key **(F)** for 2 seconds to exit the timer / counter.

## /// Lock (G, 15)

Press the "Lock" key (**G**) for 2 seconds to lock / unlock the device.  
When the device is locked, if you press any other buttons (e.g. key **E**, **F**) or turn any knobs (e.g. knob **C**, **D**), the Lock button symbol (**15**) will flash 5 times.

The device keeps the status of "unlock" or "lock" when the power is switched off and switched on again, which is important for operating mode B.

## Interfaces and outputs

### /// USB and RS 232 interfaces

The device can be operated by computer via an RS 232 or USB interface using the laboratory software labworldsoft®.

#### Note! (⚠)

Please comply with the system requirements together with the operating instructions and help section included with the software.

#### USB interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging).

Connected devices and their properties are automatically recognized. Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and also to update the firmware.

#### USB device drivers

First, download the latest driver for IKA devices with USB interface from [www.ika.com/ika/lws/download/usb-driver.zip](http://www.ika.com/ika/lws/download/usb-driver.zip) and install the driver by running the setup file. Connect the IKA device through the USB data cable to the PC. The data communication runs through a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

#### Serial RS 232 interface

##### Configuration

- › The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- › For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- › Transmission process: asynchronous character transmission in start-stop operation.
- › Transmission type: full duplex.

- › Character format: character composition according to data format in DIN 66022 for start-stop operation. 1 start bit, 7 character bits, 1 parity bit (even), 1 stop bit.
- › Transmission speed: 9,600 bit / s.
- › Data flow control: none.
- › Access procedure: data transfer from the device to the computer takes place only at the computer's request.

### Command syntax and format

#### The following applies to the command set:

- › Commands are generally sent from the computer (Master) to the lab device (Slave).
- › The lab device only sends at the computer's request. Even fault indications cannot be sent spontaneously from the lab device to the computer (automation system).
- › Commands are transmitted in capital letters.
- › Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- › Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 80 characters.
- › The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific IKA commands serve only as low level commands for communication between the lab device and the PC. With a suitable terminal or communications program these commands can be transmitted directly to the lab device. The IKA software labworldsoft®, provides a convenient tool for controlling lab device and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the IKA equipment.

NAMUR Commands	Function
IN_NAME	Read the device name
IN_PV_1	Read actual external sensor value
IN_PV_2	Read actual hotplate sensor value
IN_PV_4	Read stirring speed value
IN_PV_5	Read viscosity trend value
IN_SP_1	Read rated temperature value
IN_SP_3	Read rate set safety temperature value
IN_SP_4	Read rate speed value
OUT_SP_1	Adjust the set temperature value
OUT_SP_1 100	Set temperature value
OUT_SP_4	Adjust the set speed value
OUT_SP_4 10	Set speed value
START_1	Start the heater
STOP_1	Stop the heater
START_4	Start the motor
STOP_4	Stop the motor
RESET	Switch to normal operating mode
SET_MODE_n (n=A, B or D)	Set operating mode

NAMUR Commands	Function
OUT_SP_12@n	Setting WD safety limit temperature with set value echo
OUT_SP_42@n	Setting WD safety limit speed with set value echo
OUT_WD1@m	<p>Watchdog mode 1: If event WD1 should occur, the heating and stirring functions are switched off and ER 2 is displayed. Set watchdog time to m (20 – 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.</p>
OUT_WD2@m	<p>Watchdog mode 2: if event WD2 should occur, the speed target value is changed to the WD safety speed limit and the temperature target value is changed to the WD safety temperature limit value. The warning WD is displayed. The WD2 event can be reset with OUT_WD2@0 – this also stops the watchdog function. Set watchdog time to m (20 – 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.</p>

### “Watchdog” functions, monitoring of the serial data flow

If, once this function has been activated (see NAMUR commands), there is no retransmission of the command from the computer within the set time (“watchdog time”), the heating and stirring functions are switched off in accordance with the set “watchdog” function or are changed to the set target values.

The data transmission may be interrupted by, for example, a crash in the operating system, a power failure in the PC or an issue with the connection table between the computer and the device.

#### “Watchdog” – mode 1


If there is an interruption in data communications (longer than the set watchdog time), the heating and stirring functions are switched off and ER 2 is displayed.

#### “Watchdog” – mode 2

If there is an interruption in data communications (longer than the set watchdog time), the speed target value is changed to the WD safety speed limit and the temperature target value is changed to the WD safety temperature limit value. The warning WD is displayed.


### PC 1.1 cable

This cable is used to connect the 9-pin plug to a computer.

### USB cable A – B

This cable is used to connect the USB interface to a PC.

# Maintenance and cleaning

## /// Proper handling

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

### Cleaning

Remove the device from the mains before cleaning.

Dirt	Cleaning agent
Dyes	Isopropyl alcohol
Building materials	Water containing detergent / isopropyl alcohol
Cosmetics	Water containing detergent / isopropyl alcohol
Food	Water containing detergent
Fuels	Water containing detergent
Other materials	Please consult IKA

- › Wear protective gloves during cleaning the devices.
- › Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- › Do not allow moisture to get into the device when cleaning.
- › Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

### Ordering spare parts

When ordering spare parts, please indicate:

- › device type.
- › serial number, see product label (bottom of product).
- › line item and description of spare part, see [www.ika.com](http://www.ika.com).
- › software version (upon device starting up).

### Repairs

Please only send in for repair devices that have been cleaned and are free of hazardous substances.

Please use our „**Decontamination Certificate**“ that is available on [www.ika.com](http://www.ika.com).

If your device requires repair, please return in its original packaging.

Contact IKA if you no longer have the original packaging.

Storage packaging is not sufficient when sending the product to IKA.

Please use appropriate packaging for transportation.

## Accessories

### /// An overview

#### Magnetic stir bar

- › ø 7 mm; length max. 60 mm
- › ø 10 mm; length max. 80 mm

#### Stirring bar

- › **RS 1** set of magnetic stir bars
- › **RSE** stir bar remover

#### Beaker

- › **H 1000** beaker 1 l
- › **H 1500** beaker 1.5 l
- › **H 3000** beaker 3 l
- › **H 5000** beaker 5 l
- › **H 8000** beaker 8 l

#### Square carrier

- › **H 135.10** square carrier without handle
- › **H 135.11** square carrier with handle

#### Blocks

- › **H 135.101** block 16 × 4 ml
- › **H 135.102** block 16 × 8 ml
- › **H 135.103** block 9 × 16 ml
- › **H 135.104** block 4 × 20 ml

- › **H 135.105** block 4 × 30 ml
- › **H 135.106** block 4 × 40 ml
- › **H 135.107** block 100 ml
- › **H 135.108** block 250 ml

#### Flask carrier

- › **H 135.20** flask carrier 100 ml without handle
- › **H 135.21** flask carrier 100 ml with handle
- › **H 135.25** flask carrier 250 ml without handle
- › **H 135.26** flask carrier 250 ml with handle
- › **H 135.30** flask carrier 500 ml without handle
- › **H 135.31** flask carrier 500 ml with handle
- › **H 135.40** flask carrier 1,000 ml without handle
- › **H 135.41** flask carrier 1,000 ml with handle
- › **H 135.50** flask carrier 2,000 ml without handle
- › **H 135.51** flask carrier 2,000 ml with handle

#### Flask inlay

- › **H 135.201** flask inlay 10 ml
- › **H 135.202** flask inlay 25 ml
- › **H 135.203** flask inlay 50 ml
- › **H 135.301** flask inlay 100 ml
- › **H 135.302** flask inlay 250 ml
- › **H 135.401** flask inlay 500 ml
- › **H 135.501** flask inlay 1,000 ml

#### General accessories

- › **H 16 V** support rod
- › **H 16.3** extension
- › **H 38** holding rod
- › **H 44** boss head clamp
- › **PT1000** series temperature sensor
- › **ETS-D5** electronic contact thermometer
- › **ETS-D6** electronic contact thermometer

For further accessories see [www.ika.com](http://www.ika.com).

## Error codes

/// Troubleshooting correctly

The stirrer displays an error code; following is an overview of possible errors and solutions.

#### Proceed as follows in such cases:

- › Turn off the device by using the power button.
- › Carry out corrective measures.
- › Restart the device.

#### Error code | Causes | Effect | Solutions

##### Err. 1 – Watchdog error 1

Causes	<ul style="list-style-type: none"> <li>› PC does not transmit any data within the set watchdog time</li> <li>› connection to PC interrupted</li> </ul>
Effect	<ul style="list-style-type: none"> <li>› heating switched off</li> <li>› motor switched off</li> </ul>
Solutions	<ul style="list-style-type: none"> <li>› change watchdog time</li> <li>› transmit data from PC within set watchdog time (OUT_WDx@m)</li> <li>› check cable and plug</li> </ul>

##### Err. 2 – Watchdog error 2

Causes	<ul style="list-style-type: none"> <li>› PC does not transmit any data within the set watchdog time</li> <li>› connection to PC interrupted</li> </ul>
Effect	<ul style="list-style-type: none"> <li>› heating switched off</li> <li>› motor switched off</li> </ul>
Solutions	<ul style="list-style-type: none"> <li>› change watchdog time</li> <li>› transmit data from PC within set watchdog time (OUT_WDx@m)</li> <li>› check cable and plug</li> </ul>

---

**Err. 3 – Temperature inside unit is higher than +80 °C**

---

Causes	<ul style="list-style-type: none"><li>› heat accumulation between heating plate and housing</li><li>› permitted ambient temperature exceeded</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li><li>› change experiment</li><li>› observe maximum permissible ambient temperature</li></ul>

---

**Err. 4 – Motor control is unavailable**

---

Causes	<ul style="list-style-type: none"><li>› motor blocked or overloaded</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li><li>› motor switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› reduce load torque or use smaller magnetic rods</li><li>› reduce target speed</li></ul>

---

**Err. 5 – No temperature increase measured by temperature sensor (selected time in menu)**

---

Causes	<ul style="list-style-type: none"><li>› sensor not in medium</li><li>› volume of medium to be measured too large</li><li>› heat conductivity of medium to be measured too low</li><li>› heat conductivity of the vessel is too low</li><li>› in the case of indirect heating, the overall heat conductivity resistance is too large</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› place the sensor in the medium</li><li>› reduce the volume of the media</li><li>› use a carrier fluid with better heat conductivity properties</li><li>› replace the glass vessel with a metal pot</li><li>› increase the “Time-out” period</li></ul>

---

**Err. 6 – Sensor plug out during Heating**

---

Causes	<ul style="list-style-type: none"><li>› defective connection cable</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› replace cable</li></ul>

---

---

**Err. 13 – Hotplate safety sensor, open-circuit**

---

Causes	<ul style="list-style-type: none"><li>› target / current different of the adjustable safety circuit for minimum temperature monitoring</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› after switching on, change the SAFE TEMP to a different value; if this solves the issue, the previous value can be reset by switching the device off and on again</li></ul>

---

**Err. 14 – External temperature sensor, short-circuit**

---

Causes	<ul style="list-style-type: none"><li>› short circuit in temperature sensor plug</li><li>› short circuit in the cable or temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› check the plug</li><li>› replace the temperature sensor</li></ul>

---

**Err. 21 – Fault during heating plate safety test**

---

Causes	<ul style="list-style-type: none"><li>› safety relay does not open</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li></ul>

---

**Err. 22 – Fault during heating plate safety test**

---

Causes	<ul style="list-style-type: none"><li>› S_CHECK cannot generate H_S_TEMP</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li></ul>

---

**Err. 24 – Heating plate temperature is higher than the set safety temperature**

---

Causes	<ul style="list-style-type: none"><li>› SAFE TEMP H (Hotplate) has been set to lower than the current temperature of the heating plate</li><li>› disconnection of heating plate control temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› leave the heating plate to cool</li><li>› set the SAFE TEMP H (Hotplate) higher</li></ul>

---

---

**Err. 25 – Heating switching element monitoring fail**

---

Causes	<ul style="list-style-type: none"><li>› heater control circuit switch (TRIAC) short-circuited</li><li>› safety relay has interrupted the heating circuit</li><li>› heater or the supply line is disconnected</li><li>› disconnection of heating plate safety temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li></ul>

---

**Err. 26 – Plate temperature > plate safety temperature (more than +40 K)**

---

Causes	<ul style="list-style-type: none"><li>› irregular temperature distribution across heating plate due to sporadic heat dissipation</li><li>› defective control or safety temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li><li>› ensure regular heat dissipation when using metal blocks, etc. on the flat surface of the heating plate</li></ul>

---

**Err. 31 – Fault in the heater switch element**

---

Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› contact customer service</li></ul>

---

**Err. 44 – Heating plate safety temperature is higher than the set safety temperature**

---

Causes	<ul style="list-style-type: none"><li>› SAFE TEMP H (Hotplate) has been set to lower than the safety temperature of the heating plate</li><li>› disconnection of heating plate safety temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› leave the heating plate to cool</li><li>› set the SAFE TEMP H (Hotplate) higher</li></ul>

---

**Err. 46 – Plate safety temperature > plate temperature (more than +40 K)**

---

Causes	<ul style="list-style-type: none"><li>› irregular temperature distribution across heating plate due to sporadic heat dissipation</li><li>› defective control or safety temperature sensor</li></ul>
Effect	<ul style="list-style-type: none"><li>› heating switched off</li></ul>
Solutions	<ul style="list-style-type: none"><li>› switch device off, leave to cool and switch on again</li><li>› ensure regular heat dissipation when using metal blocks, etc. on the flat surface of the heating plate</li></ul>

**If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:**

- › contact the service department.
- › send the device for repair, including a short description of the fault.

## Warranty

### /// Regulations

In accordance with IKA warranty conditions, the warranty period is 24 months. Upon registering this device at [www.ika.com/register](http://www.ika.com/register), your product will be subject to a lifetime warranty. For claims under warranty, please contact your local dealer.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

# Technical data

## /// IKA Plate (RCT digital) in detail

### Technical data IKA Plate (RCT digital)

Number of stirring positions	1
Stirring quantity max. per stirring position (H <sub>2</sub> O)	20 l
Motor rating output	9 W
Direction of rotation	right / left
Speed display set-value	LCD
Speed display actual-value	LCD
Speed control	Turning knob
Speed range	0 / 50 – 1,500 rpm
Setting accuracy speed	10 rpm
Speed variation (no load, nominal voltage, at 1500 rpm, ambient temperature +25 °C)	± 2 %
Stirring bar length	30 – 80 mm
Self-heating of the hotplate by max. stirring (RT: +22 °C / duration:1 h)	+8 °C
Heat output	600 W
Temperature display set-value	LCD
Temperature display actual-value	LCD
Temperature unit	°C / °F
Heating temperature	RT + device self-heating to +310 °C
Heat control	Turning knob
Temperature setting range	0 – 310 °C
Heating rate (1 l water in H 1,500)	7 K / min
Temperature setting resolution of heating plate	1 K
Connection for ext. temperature sensor	PT1000 serie, ETS-D5, ETS-D6
Temperature control accuracy (500 ml water in 600 ml glass beaker, 40 mm bar, 600 rpm, +50 °C)	± 0.5 K ( with external temperature sensor PT1000) ± 0.5 K (with ETS-D5) ± 0.2 K (with ETS-D6)
Temperature control accuracy of heating plate (at +100 °C)	± 5 K

Temperature setting resolution of medium	1 K
Adjustable safety circuit	+50 °C to +370 °C (± K)
Set-up plate material	Aluminium
Set-up plate dimensions	ø 135 mm
Automatic reverse rotation	yes
Intermittent mode	yes
Viscosity trend measurement	yes
Break detection stirring bar	no
Timer	yes
pH measurement	no
Graph function	no
Programs	no
Sensor in medium detection (Error 5)	yes
Weighing function	no
Set-up plate illuminated	no
Dimensions (W × H × D)	160 × 85 × 270 mm
Weight	2.4 kg
Permissible ambient temperature	+5 to +40 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
USB / RS 232 interface	yes
Analog output	no
Voltage	220 – 230 ± 10 % / 115 ± 10 % / 100 V ± 10 %
Frequency	50 / 60 Hz
Power input	650 W
Power input standby	1.6 W

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IKA India Private Limited
IKA England LTD.