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Instructions for use

Lambda Hind III / phiX Hae III Marker

This marker was manufactured by using unmethylated Lambda-DNA and restriction endonuclease *HindIII*, as well as phiX-174 DNA and *HaeIII*.

Versatile marker for depiction of all standard DNA fragment lengths. By mixing high molecular weight Lambda (*HindIII*) fragments with short-sequence phiX-174 (*HaeIII*) fragments a marker was designed which covers a broad range from 200 bp up to 20 kb. Most fragments are distributed in an equimolar way, additionally allowing rough DNA quantification in the gel. Due to this equimolarity, however, the smallest bands of 125 bp, 118 bp and 72 bp can not be detected efficiently in standard gels. Optimal results are obtained with a running distance of approx. 10 cm.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

The following fragment sizes (in bp) are included: 23130 1), 9416, 6557, 4361 1), 2322, 2027, 1353, 1078, 872, 603, 564, 310, 281, 271, 234, 194, 118, 125, 72.

¹ Labeled fragments contain cohesive ends of bacteriophage Lambda which often hybridise, resulting in fragments/bands of higher molecular weight. We recommend heating the marker prior to use as described below, in order to separate those two fragments.

The whole Lambda phage genome contains 48.502 bps.

Content:

CP49.1 - 1x 50 µg Lambda *HindIII* / phiX *HaeIII* Marker
0100.1 - 1x sample buffer ROTI®Load DNA (with glycerol)

Application:

Depending on the application, the marker can be dissolved directly in the supplied sterile-filtered 1 x sample buffer ROTI®Load DNA (with glycerol) (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1).

The lyophilized DNA is dissolved for 15 minutes at room temperature in an appropriate volume (see below) under occasional stirring. Prior to use we recommend to dissolve the cohesive ends by incubating at 65 °C for 5 mins, followed by cooling at 4 °C or on ice.

Sample application/concentration:

Standard loading for mini to midi gels per trace is:

- with bands visible in UV-light after ethidium bromide-staining:
0.5 – 1.0 µg in 10 µl = 1 lane
- with detection after ethidium bromide-staining with signal-integrated camera systems:
0.3– 0.4 µg in 10 µl = 1 lane

Storage:

Optimal storage temperature is –20 °C. Please avoid repeated (>20 times) thawing and freezing. We recommend portioning in aliquots.


Fig.: 1 % agarose (NEEO, Art. No. 2267).

Lambda Hind III / phiX Hae III Marker
CP49.1 50 µg + sample buffer