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PAS (Periodic Acid Schiffs) - Färbekit

HP01

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HP01

Detection of Glycoconjugates

PAS staining is an approved method to reveal glycoconjugates in tissue probes, e.g. neutral mucopolysaccharides, mucoproteins, glycoproteins, glycogen.

Specific characteristic: PAS staining visualises no cell types but certain metabolic situations in cells.

Principle:

Periodic acid oxidises free 1,2-glycol groups to aldehyde groups, which subsequently are stained red by fuchsin sulfurous acid (Schiff's reagent). Hemalum solution causes counterstaining of cell nuclei.

Kit contains:

- Schiff's reagent (X900.1) **Danger** H290-H350 500 ml
- Periodic acid solution 1 % (HP00.1) **Warning** H315-H319-H373 500 ml
- Hemalum solution acidic acc. to Mayer (T865.1) **Warning** H302-H315-H319 500 ml

The staining solutions should be filtrated before use! Solutions may be bought separately.

*Note: Store **Schiff's reagent** well-sealed after opening. Discard solution if it turns reddish!*

Additional chemicals required:

- Clearing Agents: ROTI®Histol (Art. No. 6640), ROTICLEAR® (Art. No. A538) or Xylene p.a. (Art. No. 4436)
- Mounting Media: ROTI®Histokitt (Art. No. 6638), compatible with ROTI®Histol
 ROTI®Mount (Art. No. HP68), compatible with ROTICLEAR®
 ROTI®Histokitt II (Art. No. T160), compatible with Xylene

Instructions:

1. De-wax sections and rehydrate by descending alcohol series	7. Rinse with warm tap water (35 °C minimum) approx. 5 min
2. Rinse quickly with distilled water	8. Rinse quickly with distilled water
3. Hydrolyse with periodic acid solution 1 % 10 min	9. Counterstain with hemalum solution acc. to Mayer 5 min
4. Rinse with tap water 10 min	10. Blue in flowing tap water 10-15 min
5. Rinse with distilled water 2x2 min	11. Dehydrate by ascending alcohol series
6. Stain with Schiff's reagent (room temperature) 10-20 min	12. Clear with clearing agent Mount with appropriate mounting medium

For cytological preparations:

After hydrolysis with periodic acid rinse 10 min with distilled water, counterstain with hemalum solution 3-5 min, blue in running tap water 5-10 min.

Results:

PAS-positive cells: pink to violet cell nuclei: blue

Note: The colour intensity depends on the pre-treatment and the composition of the samples to be stained. It may initially be necessary to adapt the method to the respective conditions.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com

JWu 08/2023

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

PAS Staining Kit

HP01.1