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ROTI®Load DNA orange 2

(with glycerol)

6 x gel loading buffer for DNA-electrophoresis

I. Introduction:

ROTI®Load DNA orange 2 is a 6x concentrated gel loading buffer for DNA electrophoresis. Especially convenient for separation of a large range of DNA-fragments in the gel. Also well suited in order to avoid concealing of important DNA-bands by bromophenol blue. Each batch is tested for functionality in electrophoresis and for being DNase-free in order to guarantee the highest safety and reproducible results.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com

ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

II. Composition:

ROTI®Load DNA orange 2 (with glycerol) is composed as follows:

TRIS/HCl pH 7,5	60 mM
Na-acetate	30 mM
EDTA	12 mM
Glycerol	60 % (w/v)
Orange G	0.36 % (w/v)
Xylene cyanol blue	0.12 % (w/v)

The density of your DNA-solution is increased by the proportion of glycerol share in ROTI®Load which guarantees even distribution in the gel pocket base. Xylene cyanol and orange G alleviate sample application and can be used as a run marker. By adding Tris and EDTA, DNA stability is increased.

III. Application:

ROTI®Load DNA orange 2 can either be added to DNA-solution or the DNA can be dissolved directly in ROTI®Load DNA orange 2.

Dissolving DNA in 1 x ROTI®Load DNA orange 2:
Dilute ROTI®Load DNA orange 2 1:6 by pipetting 170 µl ROTI®Load DNA orange 2 into a clean, sterile test tube. Add 830 µl bidistilled, autoclaved water and mix by vortexing briefly. This solution (ROTI®Load DNA orange 2 1x) can now be used to dissolve precipitated DNA and can be applied directly onto the gel.

Adding ROTI®Load to a DNA-solution:

Volume DNA-solution	Volume ROTI®Load	Total volume Sample
5 µl	+	1 µl = 6 µl
10 µl	+	2 µl = 12 µl
15 µl	+	3 µl = 18 µl
20 µl	+	4 µl = 24 µl
25 µl	+	5 µl = 30 µl

For example, add 1 µl ROTI®Load DNA orange 2 to a DNA solution with a volume of 5 µl. Mix the sample by carefully pipetting up and down three times.

OPTIONAL:

Irrespective of which method you use, we recommend denaturing the sample for 5 mins at 65 °C and then cooling on ice, otherwise run artefacts may develop (e.g. EcoRI/HindIII-separated Lambda-DNA).

IV. Storage:

Buffer should be stored at -20 °C for long-term storage (> 3 months). Buffer currently in use can be stored for some weeks at 4 °C (or room temperature).

V. Quality assurance:

Every ROTI®Load DNA orange 2 batch is tested for being DNase-free and for functionality in electrophoresis.

For this purpose, 1 µg linearized plasmid-DNA is incubated for 48 h at 37 °C in 1x DNA-buffer and the integrity of the DNA electrophoretically compared with a fresh sample.

ROTI®Load DNA orange 2

(with glycerol) HP05.1 5x1.8 ml