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Pen/Strep-PreMix

CELLPURE®

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Pen/Strep PreMix

CELLPURE®

Antibiotic mixture for elimination and prevention of bacterial contamination in cell culture

I. Introduction:

Antibiotic mixture made of penicillin, streptomycin and stabilising sodium chloride for application in cell culture. The premixed powder is extremely stable during storage at room temperature and, due to the excellent solubility of the reagents, a 100 x stock solution can be easily prepared.

II. Content:

Roth Pen/Strep PreMix is composed as follows (amount per vial):

Penicillin G Na-salt	312.5 mg
Streptomycin	500 mg
NaCl	450 mg

III. Mechanism:

Penicillin inhibits cell wall biosynthesis in Gram-positive bacteria. Streptomycin produces the antibiotic effect in Gram-negative bacteria through inhibition of translation.

Cytotoxic concentration is

Penicillin: 10.000 U / ml

Streptomycin: 20 mg / ml

IV. Application:

The powder mixture of one vial may be solubilised with 49.2 ml water to gain exactly 50 ml 100 x stock solution. The powder solubilises rapidly and, after tilting the vial several times, results in a homogenous solution. We recommend to use pyrogen-free water for solubilisation (e.g. Art. No. 3255.1). The dilution should be sterilised through 0.2 µm filters and frozen in aliquots of 1 - 10 ml.

Concentrations in the 100 x stock solution:

Penicillin G Na-salt 10.000 I.U / ml

Streptomycin 10 mg / ml

NaCl 0.9 % (154 mM)

Amount of endotoxins is ≤210 I.U. per vial and ≤4.2 I.U./ml in the 100 x stock solution when pyrogen-free water is used. After solubilisation, pH value is approx. 5.5.

For cell culture, we recommend to dilute the stock solution by 1:100 in cell culture medium. Antibiotics solution should be added to the cell culture medium prior to its application to the cells.

Dilution by 1:100 results in the following end concentrations:

Penicillin G Na-salt 100 I.U / ml

Streptomycin 100 µg / ml

NaCl 0.009 % (1.54 mM)

The cell culture medium should be changed after 3 days at most and fresh medium including antibiotics should be added to the cells. Stock solution (2 x 50 ml) is sufficient for 10 l medium.

IV. Storage and Stability:

Powder: store at +4°C.

Stock solution: store in aliquots at -20 °C.

Stock solution should be frozen and thawed only once.

Antibiotics in solution are stable for approx. 1 year at -20 °C and for approx. 3 days at 37 °C.

Danger H302-H317-H334-H361d P260-P270-P302+P352-P309+P311

Pen/Strep-PreMix

HP10.1 2 vials for 50 ml stock solution each

Carl Roth GmbH + Co. KG

Schoenperlenstraße 3-5 • 76185 Karlsruhe

P.O. Box 100121 • 76231 Karlsruhe

Phone: +49 (0) 721/ 5606-0

Fax: +49 (0) 721/ 5606-149

info@carlroth.com • www.carlroth.com

ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.