ROTI®PinkDNA

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Instructions for use

ROTI®PinkDNA

Synthetic co-precipitant for efficient precipitation of nucleic acids. Ready-to-use.

- **With coloured dye for easily visible pellets.**
- Simply add to nucleic acid precipitation. Suitable for precipitation of all nucleic acids.
- Recovery rates of up to 99 %.
- Recovered DNA may be immediately applied to all standard applications.
- Not compatible with photometric concentration measurement at 260-280 nm.

Precipitation additive for recovery of DNA from all solutions, even from solutions with very low concentrations of nucleic acids. ROTI®PinkDNA is produced synthetically and is free of nucleic acids. DNA or dsRNA precipitated with ROTI®PinkDNA is highly pure and may be applied directly to all standard applications like transformation, restriction digest, cloning, sequencing, *in vitro* transkription etc.

Suitable for all nucleic acid fragments, most efficient for fragments ≥ 25 bp.

Content:


Linear Polyacrylamide: 5 mg/ml
Pink coloured dye

Application - Ethanolic Precipitation

[Examples given for 500 μ l DNA solution]

1. Add 1/10 volume sodium acetate (3 M, pH 5.2) to DNA solution [50 μ l]
2. Vortex
3. Add 6 μ l (30 μ g) co-precipitant ROTI®PinkDNA
4. Vortex
5. Add 2 volumina ethanol_{abs.} (min. 1.5 volumina, max. 2.5 volumina)* [1 ml]
6. Incubate at roomtemperature for 10 min.
7. Centrifuge at min. 15.000 g (13.000 rpm in Mini-Centrifuge) for 10 min.
8. Carefully remove supernatant
9. Wash pellet with ca. 1 vol. 70 % ethanol [500 μ l]
10. Centrifuge at min. 15.000 g (13.000 rpm in Mini-Centrifuge) for 10 min.
11. Carefully remove supernatant
12. Shortly centrifuge pellet, in order to collect residual ethanol in the tip of the reaction tube (ca. 10 sec.)
13. Carefully remove residual supernatant
14. Dry pellet at roomtemperature (ca. 3 min.)
15. Resuspend in suitable buffer (sterile water, 10 mM Tris-buffer, TE-buffer)

* Alternatively, 1 volumen isopropanol instead of ethanol may be used. Please note: the pellet does then not adhere firmly to the side of the reaction tube.

Please use additionally:

Sodium acetate 3 M in water, pH 5.2
(Roth Art. No. 6779)

Ethanol_{abs.} (Roth Art. No. 9065)

Ethanol 70 % (Roth Art. No. T868)

Resuspension buffer:

sterile water *or*

10 mM Tris pH 8.0 *or*

1 x TE-buffer (ROTI®Stock 100 x TE-buffer,
Roth Art. No. 1052)

Storage:

Standard storage for frequent use should be done at +4 °C.

For long-term storage we recommend -20 °C.

ROTI®PinkDNA

HP54.1	0.5 ml
HP54.2	1.5 ml
HP54.3	3.0 ml

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.