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HP69

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ROTI® Nucleic Acid-free

HP69

Ready-to-use solution for efficient removal of nucleic acid contaminations from surfaces

I. Introduction

ROTI® Nucleic Acid-free was specifically designed for extensive removal of nucleic acid contaminations. Being a mild reagent, ROTI® Nucleic Acid-free may be used on tables, laboratory devices, as well as on glass and plastic wares. By incubation or by wiping, laboratory benches, heating blocks, thermocyclers, pipettor shafts, reaction tubes a.s.o. can easily be cleaned. Contamination with foreign DNA or RNA, which may heavily interfere with sensitive assays like PCR, reverse transcription or sequencing, is hence efficiently eliminated.

II. Application

Wiping (table surface, thermal cycler etc.)

1. Apply undiluted ROTI® Nucleic Acid-free to the surface that is to be cleaned.
2. Incubate for approx. 1 min.
3. Wipe the surface covered with ROTI® Nucleic Acid-free with clean wipes and rub down.
4. Sensitive surfaces: wash with clear water free of nucleic acid and rub down.

Incubation (reaction tubes, glass ware, pipettor shafts etc.)

4. Place items for cleaning into undiluted ROTI® Nucleic Acid-free. All items should be completely covered by the solution.
5. Incubate at room temperature over night.
6. Rinse with nucleic acid-free water and rub down or let dry. Items made of alcohol-resistant material may additionally be rinsed with 70 % ethanol and left to air dry.

III. Additionally Recommended Reagents

Nucleic acid free water:

Water ROTIPURAN®, highly pure water for sample preparation in organic trace analysis. Optical absorption at 254 nm: ≤ 0.005 Art. No. HN57
Water, sterile, pyrogen-free, hypotonic, acc. to regs. of the Ph. Eur. for *Aqua purificata* Art. No. 3255
Water, DEPC-treated, sterile and autoclaved, for molecular biology Art. No. T143

Paper wipes:

ROTIZELL®-Tissue-Wipes Art. No. 0087

 **Warning** H315-H319
P280-P302+P352-P305+P351+P338

ROTI® Nucleic Acid-free

Spray (500 ml)	HP69.1
1 Litre replacement solution	HP69.2
2.5 Litre replacement solution	HP69.3

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.