

# Instructions for use



## Gel leveling table

**N854.1**

**WARNING:**

**Please read the entire operator's manual thoroughly before operating this unit.**

Our gel leveling table comply with the statutory CE safety directives:

Our gel leveling table is designed to give long service and reproducible results in your laboratory. Please take a few minutes to read these instructions.

First check that the apparatus has been received complete and undamaged following shipment. Refer to the packing list and check that all components and accessories are present. Any faults or losses must be notified to us immediately.

**PLEASE RETAIN ALL PACKAGING MATERIALS UNTIL THE WARRANTY PERIOD HAS EXPIRED.**

### 1. PACKING LIST

No. Items	Description
1	30 x 20 cm Leveling Table
1	Operating Manual

## 2. SPECIFICATIONS

### A. Construction

- White Acrylic Construction.
- All acrylic joints chemically bonded.
- Easy to use Leveling Feet.
- Level Bubble.

### B. Environmental Conditions

- This apparatus is intended for indoor use only.
- This apparatus can be operated safely at an altitude of 2,000 m.
- The normal operating temperature range is between 4 °C and 65 °C.
- Maximum relative humidity 80 % for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C.

## 3. USING THE GEL LEVELING TABLE

### A. Safety Precautions

READ the instructions before using the apparatus.

### B. Cleaning

Rinse the Gel Leveling Table with DISTILLED WATER ONLY.

#### **IMPORTANT:**

Acrylic plastic is NOT resistant to aromatic or halogenated hydrocarbons, ketones, esters, alcohol's (over 25 %) or acids (over 25 %). NEVER USE ORGANIC SOLVENTS. Dry the gel levelling table then dry with tissues.

### C. Operating instructions

- Place the unit on a flat surface.
- Adjust the table by using the leveling feet until the bubble is centered in the centrally positioned bubble level.
- Once level, any electrophoresis unit gel casting apparatus can be place on the table for level gel pouring.

**Gel leveling table**

**N854.1**

#### **Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com

sse 06/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.