



ROTI®Lumin

Chemiluminescent substrate for Horseradish Peroxidase (HRP)-labelled reporter molecules

I. Description and Applications

ROTI®Lumin is a luminol-based chemiluminescent substrate designed for use with peroxidase-labelled (HRP) reporter molecules.

ROTI®Lumin provides good sensitivity over chromogenic substrates in both blotting and microwell assays.

Positive reaction sites are rapidly detected with high sensitivity and minimal background. The use of ROTI®Lumin allows for multiple stripping and reprobing of blots. In blotting applications, permanent results are recorded on X-ray film. In microwell assays, positive reactions are rapidly detected and read in a luminometer.

ROTI®Lumin can be used in applications such as ELISA, Western blotting, Southern blotting, dot blotting und colony hybridizations.

II. Principle

In the presence of hydrogen peroxide, HRP converts luminol to an excited intermediate dianion. This dianion emits light on return to its ground state. After reaction with HRP-conjugate, the light emission from ROTI®Lumin reaches maximum intensity within 5 minutes and is sustained for approximately 1-2 hours.

III. Membranes

ROTI®Lumin can be used with nitrocellulose (NC), nylon, and PVDF membranes.

ROTI®PVDF (T830.1), ROTI®Nylon plus (K058.1) and ROTI®NC membranes (HP40.1) are recommended for use with ROTI®Lumin.

IV. Blocking reagents

Milk or casein-based blocking solutions are recommended for use with ROTI®Lumin. BSA or serum-based blocking agents may cause elevated background. For Western-Blotting applications we recommend ROTI®Block (A151.1), an optimised polymer-based blocking reagent.

V. Application

- Perform standard blotting and/or immunoassay procedures.
- Following incubation with HRP antibody or HRP streptavidin, perform at least three washes.
- Mix ROTI®Lumin 1 (P079.1) and ROTI®Lumin 2 (P080.1) in equal volumes. Warm this ROTI®Lumin working solution to room temperature before use. (Stock solutions should be put back to cooling at once.) ROTI®Lumin working solution is stable for up to one hour at room temperature or up to 24 hours when stored at 2 °C - 8 °C. Working solution needs not be protected from light.

ELISA Procedure

Add 100 µl / well ROTI®Lumin working solution.

Read in a luminometer with 1 second integration time per well. ROTI®Lumin provides consistent results when read 5-45 minutes after addition of substrate.

Western Blotting Procedure

Use 1 ml substrate / 10 cm² membrane

- Immerse membrane in ROTI®Lumin working solution for one minute. This incubation may be performed in the light.
- Remove membrane from ROTI®Lumin working solution, touch the edge of the membrane to a piece of filter paper to remove excess substrate.
- Place membrane between two sheets of plastic foil.
- Expose membrane to X-ray film 1-10 minutes. Light emission begins immediately upon removal from substrate and continues for approximately 2 hours.
- Develop film.

VI. Repeated use of a blot: Removing bound antibodies and probes

For Western-Blots:

- Wash membrane in TBST-buffer
- Incubate in ROTI®Free Stripping Buffer for 30 minutes at 56 °C (e.g. in a water bath with shaking device) in fume hood.
- Wash twice in TBST for 20 minutes each.
- Blocking of the membrane (see IV.)
- Incubation in primary antibody and further detection according to protocol.

Efficiency of stripping may be controlled with brief chemoluminescence staining after washing in TBST. A slight reduction of signals of membrane-bound proteins and an increase in background may occur after repeated stripping.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

Please note:

By using ROTI®Block the saturation of nonspecific binding sites will remain stable. No additional blocking is required.

For Southern-Blots:

- Rinse membrane for 5 minutes in Washing buffer.
- Cover membrane with stripping buffer and incubate under gentle shaking at 55 °C for 20 min.
- Rinse membrane twice for 5 minutes each in 2x SSC.
- Prehybridize membrane before reprobing.

VII. Additionally Required Reagents

Washing Buffer: (example)

100 mM Tris, pH 7.5; 150 mM NaCl, 0.5 % Tween 20.

Stripping Buffer:

For Western-Blots: ROTI®Free Stripping Buffer (Art. No. 0083.1, 3337, 3319)

For Southern-Blots: (prewarm to 55 °C) 0.2 N NaOH, 0.1 % SDS

SSC: ROTI®Stock 20x SSC (Art. No. 1054.1)

TBST: ROTI®Stock 10x TBST (Art. No. 1061.1)

PBST: ROTI®Stock 10x PBST (Art. No. 1059.1)

Blocking: ROTI®Block (Art. No. A151.1)

Casein (Art. No. 8569.1)

Powdered milk (Art. No. T145.1)

VIII. Trouble Shooting

Blotting

Excess signal or high background	<ul style="list-style-type: none">• Decrease film exposure time• Decrease HRP conjugate concentration• Increase blocking times• Increase washing times• Load less protein/DNA onto gel
No signal	<ul style="list-style-type: none">• Verify transfer by staining protein gel with Coomassie blue or DNA gel with ethidium bromide• Verify protein transfer by staining membrane with Ponceau S or Amido black• Make sure HRP secondary antibody is specific for the primary antibody• “Blitzed” signal. Reduce concentration of primary and/or secondary antibody
Weak signal	<ul style="list-style-type: none">• Increase film exposure time.• Place the membrane on top of a Blotting paper wetted with ROTI®Lumin during exposition• Increase conjugate concentration• Increase conjugate incubation time• Load more protein/DNA onto gel

ELISA

Excess signal or high background	<ul style="list-style-type: none">• Decrease HRP conjugate concentration• Reduce conjugate incubation times• Increase blocking times• Increase washing times• Decrease the amount of protein coated to plate
No signal	<ul style="list-style-type: none">• Verify the luminometer is working correctly• Make sure HRP secondary antibody is specific for the primary antibody

Weak signal	<ul style="list-style-type: none">• Increase conjugate incubation time• Increase conjugate concentration• Increase the amount of protein coated to the plate
-------------	--

IX. Storage

Always store ROTI®Lumin stock solutions at approx. 4 °C (2-8 °C). Protect from light. Mixed Working solution may be stored at room temperature for approx. 1 h or at 4 °C for approx 24 hs. Working solution needs not be protected from light.

X. Shelf life

ROTI®Lumin stock solutions expire 1 year after production when stored according to instructions.

XI. Content

One MINI-Kit contains:

25 ml ROTI®Lumin sol. 1 (P079)
25 ml ROTI®Lumin sol. 2 (P080)
Sufficient for 500 cm² membrane.

One Kit contains:

120 ml ROTI®Lumin sol. 1 (P079)
120 ml ROTI®Lumin sol. 2 (P080)
Sufficient for 2.400 cm² membrane.

Contents of this Kit may not be bought separately.

ROTI®Lumin	P078.2	1 Mini-Kit
	P078.1	1 Kit