Kryoetiketten

Cryo Labels

Manual | Anleitung

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Right handling of the Cryo Labels


Temperature Resistance

The cryogenic labels withstand temperatures from minus 196 °C to plus 121 °C and are suitable for long-term storage in liquid nitrogen and transport on dry ice. Even after repeated freezing and thawing, the cryogenic labels have shown good adhesion. If a label partly peeling off, it is possible to fix it again by pressing. The cryo-labels can be attached to a frozen tube. You should nevertheless dry the surface thoroughly first.

Applying the Labels

Make sure that the surface of your vessel or tube is not damp before applying the label. To prevent the glue from drying out, apply the label at room temperature immediately after removing it from the backing film. For getting the best result, labels should be fixed at least 30 minutes before application at room temperature (approx. 20°C). If possible, stick the labels slightly overlapping.

Print Information

Please make sure that the label sheets do not slip when printing, to avoid any damage. Use the enclosed adjustment template to align your printer.

Chemical Resistance

Materials and Methods

Cryo Labels, format 28×12 mm/ white, marked with:

- Laser printing (HP Color Laser Jet CM 1312),
- Ball-point pen (no-name product),
- Pencil (Rheita, HB)

Labels were placed on 2 ml tubes and submersed in the respective solvent for 24 h. Shorter exposition times may increase resistance.

Overview Chemical Resistance at Room Temperature

Solvents	Adhesion	Laser printing		Ball-point pen		Pencil	
		Legibility	abrasion resistance	Legibility	abrasion resistance	Legibility	abrasion resistance
Acetic acid 25 %	G	G	G	G	G	G	L
Acetone	N	L	L	L	L	G	L
Acetonitrile	L	G	G	L	G	G	L
Ammoniumhydroxide (1M)	G	G	G	G	G	G	L
Ammoniumhydroxide (25%)	G	G	G	L	G	G	L
Ethanol	G	G	G	G	G	G	L
Ethanol 70 %	G	G	G	G	G	G	L
Formic acid (25%)	G	G	G	L	G	G	L
Hexane	N	G	G	G	G	G	L
Hydrochloric acid (10%)	G	G	G	G	G	G	L
Hydrogenperoxide (30%)	G	G	G	G	G	G	L
1-Propanol	L	G	G	G	G	G	L
Methanol	G	G	G	G	G	G	L
Methylenchloride	N	L	L	N	N	G	L
Phosphoric acid (25%)	G	G	G	G	G	G	L
Potassium Hydroxide (KOH, 1M)	L	G	G	G	G	G	L
Potassium Hydroxide (KOH, 32%)	L	L	L	L	L	L	L
Sodium hydroxide (1M)	G	G	G	G	G	G	L
Tetrahydrofurane (THF)	N	N	N	G	G	G	L
Toluene	N	L	N	G	G	G	L
Water distilled deionized	G	G	G	G	G	G	L

G	Good chemical resistance
L	Limited chemical resistance
N	No chemical resistance, use not recommended

Using the Adjustment Template to Set the Printer

To adjust the exact print alignment, we recommend one of the following procedures, depending on whether you use the labelling program Barcode Forge scienova or the Microsoft Word templates available in our Onlineshop.

1. Use of the calibration template with the labelling software Barcode Forge scienova

You can purchase the program Barcode Forge scienova in the scienova shop or download a free test version (runtime 30 days).

- Start the program and select the format corresponding to your label sheets under "Label" → "Format"
- Under "Label" → "Etikett", add a frame line with a thickness of 0.1 mm or more with a margin of 0 mm. Select a corner radius that corresponds your label format.
- Now go to "Drucken" → "Projekt drucken" and select your printer. In the printer settings under properties, select the desired paper tray.
- Advice: We recommend that you adjust your default print settings in the Microsoft Windows Control Panel. So you do not have to make the print settings again before print operation.
- Click "Weiter" and to print the label under [...nicht mit Daten befüllen, sondern...], click on the first one (top left).
- Insert the supplied calibration template into the printer and print the labels. Attention: The positioning aids in the corresponding compartment must be set as precisely as possible to A4 format.
- If the printed frame differs from the frame line of the template, measure the horizontal and vertical displacement.
- Under "Drucken" → "Projekt drucken" → "Drucker auswählen" you can configure a corresponding adjustment (in the tenth-millimeter range).
- Print the next label on the sheet. Continue to pay attention to the accuracy of the insertion.
- Repeat the process until you see no or only a slight deviation. To check accuracy, print one of the lower labels.

Once you have made the appropriate settings for your printer, you can now create the label and print it reliably. Place the elements on the labels at least 1-1.5 mm from the edge.

2. Use of the Calibration Template with Microsoft Word template

You can download the Microsoft Word templates for the labels on the product pages of the scienova webshop.

- First open the template document in Microsoft Word.
- Enter a sample text and/ or illustration in one of the fields. Place the elements on the labels at least 1-1.5 mm from the edge.
- Advice: You can show the grid lines in the table properties . Click on the Table Layout tab and then under Settings on "show Gridlines". This will simplify the orientation and the optimal alignment of your printer to the template.
- Insert the supplied calibration template into the printer and print the labels. Attention: The positioning aids in the corresponding compartment must be set as precisely as possible to A4 format.
- If the printed frame differs from the frame line of the template, measure the horizontal and vertical displacement. This allows you to set a small offset in the settings menu of your printer (directly on the device).
- Repeat the process with the next table field until you have the best accuracy you can achieve.
- In Microsoft Word, you can move the template to the side and align it even more precisely: In the table tools under "Layout" → "Table Properties", you can easily move the entire table on the sheet. Test the different alignment variants [left, centered and right] and approach the optimal alignment.

scienova labelling software

If you print labels frequently and need to process a lot of information, our label printing software Barcode Forge scienova makes your work easier. You can easily export existing data sets, integrate graphics and generate a variety of barcode types.

We provide a free 30-day demo version, which expires after this period and becomes inactive. You can find the demo version of Barcode Forge scienova in our webshop www.scienova.com.

Once you have decided to purchase the software, you have a license without a time limit. You buy once and benefit from unlimited free updates.

Scope of functions Barcode Forge scienova

- Freely formatable text fields
- Field function for inserting date or consecutive numbers
- Inserting graphics or images
- Printing barcodes:
 - Code 128, Code 128 B, Code 128 C / 2 of 5 Interleaved / Code 39, Code 39 Extended (Full ASCII) / Code 11
 - Codabar (bzw. NW-7) / EAN-13 (GTIN-13), EAN-8 (GTIN-8), UPC-A / ISBN-13
 - Code 93, Code 93 Extended (Full ASCII) / PZN, PZN8
 - QR Code / DataMatrix / PDF 417 / Aztec
- Data import is supported:
 - via Microsoft Excel or Microsoft Access,
 - ODBC and SQL queries
 - via CSV

Please do not hesitate to contact us if you have any questions or suggestions. You can reach our support team Mo. - Fr. 9:00 - 16:00 o'clock under the e-mail address: info@scienova.com or by telephone at +49(0)3641 352 678.

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