HUPFER

Hupfer Metallwerke GmbH & Co. KG

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EN | Instruction manual

Instrument baskets, sieve trays and instrument fixations

Scope

This instruction manual is for the following medical products of Class 1 according to regulation (EU) 2017/745, manufactured by Hupfer Metallwerke GmbH & Co. KG.

Product category	Base UDI
Instrument basket	4044593SiebkoerbeNU
Sieve tray	4044593SiebschalenZE
Instrument fixation	4044593InstrumentenfixWY
Other designations	Fine mesh basket, adjustable fixing clamp, composite sieve tray, washing basket, suspended shelf, cover, ENT rack, instrument basket.

Safety information

Intended use

The product is intended for the transportation, treatment including sterilisation, storage and provision of medical instruments. The product is intended for reuse and treatment including sterilisation, by the operator.

The product may be used by medical professionals and laypersons after instruction by medical professionals. The maximum permitted payload is 15 kg according to ANSI/AAMI ST 77:2013. The maximum permitted temperature is 150° C.

Read and observe the operating instructions.

Observe the information in the treatment instructions.

The product is not sterile when delivered and must undergo a treatment cycle prior to use.

If the product has been used for patients with suspected Creutzfeldt-Jakob disease, it must be disposed of according to national regulations. In Germany, follow the instructions of the Robert Koch Institute.

The product is heated during cleaning, disinfection and sterilisation. There is a risk of burns on contact with the hot product.

- Wear heat-insulating protective gloves.
- Let the product cool down sufficiently.

The product may transport contaminated instruments. There is a risk of infection when coming into contact with contaminated instruments.

- Wear suitable protective gloves.
- Handle the product with care.

Damaged products have sharp edges. There is a risk of injury on sharp edges.

- Check the product regularly.
- Dispose of damaged products.

Loose instruments and inserts may fall out during transport and cause injury.

- Please observe the assembly diagram.
- Seal the product if possible.

Performance characteristics

All materials used are intended for treatment including sterilisation and withstand temperatures up to 150° C.

Stainless steel and silicone components are intended for treatment and reuse. If no damage is found during visual inspection, the product can be used after treatment.

Instructions for use

The product is intended for reuse and treatment including sterilisation, by the operator. Treatment includes inspecting and checking the product.

The following are indicators of damage to the product. Damaged products must not be used.

- Components made of stainless steel
 - Damaged overall structure, e. g. frame, perforated plates, sieve structure
 - Broken wire supports
 - Loose weld spots
 - Cracks
 - Corrosion
 - Components made of silicone
 - Sticky, swollen surface
 - Cracks
- Replace silicone components with the above damage to continue using the product.

Loading must be specified by the operator in a loading plan.

Store and transport the product dry and protected from sunlight.

Recommendation for use

Do not put the product in isotonic solution. Prolonged contact with isotonic solution will lead to corrosion. Only load the product with pre-cleaned instruments. Do not transport waste, residues of skin disinfectants, saline solution, etc. with the product.

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