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




Tootja

REF Referents = toote number

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## Operation Instructions: Disposable micro pipettes with ring cap 'ringcaps®'

1+2+3+4+5µL, 5+10µL, 10µL, 20µL, 25µL, 20+40µL, 40µL, 44,7µL, 50µL, 50+100µL, 100µL, 100+200µL  
(Sodium heparinised, without coating)

### Instructions for use

Capillaries for the transfer and quantitative analysis of defined volumes of capillary blood.

Suitable for medical professionals only

Application temperature: +2 – +40°C

National guidelines for occupational safety and quality assurance must be observed.

Accuracy [A]: 0.25 % (calculation in accordance with DIN ISO 7550)

Precision [P]: 0.5 % (calculation in accordance with DIN ISO 7550)

Accuracy [A]: 0.3 % (calculation in accordance with DIN ISO 7550)

(only for item no. 9600105)

Precision [P]: 0.6 % (calculation in accordance with DIN ISO 7550)

(only for item no. 9600105)

### Purpose

✓ **Intended purpose:** Disposable micro pipettes for transferring defined volumes of capillary whole blood for quantitative analysis. The blood is entered into the capillary directly at the analysis device.

✗ **Prohibitions of use:** Disposable micro pipettes are NOT intended for the following applications:

Sampling of other liquids than capillary whole blood, storage of samples, measurements within the capillary itself.

## Hazards and safety

The precautionary measures required for the use of laboratory reagents and body fluids as well as microbiological samples must be observed. Handling may only be carried out by qualified personnel. National and internal laboratory guidelines for occupational safety and infection control must be observed. Suitable protective clothing, protective goggles and disposable gloves must be worn during work.

Attention must be paid to effective protection against infection in accordance with the laboratory guidelines.

The user is committed to report any serious incidents in connection with the product to the competent authority of the Member State in which the user and/or the patient is registered and to Hirschmann Laborgeräte GmbH & Co. KG.

The capillaries comply with the requirements of Regulation (EU) 2017/746 regarding in vitro diagnostic medical devices.

## Disposal

The national statutory regulations must always be observed when disposing of the capillaries. Used capillaries contaminated with potentially infectious substances of human origin may pose a risk of injuries, infections or microbiological hazards. Potentially contaminated capillaries must be disposed of in such a way that injuries, infections or microbiological hazards to third parties are excluded. Unused capillaries can be disposed of with the household waste. Unused capillaries must be disposed of in such a way that injuries to third parties are excluded.

No.	Description
9600105	RINGCAPS 1+2+3+4+5 µl CC, BOX = 250 PIECES
9600110	RINGCAPS 10µl CC, BOX = 250 PIECES
9600120	RINGCAPS 20 µl CC, BOX = 250 PIECES
9600122	RINGCAPS 200 µl CC, BOX = 100 PIECES
9600125	RINGCAPS 25 µl CC, BOX = 250 PIECES
9600140	RINGCAPS 40 µl CC, BOX = 250 PIECES
9600144	RINGCAPS 44,7 µl CC, BOX = 250 PIECES
9600150	RINGCAPS 50 µl CC, BOX = 250 PIECES
9600199	RINGCAPS 100 µl CC, BOX = 250 PIECES
9600210	RINGCAPS 5+10 µl CC, BOX = 250 PIECES
9600222	RINGCAPS 100+200 µl CC, BOX = 100 PIECES
9600240	RINGCAPS 20+40 µl CC, BOX = 250 PIECES
9600299	RINGCAPS 50+100 µl CC, BOX = 250 PIECES
9602144	RINGCAPS 44.7 µl, NA-HEP CC, BOX = 250 PIECES

**IVD** In vitro diagnostic medical devices

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