KCL- INFORMATIONSBROSCHÜRE

**Art. 890
897
898**

- Ⓧ Informationsbroschüre
- Ⓛⓖ Информационна брошюра
- ⓁⓏ Informační brožura
- ⓁⓀ Informationsbrochure
- Ⓛⓔ Folleto informativo
- ⓁⓔⓈⓐ Infobrošūr
- ⓁⓎ Tiedotuslehtisen
- Ⓛⓕ Brochure d'information
- Ⓛⓑ Information brochure
- ⓁⓖⓎ Ενημερωτικό φυλλάδιο
- ⓁⓎ Információs füzet
- Ⓛⓔ Informazioni brochure
- ⓁⓇ Informacije brošura
- Ⓛⓐ Informacija brošiūra
- ⓁⓋ Informācijas brošūra
- ⓁⓕⓁ Informatie brochure
- Ⓛⓓ Informasjonsbrosjyre
- ⓁⓔⓁ Broszura informacyjna
- Ⓛⓕ Folheto informativo
- ⓁⓇⓐ Informații broșură
- ⓁⓎ Информационная брошюра
- ⓁⓀ Informačná brožúra
- Ⓛⓔ Informacijska Brošura
- Ⓛⓔ Informationsbroschyr
- ⓁⓇ Bilgi Broşürü

KCL GmbH
Industriepark Rhön
Am Kreuzacker 9
36124 Eichenzell
Germany
Hotline: ++49-6659-87300
Fax: ++49-6659-87155
Web: www.kcl.de
<<http://www.kcl.de>>

by Honeywell

IIICBUT890

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by Honeywell

Information brochure in section 1.4 of Annex II of Directive 89/686/EEC

Personal chemical resistant protective gloves Cat III. Article 890, 897, 898

Protective Gloves conform to the EN 420. Approved to EN 374 and EN 388 according to Article 10 of Directive 89/686/EEC:

IFA, Alte Heerstraße 111, D-53757 St. Augustin, Identification number: 0121

SGS United Kingdom Ltd, Systems & Services Certification 2028 Worel Parkway, Weston-super-Mare, BS22 6WA, UK, Identification number: 0120

Surveillance after 11 B of Directive 89/686/EEC, chemical risks according to EN 374: same location as type test.

The identification number for testing and certification centre 0121 for EC type testing and monitoring of quality assurance measures relates solely to the content of EN374-1:2003 and the PPE Directive 89/686/EEC.

Protection against bacteriological contamination, EN 374

KCL chemical protective gloves have been penetration-tested to the highest performance level (3) set out in EN 374-2. This quality limit equates to an AQL < 0.65.

Full-protection protective gloves against chemicals, EN 374

The protective index refers to the permeation time determined during uninterrupted contact with the test chemical under stable laboratory conditions. EN 374-3 – Permeation

A glove is resistant to chemicals, when a protection index is achieved at least Level 2 in three of the chemicals listed below. The chemicals which have passed the test are marked on the gloves with the letters A-L. (KB=classification letters.)

KB	Chemical	CAS no.:
A	Methanol	67-56-1
B	Acetone	67-64-1
C	Acetonitrile	75-05-8
D	Dichloromethane	75-09-2

KB	Chemical	CAS no.:
F	Toluene	108-88-3
G	Diethylamine	109-89-7
H	Tetrahydrofuran	109-99-9
I	Ethylacetate	141-78-6

Penetration time in min.	Protection index
> 10	1
> 30	2
> 60	3
> 120	4
> 240	5
> 480	6

By Notified Body examined values

Article	Name	Size – Gr.	EN388	classification letters	Testing institute	Test chemical protection index									
						A	B	C	D	F	G	H	I		
890	Vitoject*	8, 9, 10, 11	3101	D F G	0121	4	0	1	4	6	4	1	1	4	
897	Butoject*	7, 8, 9, 10, 11	2010	B C I	0120	6	6	6	0	1	1	1	4		
898	Butoject	8, 9, 10, 11	0010	B C I	0121	6	6	6	1	1	1	1	4		

The protection index of the test chemicals was determined under stable laboratory conditions in KCL's own laboratory.

Mechanical risks, EN 388

1st digit Abrasion resistance

(min. 0; max. 4)

3rd digit Tear propagation resistance

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*Article 898 conforms to EN 421, radioactive contamination

Warning!

Resistance to the chemicals listed was determined under laboratory conditions and may be adversely affected by changes to physical properties such as temperature, abrasion, stretching, etc. When using highly corrosive chemicals, degradation is the most important factor in the choice of protective gloves.

These protective gloves provide no protection against extreme cold (< -5°C), heat (> 50°C), electricity.

Do not use in the vicinity of moving machine parts. Risk of being drawn into the machine.

Allergy Notice: All gloves may contain traces of Mercapto benzothiazole, Article 898, 897 in addition traces of Thiuram, Dithiocarbamate.

Article 898 may also contain traces of PAH, as well as Sulphenamide.

Storage /Transport: Flat, dry, dark, with no additional weight load in its original packaging, at a temperature of 5 °C - 25 °C.

Protect from sunlight and ozone.

Use: The protective gloves you wear must be of the correct size. Note that using undergloves may result in some usage restrictions.

Check the protective gloves for damage before you use them. Damaged protective gloves must not be used under any circumstances.

Prevent the penetration of pollutants over the edge of the glove.

Prevent carryover of contaminants located on the glove and the cross-contamination when removing gloves.

Cleaning: If possible under running cold water. If the gloves have been exposed to chemicals, follow the disposal instructions of the manufacturer of the chemicals. The specified gloves are not washable.

Expiry date: With proper storage, no reduction in performance levels within 36 months.

Disposal: If not contaminated with chemicals, gloves can be disposed of in the household waste.

If the gloves have been exposed to chemicals, follow the disposal instructions of the manufacturer of the chemicals.

For further information please contact.

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