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11/15

Instruction manual Digital probe thermometer Kat. Nr. 30.1018

- 1. Before you start using it**
 - Please make sure to read the instruction manual carefully.
 - Following and respecting the instructions in your manual will prevent damage to your instrument and loss of your statutory rights arising from defects due to incorrect use.
 - We shall not be liable for any damage occurring as a result of non following of these instructions.
 - Likewise, we take no responsibility for any incorrect readings and for any consequences which may result from them.
 - Please take particular note of the safety advice!
 - Please keep this instruction manual for future reference.

2. Field of operation

- Current temperature / MAX/MIN/HOLD / probe 125 mm / also for food checks (according to HACCP)

For your safety

- The device should only be used as described in your instruction manual.
- Keep this instrument and the battery out of the reach of children.
- Batteries can be fatal if swallowed. If a battery has been swallowed, get medical assistance immediately.
- Batteries must not be thrown in fire, short-circuited, taken apart or recharged. Danger of explosion!
- Batteries contain harmful acids. Low batteries should be changed as soon as possible to prevent damage caused by a leaking.
- Remove the battery if you do not use the instrument for a longer period of time.
- Do not place your product near extreme temperatures, vibrations or shocks. Protect it from moisture.
- Only the probe is heat resistant up to 200 °C.
- Never hold the probe directly over fire.
- Do not use the instrument in the oven or microwave.
- Unauthorized repairs, modifications or changes to the product are prohibited.
- Clean your instrument with a soft damp cloth. Do not use solvents or scouring agents.
- Do not immerse the display element into water. Water can penetrate into the precision sensor and cause malfunction. Not suitable for dishwasher.
- Be careful when handling the probe. When the device is not used, put the protective cap over the probe.

3. Getting started / Operation

- Open the battery compartment and remove the battery insulation strip. The device is ready for use.
- Remove the protective cap from the probe.
- Press the ON/OFF button, to turn on the device.
- The display shows the currently measured temperature. The sensor is located in the point of the probe.
- To measure the temperature, insert the probe at least 2,5 cm deep in the object.
- The temperature unit can be changed from °C to °F by the °C/°F button.
- Press the MODE button once in normal mode and the hold-function is activated (on the display appears HOLD).
- Press the MODE button again and the maximum temperature is shown since the last reset (on the display appears MAX).
- Press the MODE button again and the minimum temperature is shown since the last reset (on the display appears MIN).
- Press the MODE button again to return to normal mode.
- To clear the maximum and minimum memory, press and hold the MODE button for 3 seconds in the corresponding mode while MAX or MIN is displayed (indicated by ---). Press the MODE button again to return to normal mode.
- Press and hold MODE button for 3 seconds in normal mode, to deactivate or activate the auto-matic off-function.
- X appears or disappears on the display.
- Change the battery when the functions of the device become weak.
- Open the battery compartment and insert a new battery LR44 polarity as illustrated (+pole above).

4. Troubleshooting

Problems	Solutions
No display	→ Ensure battery polarity is correct (+ pole above) → Change the battery → Switch on the instrument (ON)
Incorrect display	→ Check the position of the sensor → Change the battery

5. Waste disposal

Never dispose of empty batteries and rechargeable batteries in household waste. As a consumer, you are legally required to take them to your retail store or to an appropriate collection site depending on national or local regulations in order to protect the environment.

The symbols for the heavy metals contained are:
Cd=cadmium, Hg = mercury, Pb=lead

This instrument is labelled in accordance with the EU Waste Electrical and Electronic Equipment Directive (WEEE).

Please do not dispose of this instrument in household waste. The user is obligated to take end-of-life devices to a designated collection point for the disposal of electrical and electronic equipment, in order to ensure environmentally-compatible disposal.

6. Specifications

Measuring range	-40°C to +200°C (-40°F to 392°F)
Operating temperature	-20°C to +50°C (-4°F to 122°F)
Precision	±0,8°C (±1,4°F) from 0°C ... +100°C (32°F ... 212°F), ±1,0°C (±1,8°F) from -20...0°C (-4°F...32°F) otherwise ±1,5°C (±2,7°F)
Resolution	0,1 °C (0,1 °F)
Power consumption	LR 44 button cell (included)
Auto switching-off function	after 5 minutes
Dimensions	20 x 16 x 205 mm
Weight	29 g (instrument only)

This product fulfills the guidelines according to EN 13485.

Suitability: S, T (Storage, Transport)

Location: E

Accuracy class: 1

Measuring range: -40°C...+200°C

In accordance with EN 13485, this instrument is subject to regular inspections as per EN 13486 (recommendation: yearly)

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TFA Dostmann GmbH & Co.KG, Zum Ottersberg 12, D-97877 Wertheim